

## Covid 19 Rapid Antigen Test, Clongene Covid 19 Antigen Rapid Test Diagnostic Colloidal Gold Antigen Test

Our Product Introduction

for more products please visit us on [oxygenconcentratormedical.com](http://oxygenconcentratormedical.com)

### Basic Information

- Place of Origin: China
- Brand Name: DVOT
- Certification: CE/BfArM/ PEI /TGA
- Model Number: Swab / Saliva
- Minimum Order Quantity: Negotiable
- Price: Negotiable
- Packaging Details: Each boxes 25tests packing /17.5 kg per carton. 64\*34\*55cm 1250 test per Carton.
- Delivery Time: 5-8 Working days
- Payment Terms: T/T, Western Union
- Supply Ability: 5000000Test



### Product Specification

- Product Name: Diagnostic Kit For COVID-19 Antigen Rapaid Test
- Format: Cassette
- Safety Standard: EN149-2001+A1-2009
- Warranty: 12months
- Keyword: Antigen Rapid Test
- Specimen: Whole Blood/ Serum/ Plasma
- Result Time: Read In 15 Minutes
- Accuracy: Over 99.4%
- After Sale Service: Online Technical Support
- Storage: 1~30° Normal Temperature
- Highlight: Clongene Covid 19 Antigen Rapid Test, Colloidal Gold Covid 19 Antigen Rapid Test, Diagnostic Colloidal Gold Antigen Test Kit



### More Images



## Product Description

### Immunochromatography Lateral Flow Test Kit ,SARS-CoV-2 Lateral Flow Test Kit Professional Use

#### Intend Use

This rapid test kit is intended for the qualitative detection of SARS-CoV-2 infection from patients. It is for professional use only. It is an aid in the diagnosis of the patients with suspected SARS-CoV-2 infection in conjunction with clinical presentation and results of other laboratory tests. Results from this test kit should not be used as the sole basis for diagnosis.






### Product Advantage

- Simplicity: It is extremely easy to use with simple to understand instructions.
- Fast: After the extraction of the sample, it will take just 10-15 mins for the result to be reflected on the testing device.
- Reliable: CE, ISO13485 ,TGA Certificated, Stric quality control, High quality
- Extraction tube is fully enclosed for disposal

 **CE Marked**

 **Instant result at 15 minutes**

 **Easy to collect samples**

 **Results are clearly visible**

 **No equipment required**

 **Suitable for large-scale rapid screening**

### Main Components

- 25 Test Cassettes
- 25 Sample tubes
- 1 Sample extraction buffers
- 25 Swabs
- 1 Package Insert





## PERFORMANCE CHARACTERISTIC

### 1. Clinical Verification

The performance of Test Card was established with 243 sample collected from symptomatic patients, who with symptoms onset within 7 days.

## Performance Characteristics

### Clinical Performance

#### Nasopharyngeal Swab

770 nasopharyngeal swabs were collected from individual symptomatic patients by using CLUNGENE® COVID-19 Antigen Rapid Test. The RT-PCR cycle threshold (Ct) is the relevant signal value. Lower Ct value indicate higher viral load. The sensitivity was calculated for the different Ct value range (Ct value ≤ 33 and Ct value ≤ 37).

Summary data as below:

COVID-19 Antigen		RT-PCR (Ct value ≤ 33)		Total
		Positive	Negative	
CLUNGENE®	Positive	145	2	147
	Negative	3	593	596
Total		148	595	743

PPA (Ct ≤ 33): 98.0% (145/148), (95%CI: 94.2%~99.3%)  
NPA: 99.7% (593/595), (95%CI: 98.8%~99.9%)

COVID-19 Antigen		RT-PCR (Ct value ≤ 37)		Total
		Positive	Negative	
CLUNGENE®	Positive	161	2	163
	Negative	14	593	607
Total		175	595	770

PPA (Ct ≤ 37): 92.0% (161/175), (95%CI: 87.0%~95.2%)  
NPA: 99.7% (593/595), (95%CI: 98.8%~99.9%)

PPA - Positive Percent Agreement (Sensitivity)

NPA - Negative Percent Agreement (Specificity)

#### Nasal Swab

617 nasal swabs were collected from individual symptomatic patients by using CLUNGENE® COVID-19 Antigen Rapid Test. The sensitivity was calculated for the different Ct value range (Ct value ≤ 33 and Ct value ≤ 37).

Summary data as below:

COVID-19 Antigen		RT-PCR (Ct value ≤ 33)		Total
		Positive	Negative	
CLUNGENE®	Positive	132	3	135
	Negative	4	462	466
Total		136	465	601

PPA (Ct ≤ 33): 97.1% (132/136), (95%CI: 92.7%~98.9%)  
NPA: 99.4% (462/465), (95%CI: 98.1%~99.8%)

COVID-19 Antigen		RT-PCR (Ct value ≤ 37)		Total
		Positive	Negative	
CLUNGENE®	Positive	139	3	142
	Negative	13	462	475
Total		152	465	617

PPA (Ct ≤ 37): 91.4% (139/152), (95%CI: 85.9%~94.9%)  
NPA: 99.4% (462/465), (95%CI: 98.1%~99.8%)

PPA - Positive Percent Agreement (Sensitivity)

NPA - Negative Percent Agreement (Specificity)

### Limit of Detection (Analytical Sensitivity)

The study used cultured SARS-CoV-2 virus which is heat inactivated and spiked into nasopharyngeal swab specimen. The Limit of Detection (LoD) is  $5.7 \times 10^2$  TCID<sub>50</sub>/mL.

### Cross Reactivity (Analytical Specificity)

32 commensal and pathogenic microorganisms that may be present in the nasal cavity have been evaluated and no cross-reactivity was observed.

### Interference

17 potential interference substances with different concentration were evaluated and found no affect to the test performance.

### High-dose Hook Effect

The COVID-19 Antigen Rapid Test was tested up to  $1.0 \times 10^{6.67}$  TCID<sub>50</sub>/mL of inactivated SARS-CoV-2 and no high-dose hook effect was observed.

## Certificates

## Our Certifications



## LAB



## Factory Workshop



## Why Choose Us?

10

FACTORY DIRECT SUPPLY,  
Professional manufacturer



We have many certificates  
CE ,ISO,TUV,TGA



Industry and trade  
Integration,more discount!



Overseas Warehouse in  
Europe,Fast shipment



## Packing & Delivery



Shenzhen Emeng Health Technology Co., LTD

+86 13760822077

kaylayi0612@foxmail.com

oxygenconcentratormedical.com

Room D, 7th Floor, Guang Long Building, 162 Ping Xin North Road, PingHu Street, LongGang District,  
ShenZhen