

Colloidal Gold Antigen Test, Influenza A/B SARS-CoV-2 Rapid Antigen Test Detection Kit Cassette

Our Product Introduction

for more products please visit us on oxygenconcentratormedical.com

Basic Information

- Place of Origin: China
- Brand Name: DVOT
- Certification: CE/BfArM/ PEI /TGA
- Model Number: Swab / Saliva
- Minimum Order Quantity: 1250Test
- Price: 0.8~1USD per Test
- Packaging Details: Each boxes 25tests packing /17.5 kg per carton. 64*34*55cm 1250 test per Carton.
- Delivery Time: 5-8 Working days
- Payment Terms: T/T, Western Union
- Supply Ability: 1000000pcs




Product Specification

- Product Name: Influenza A + B Rapid Detection Kit
- Formats: Cassette
- Specifications: Cassette(3mm, 4mm) Or According To The Customers' Requirements
- Reaction Time: 10-15 Minutes
- Shelf Life: 24 Months At Room Temperature 4- 30 Degrees
- Others: The Kits Can Be Made According To The Customers' Artwork Or Design
- Time Of Use: One Time Off Use Only
- Highlight: SARS-CoV-2 Rapid Antigen Test Cassette, SARS-CoV-2 Rapid Antigen Test, Influenza A/B Rapid Test Cassette



More Images



Product Description

Rapid answers, combined.

The SARS-CoV-2 & Flu A/B Rapid Antigen Test helps healthcare professionals at the point of care to quickly detect and differentiate infections with any of the three respiratory viruses and aid in pursuing appropriate courses of action, including patient management decisions. Also, it enables scaling up test capacities during flu season to meet the need for extensive testing during peak periods.

Our Product Introduction

to



Key benefits

- Reliable test performance
- Rapid answers in only 15 minutes
- Support in the differential detection of SARS-CoV-2, Influenza virus A, and Influenza virus B
- Easy handling, only 4 steps needed
- Fully integrated test, only one application area and single sample are needed
- Optimized composition to detect antigens against several dominant variants of the SARS-CoV-2
- Access to testing in areas where laboratory testing is not available
- Data matrix code on each test device for facilitated result sharing

Performing the test

Simple to perform and fast results within 15 – 30 min*

Collection of nasopharyngeal sample

Insert a sterile swab into the nostril of the patient. Rotate the swab 3 – 4 times against the nasopharyngeal surface, then leave the swab in place for 10 seconds to absorb secretions. Withdraw the swab from the nasal cavity.

Sample preparation

- Insert the swab into an extraction buffer tube. While squeezing the buffer tube, stir the swab at least 10 times.
- Remove the swab while squeezing the sides of the tube to extract the liquid from the swab.
- Press the nozzle cap tightly onto the tube.

WARNING! Failure to squeeze the tube can lead to incorrect results due to excess buffer in the swab.

Performing the test

Place the test device on a flat surface and apply 4 drops of extracted sample at a 90° angle to the specimen well of the test device.


Read the results

Read the test result at 15 to 30 minutes. Do not touch or move the test device until the result can be read.

WARNING! Test results that are read before 15 minutes or after 30 minutes may be incorrect.

* Before using the SARS-CoV-2 & Flu A/B Rapid Antigen Test, carefully read the instructions for use that are provided in the kit.

Results interpretation



Test kit information

Each kit contains 25 individually packaged, ready-to-use tests.

The following components are needed for a test and included in the kit:

- Test devices
- Nozzle cap
- Extraction buffer tube
- Buffer tube holder
- Sterile nasopharyngeal swabs
- 1 set of positive and negative quality control swabs
- Instructions for use
- Quick Reference Guide

10+ Years Experience

- OEM/ODM Available
- Factory Direct Price
- CE&ISO Marked
- Intimate service



Shenzhen Emeng Health Technology Co., LTD

+86 13760822077

kaylayi0612@foxmail.com

oxygenconcentratormedical.com

Room D, 7th Floor, Guang Long Building, 162 Ping Xin North Road, PingHu Street, LongGang District, ShenZhen