

Medical Grade Disposable fFN Rapid Test Kit Fetal Fibronectin Rapid Detection, Colloidal Gold Antigen Rapid Test

Our Product Introduction

for more products please visit us on oxygenconcentratormedical.com


Basic Information

- Place of Origin: China
- Brand Name: DVOT
- Certification: CE/BfArM/ PEI
- Model Number: fFN
- Minimum Order Quantity: 1250Test
- Price: 0.95USD per Test
- Packaging Details: Each boxes 25tests
- Delivery Time: 5-8 Working days
- Payment Terms: T/T, Western Union
- Supply Ability: 1000000pcs



Product Specification

- Product Name: FFN Rapid Test Kit
- Formats: Cassette
- Reaction Time: 10-15 Minutes
- Shelf Life: 24 Months At Room Temperature 4- 30 Degrees
- Others: The Kits Can Be Made According To The Customers' Artwork Or Design
- Time Of Use: Disposable
- Highlight: **Medical fFN Rapid Test Kit,
Disposable fFN Rapid Test Kit,
Medical Fetal Fibronectin Rapid Detection**



Product Description


Medical Grade Rapid Detection Of Fetal Fibronectin (fFN)

The fetal fibronectin (fFN) rapid test (private parts secretion) is a visually interpreted, qualitative immunochromatographic test device for detection of fFN in private parts secretions during pregnancy, which is a special protein that literally holds your baby in place in the womb. The test is intended for professional use to help diagnose the rupture of fetal membranes (ROM) in pregnant women.



Detection Principle

The fFN (private parts secretion) has been designed to detect fFN through visual interpretation of color development in the internal strip. The membrane was immobilized with anti-fFN antibodies on the test region. During the test, the specimen is allowed to react with colored anti-fFN antibodies colloidal gold conjugates, which were precoated on the sample pad of the test. The mixture then moves on the membrane by a capillary action, and interact with reagents on the membrane. If there were enough fFN in specimens, a colored band will form at the T region of the membrane. Presence of colored band indicates a positive result, while its absence indicates a negative result. Appearance of a colored band at the control region serves as a procedural control. This indicates that proper volume of specimen has been added and membrane wicking has occurred.



SPECIMEN COLLECTION AND STORAGE

The Fetal fibronectin (fFN) rapid test (private parts secretion) is intended only for use with women private parts secretions specimens.

The specimen is cervicoprivate parts secretion that is extracted into the Specimen Extraction Solution provided. A private parts secretion sample is obtained using a sterile polyester swab from the posterior fornix of the private parts during a sterile speculum examination or, if no private parts fluid is visible, the sample may be taken from the cervix. Take not to touch anything with the swab before taking the sample. The swab should be left in the private parts or cervix for approximately 10~15 seconds to allow it to absorb the secretion samples.

Open the Specimen Extraction Solution tube and put it in a vertical position. The specimen is extracted immediately from the swab by swirling the swab vigorously in the extraction solution for approximately 10 seconds. Specimens should be tested as soon as possible after extraction but in any case no more than 4 hours after specimen collection and extraction. If a specimen can not be tested within this time it should be frozen. After thawing, the specimens can be tested as described below.

Perform the testing immediately after the specimen collection. Do not leave the specimens at room temperature for prolonged periods. Specimens may be stored at 2- 8°C for up to 72 hours.

Bring specimens to room temperature prior to testing.

Pack the specimens in compliance with applicable regulations for transportation of etiological agents, in case they need to be shipped.

PROCEDURE

Bring tests, specimens, buffer and/or controls to room temperature (15-30°C) before use.

Remove the test from its sealed pouch, and place it on a clean, level surface. Label the test with patient or control identification. To obtain a best result, the assay should be performed within one hour.

Hold the label area and put the strip into the extracted buffer towards the arrow label. Leave the strip into the buffer until you see the liquid front enter the result area. Remove the strip from the solution and place it in a horizontal position.

As the test begins to work, you will see color move across the membrane.

Wait for the colored band to appear. The result should be read at 10 minutes. Do not interpret the result after 20 minutes.

INTERPRETATION OF RESULTS

NOTE:

The intensity of the color in test region (T) may vary depending on the concentration of aimed substances present in the specimen. Therefore, any shade of color in the test region should be considered positive. Besides, the substances level cannot be determined by this qualitative test.

Insufficient specimen volume, incorrect operation procedure, or performing expired tests are the most likely reasons for control band failure.

LIMITATIONS OF THE TEST

The Fetal fibronectin (fFN) rapid test (private parts secretion) is for professional in vitro diagnostic use, and should be used for the qualitative detection of fFN only.

As with all diagnostic tests, a definitive clinical diagnosis should not be based on the results of a single test, but should only be made by the physician after all clinical and laboratory findings have been evaluated.

If the test result is negative and clinical symptoms persist, additional testing using other clinical methods is recommended. A negative result does not at any time preclude the possibility

Preservation and stability of products

The kit should be stored at 2-30°C until the expiry date printed on the sealed pouch.

The test must remain in the sealed pouch until use.

Do not freeze.

Cares should be taken to protect components in this kit from contamination. Do not use if there is evidence of microbial contamination or precipitation. Biological contamination of dispensing equipments, containers or reagents can lead to false results.

KIT COMPONENTS

Individually packed test strips Each strip contains colored conjugates and reactive reagents pre-spread at the corresponding regions.

Specimens collection swab For specimens collection use.

Specimens dilution tube with buffer 0.1 M Phosphate buffered saline (PBS) and preservative

Package insert For operation instruction.

MATERIALS REQUIRED BUT NOT PROVIDED

Timer For timing use.



Shenzhen Emeng Health Technology Co., LTD



+86 13760822077



kaylayi0612@foxmail.com



oxygenconcentratormedical.com

Room D, 7th Floor, Guang Long Building, 162 Ping Xin North Road, PingHu Street, LongGang District,
ShenZhen