

HIV 1/2 Antibody Test Kit Rapid Antigen Test Kit for Home Use 99% Accurate 15-Minute Results

Our Product Introduction

For more products please visit us on oxygenconcentratormedical.com

Basic Information

- Place of Origin: china
- Brand Name: DVOT
- Certification: CE/ISO
- Model Number: HIV 1/2
- Minimum Order Quantity: 360 test
- Price: 0.8-0.9USD
- Packaging Details: 360T/PER CTN Weight :11KGS CTN SIZE :46*44*44CM
- Delivery Time: 15-20 work days
- Payment Terms: D/P, T/T, Western Union, MoneyGram
- Supply Ability: 1000000pcs



Product Specification

- Test Time: 15 Minutes
- Keyword: Antigen Rapid Test
- Result Time: 15 Minutes To Get
- Accuracy: 99%
- Time Of Use: Disposable
- Type: Rapid Test Kit
- Certifications: CE ISO
- Highlight: ce clungene antigen rapid test, ce medical diagnostic test kits, disposable clungene antigen rapid test



More Images



Product Description

HIV 1/2 Antibody Test Kit, Rapid Antigen Test Kit for Home Use 99% Accurate 15-Minute Results

Product Description:

Our test kits are certified by both CE and ISO, ensuring that you can trust the results you receive. With our Antigen Rapid Self Test Kit, you can easily test diseases, for example HIV.

Our HIV1-2 third generation Test Kits are designed for easy use and can be administered . The test time is just 15 minutes, making it a quick and convenient option for busy professionals and individuals alike.

Our Antigen Rapid Test has been proven to be 99% accurate, giving you confidence in the results you receive. With our test kits, you can quickly and easily identify diseases and take appropriate measures to protect yourself and those around you.



Technical Parameters:

Product Attribute	Description
Type	Rapid Test Kit
Keyword	Antigen Rapid Test
Accuracy	99%
Time Of Use	Disposable
Certifications	CE, ISO
Test Time	15 Minutes

Result Time	15 Minutes To Get
Other Keywords	Antibody Rapid Test Kit, COVID19 Infection Test Kit, Antigen Home Test Kit

Advantages:

Easy to Use
Reliable
Total Privacy
Accuracy

Applications:

One of the most common applications of this Antigen Rapid Self Test Kit is for testing HIV1-2, which is a crucial test for anyone who may have been exposed to the virus. The test kit is easy to use and provides an accuracy rate of 99%, making it a reliable option for professionals as well as individuals who want to test themselves at home.

In HIV1-2 testing, providing accurate results in a short amount of time. This makes it a valuable tool for healthcare professionals who need to quickly diagnose patients and develop treatment plans.

The HIV1-2 third generation Antigen Rapid Test Kit is manufactured in China, ensuring high standards of quality and reliability. Its small size and disposable nature make it a convenient option for various scenarios, including home testing kits and medical facilities. The kit features easy-to-use instructions and requires no additional equipment, making it a user-friendly option for anyone who needs to perform diagnostic tests.

Support and Services:

Our Medical Diagnostic Test Kits come with comprehensive technical support and services to ensure accurate and reliable results.

Steps to use:

Preparation stage:

Verify that the product packaging is intact and check the date of manufacture and lot number.

Prepare the required sample collection tools, such as a disposable blood collection needle (for whole blood samples), a matching transfer pipette or capillary pipette.

For serum/plasma sample testing, ensure that sufficient serum/plasma (approximately 20 ml) is collected.

Serum/Plasma Sample Test:

Remove the test strip and ensure that it is dry and undamaged.

Using the serum/plasma in the serum/plasma collection container, slowly drop approximately 20 ml of serum/plasma along the designated area on the test strip (usually the test line area).

Wait 15 minutes for the test results to appear naturally.

Whole blood sample test:

Wash and sterilize your finger to ensure that the test area is clean.

Using a disposable blood collection needle, prick the edge of the finger to collect a whole blood sample of approximately 25 µL.

Using a matched capillary or transfer pipette, add the collected whole blood sample to the designated area on the test strip.

It may be necessary to add 1 to 2 drops of buffer solution to facilitate the reaction, as indicated in the instructions.

Wait for the reaction:

Whether it is a serum/plasma sample or a whole blood sample, wait for the designated amount of time (usually 15 minutes) to allow the test results to appear naturally.



Reading the results:

Observe the result indicator area on the test strip.

If a line similar in color to the Control line appears, and a line of the same color appears at the Test line, the result is usually "Positive. Positive.

If the color line appears only at the Control line, but not at the Test line, the result is usually Negative.

If there are no lines at the Control line, or if the result is unclear and difficult to determine, then the result may be Invalid and the test needs to be repeated.





+86 13760822077



kaylayi0612@foxmail.com



oxygenconcentratormedical.com

Room D, 7th Floor, Guang Long Building, 162 Ping Xin North Road, PingHu Street, LongGang District,
ShenZhen