

15 Minutes Test Time Malaria Pf/Pan Antigen Rapid Test Medical Diagnostic Test Kits CE and ISO Certified

Our Product Introduction

For more products please visit us on oxygenconcentratormedical.com


Basic Information

- Place of Origin: china
- Brand Name: DVOT
- Certification: CE/ISO
- Model Number: Malaria Pf/Pan Rapid Test Kit
- Minimum Order Quantity: 360 test
- Price: 1-1.2USD
- Packaging Details: 360T/PER CTN Weight :11KGS CTN SIZE :46*44*44CM
- Delivery Time: 15-20 work days
- Payment Terms: D/P, T/T, Western Union, MoneyGram
- Supply Ability: 1000000pcs



Product Specification

- Time Of Use: Disposable
- Result Time: 15 Minutes To Get
- Type: Rapid Test Kit
- Keyword: Antigen Rapid Test
- Certifications: CE ISO
- Test Time: 15 Minutes
- Accuracy: 98%
- Highlight: medical diagnostic test kits ce, ce clungene antigen rapid test, disposable clungene antigen rapid test



More Images



Product Description

15 Minutes Test Time Malaria Pf/Pan Antigen Rapid Test Medical Diagnostic Test Kits CE and ISO Certified

Product Description:

With an accuracy rate of 98%, Antigen Rapid Test Kits are a reliable and effective way to diagnose patients quickly and accurately. These kits typically take only 15 minutes to complete, making them an ideal choice for busy healthcare facilities and clinics. Additionally, Antigen Rapid Test Kits are designed to be disposable, which means that they can be used once and then discarded, reducing the risk of contamination and infection in healthcare settings.

In addition to Antigen Rapid Test Kits, there are a variety of other Medical Diagnostic Test Kits available on the market today. For example, Antibody Rapid Test Kits are designed to detect the presence of specific antibodies in a patient's body, which can be an indicator of a recent or current infection. These kits are often used to diagnose illnesses like Malaria Pf/Pan.

Overall, Medical Diagnostic Test Kits are a crucial tool for healthcare providers and patients alike. Whether you need to diagnose the flu with Influenza Test Kits, detect the presence of antibodies with Antibody Rapid Test Kits, or predict ovulation with a Digital LH Test Kit, there is a test kit available to meet your needs. With their fast and accurate results, these test kits are helping to improve patient outcomes and reduce the spread of illness and infection in healthcare settings.

Technical Parameters:

Test Time	15 Minutes
Accuracy	98%
Time Of Use	Disposable
Type	Rapid Test Kit
Certifications	CE, ISO
Result Time	15 Minutes To Get
Keyword	Antigen Rapid Test

Applications:

The DOVT Malaria Pf/Pan Rapid Test Kit is a disposable antigen home test kit, making it easy to use and dispose of. The kit is designed to provide results in just 15 minutes, making it a quick and efficient way to diagnose malaria. The rapid test kit is suitable for use in various scenarios, such as medical clinics, hospitals, and even in remote areas where laboratory testing is not easily accessible.

Aside from malaria, this test kit can also be used for other infectious disease rapid test kits, such as COVID-19 infection test kits. The DOVT Malaria Pf/Pan Rapid Test Kit is a great tool to have in the medical field as it can be used to diagnose a range of infectious diseases quickly and accurately.

Overall, the DOVT Malaria Pf/Pan Rapid Test Kit is a valuable tool for medical professionals. Its accuracy, convenience, and speed make it a reliable option for diagnosing infectious diseases. Whether it is used in a medical clinic, hospital, or in remote areas where laboratory testing is not easily accessible, this rapid test kit is an excellent addition to any medical professional's toolkit.

Support and Services:

Our medical diagnostic test kits are designed to provide accurate and reliable results for a variety of medical conditions. We offer technical support and services to ensure that our customers are able to use our products effectively and efficiently.

Steps to use the product:

Prepare the tools:

- Make sure you have matching capillary pipettes (capillary pipette) and clean blood collection tools (e.g., blood collection needles, alcohol pads, etc.).

Collect a blood sample:

- Clean the tip of your finger (or other blood collection site) with an alcohol pad and allow the alcohol to evaporate before collecting blood.
- Collect an appropriate amount of blood sample (usually 25 microliters) by gently pricking the tip of the finger with a blood collection needle.
- Aspirate the blood sample using a matched capillary pipette.

Perform the test:

- Add 2 drops of the just drawn blood to the test card
- Then add the buffer solution

Read Results:

- At the end of the waiting time, observe the results area on the test card or container.

How to read test results:

Observe the results area:


- Look for the results interpretation area or color comparison chart according to the product description.
- Observe the results interpretation area on the test card or container for color changes, lines appearing, or other recognizable indications.

Compare Results:

- A "Positive" result indicates the presence of the substance or condition to be tested for in the blood sample and may require consultation with a physician on a case-by-case basis.
- If the result is "Negative", the substance or condition to be tested for is not detected in the blood sample.
- If the result is "Invalid", the test may need to be repeated due to improper handling, sample contamination or invalid test card/reagent.

Caution:

- Ensure that the time requirements in the product instructions are followed when reading the results.
- If there is any doubt or uncertainty about the results, it is advisable to consult a physician or relevant healthcare provider for professional advice.
- Blood collection kits and test cards/containers should be disposed of properly after use to avoid contamination of the environment.





Shenzhen Emeng Health Technology Co., LTD

+86 13760822077

kaylayi0612@foxmail.com

oxygenconcentratormedical.com

Room D, 7th Floor, Guang Long Building, 162 Ping Xin North Road, PingHu Street, LongGang District, ShenZhen