

## Accuracy Disposable Cholera O1/O139 Antigen Rapid Test Kits at home, Easy to Use

### Our Product Introduction

For more products please visit us on [oxygenconcentratormedical.com](http://oxygenconcentratormedical.com)

#### Basic Information

- Place of Origin: china
- Brand Name: DVOT
- Certification: CE/ISO
- Model Number: O1 Antigen/ O139 Antigen
- Minimum Order Quantity: 360 test
- Price: 0.9-1 USD
- Packaging Details: 360T/PER CTN Weight :11KGS CTN SIZE :46\*44\*44CM
- Delivery Time: 15-20 work days
- Payment Terms: D/P, T/T, Western Union, MoneyGram
- Supply Ability: 1000000pcs




#### Product Specification

- Result Time: 15 Minutes To Get
- Accuracy: 98%
- Keyword: Antigen Rapid Test
- Type: Rapid Test Kit
- Certifications: CE ISO
- Test Time: 15 Minutes
- Time Of Use: Disposable
- Highlight: disposable clungene antigen rapid test, disposable medical diagnostic test kits, ce clungene antigen rapid test



#### More Images



## Product Description

### Accuracy Disposable Cholera O1/O139 Antigen Rapid Test Kits at home,Easy to Use

#### Product Description:

Our Antibiotic Test Kit is CE and ISO certified, ensuring that you get the highest quality of testing available on the market. They are disposable, making them an affordable and convenient option for medical professionals. Our test kits are specifically designed to provide rapid results, allowing medical professionals to make quick diagnoses and provide timely treatment.

Our Infectious Disease Test Kit is perfect for diagnosing a wide range of infections, including Cholera O1/O139. Our Antigen Rapid Test is highly accurate, with a 98% accuracy rate, making it an effective tool for detecting infections in a timely manner. This can help in the early detection of the disease and prevent its spread.

Our Cholera O1/O139 Test Kits are designed to be user-friendly and provide quick results. They are easy to use, even for those without medical training, and can be used in a variety of settings, including hospitals, clinics, and even at home. With our test kits, individuals can get quick and reliable results, allowing them to take the necessary precautions to prevent the spread of the virus.

Our Rapid Test Kits are a cost-effective solution for medical professionals and individuals alike. They are designed to be easy to use, affordable, and reliable, making them an ideal choice for diagnosing various infections. With our test kits, you can rest assured that you will get accurate results every time.

#### Technical Parameters:

Product Name	Antigen Rapid Test Kit
Type	Rapid Test Kit
Time of Use	Disposable
Accuracy	98%
Result Time	15 Minutes to Get
Test Time	15 Minutes
Certifications	CE, ISO
Keywords	Antigen Rapid Test, Antigen Home Test Kit, Antibody Rapid Test Kit, Home Testing Kits

#### Applications:

The O1 Antigen/ O139 Antigen Rapid Self Test Kit is used for the detection of *Vibrio cholerae* in human feces. This test is commonly used in areas where cholera is endemic, and it is essential for early diagnosis and treatment. The test kit is easy to use and provides accurate results in a short amount of time, making it ideal for use in remote areas with limited access to medical facilities.

The Antibody Test Kit is used to detect the presence of antibodies in the blood, which can indicate a recent or past infection with a particular virus or disease. This test is useful for monitoring the spread of infectious diseases, identifying individuals who may have immunity to a particular virus, and identifying individuals who may have been exposed to a virus but are asymptomatic.

The Antigen Home Test Kit is designed for home use, allowing individuals to test themselves for the presence of a particular virus or disease. This test is ideal for individuals who are experiencing symptoms but are unable to visit a medical facility for testing. The test is easy to use and provides accurate results in a short amount of time, allowing individuals to take appropriate measures to protect themselves and others.

Overall, the DVOT medical diagnostic test kits are essential tools for early diagnosis and treatment of infectious diseases. They are easy to use, provide accurate results, and are ideal for use in remote areas or for individuals who are unable to visit a medical facility. With their certifications and high accuracy, these test kits provide a reliable and effective solution for detecting and monitoring the spread of infectious diseases.

#### Advantages:

- Easy to Use
- Reliable

- Total Privacy
- Accuracy
- 2 In One

## Support and Services:

Our Medical Diagnostic Test Kits come with comprehensive technical support and services to ensure the highest level of accuracy and reliability in results.

## Steps for product use:

### Prepare the sample:

- Make sure you have the type of sample you need and have prepared it according to the product instructions.
- Check that the sample container is clean and free of contamination.

### Place the sample into the container:

- Follow the first step instructions on the front of the box and carefully place the appropriate amount of sample into the designated container.

### Stir & Shake:

- First, stir 15 times to ensure the sample is well mixed.
- Then, shake 10 times; this step helps activate the test reagents or better mix the sample with the reagents.

### Drip the sample:

- Next, place three drops of the sample on the test card or at the designated spot in the designated area, as instructed in step 3 on the front of the box.

### Wait and read the results:

- Wait for the test card or reagent to react according to the wait time on the product instructions.
- After the designated reading time, view the test results.

### How to see test results:

Observe the results area:


- Follow the instructions on the back of the box to locate the results interpretation area or color comparison chart.
- Look carefully at the results display area on the test card or reagent for color changes, lines appearing, or other recognizable indications.

### Compare results:

- Compare the observed results with the result interpretation on the back of the box or on the instructions.
- If the result shows "Positive", the substance or condition to be detected is present in the sample.
- If the result is "Negative", the substance or condition to be detected is not present in the sample.
- If the result is "Invalid", it may be caused by improper operation, sample contamination or test card failure, etc., and the test needs to be repeated.

### Caution:

- When reading the results, please make sure to follow the time requirements on the product instructions, reading the results too early or too late may result in inaccuracies.
- If there is any doubt or uncertainty about the results, it is recommended to consult a physician or relevant healthcare provider for professional advice.





**Shenzhen Emeng Health Technology Co., LTD**

+86 13760822077

kaylayi0612@foxmail.com

oxygenconcentratormedical.com

Room D, 7th Floor, Guang Long Building, 162 Ping Xin North Road, PingHu Street, LongGang District, ShenZhen