

Disposable HIV 1/2 Antigen Rapid Test Kit Accurate Results in 15 Minutes

Our Product Introduction

For more products please visit us on oxygenconcentratormedical.com

Basic Information

- Place of Origin: china
- Brand Name: DVOT
- Certification: CE/ISO
- Model Number: HIV 1/2
- Minimum Order Quantity: 360 test
- Price: 0.8-0.9USD
- Packaging Details: 360T/PER CTN Weight :11KGS CTN SIZE :46*44*44CM
- Delivery Time: 15-20 work days
- Payment Terms: D/P, T/T, Western Union, MoneyGram
- Supply Ability: 1000000pcs



Product Specification

- Keyword: Antigen Rapid Test
- Type: Rapid Test Kit
- Time Of Use: Disposable
- Test Time: 15 Minutes
- Result Time: 15 Minutes To Get
- Accuracy: 99%
- Certifications: CE ISO
- Highlight: disposable clungene antigen rapid test, disposable medical diagnostic test kits, ce clungene antigen rapid test



More Images



Product Description

Disposable HIV 1/2 Antigen Rapid Test Kit Accurate Results in 15 Minutes

Product Description:

The Medical Diagnostic Test Kits product is designed to be used once and disposed of, making it a convenient and cost-effective option for many people. With a test time of just 15 minutes, you can quickly and easily get the results you need without having to wait for hours or days. Whether you are testing for a medical condition or simply want to ensure your health and safety, this product is an excellent choice.

The Medical Diagnostic Test Kits product is certified by both CE and ISO, ensuring that it meets the highest standards of quality and safety. This means that you can trust this product to provide accurate and reliable results every time you use it. Whether you are a healthcare professional or simply someone who wants to take charge of their health, this product is an excellent choice.

Technical Parameters:

Accuracy	99%
Certifications	CE, ISO
Keyword	Antigen Rapid Test
Time Of Use	Disposable
Test Time	15 Minutes
Result Time	15 Minutes To Get
Type	Rapid Test Kit

Applications:

The DVOT HIV 1/2 Antigen Rapid Test Kit is ideal for use in various medical settings, including hospitals, clinics, and laboratories. It is also a popular choice for home testing kits, as it is easy to use and provides fast and reliable results. This makes it an excellent choice for individuals who want to test themselves in the comfort of their own home without having to visit a doctor or a medical facility.

The DVOT HIV 1/2 Antigen Rapid Test Kit is suitable for use in a wide range of scenarios, including HIV screening programs, blood banks, and public health clinics. It is also useful for healthcare professionals who need to quickly diagnose HIV infections in their patients. This rapid test kit is also ideal for use in research laboratories that study HIV and other infectious diseases.

The DVOT HIV 1/2 Antigen Rapid Test Kit is a type of Rapid Test Kit, which means that it uses innovative technology to deliver fast and accurate results. This makes it an excellent choice for healthcare professionals who need to quickly diagnose HIV infections in their patients. It is also a popular choice for individuals who want to test themselves for HIV without having to wait for days or weeks for their results.

Support and Services:

Our Medical Diagnostic Test Kits come with comprehensive technical support and services to ensure accurate and reliable results.

Steps to use:

Preparation stage:

Verify that the product packaging is intact and check the date of manufacture and lot number.

Prepare the required sample collection tools, such as a disposable blood collection needle (for whole blood samples), a matching transfer pipette or capillary pipette.

For serum/plasma sample testing, ensure that sufficient serum/plasma (approximately 20 µl) is collected.

Serum/Plasma Sample Test:

Remove the test strip and ensure that it is dry and undamaged.

Using the serum/plasma in the serum/plasma collection container, slowly drop approximately 20 µl of serum/plasma along the designated area on the test strip (usually the test line area).

Wait 15 minutes for the test results to appear naturally.

Whole blood sample test:

Wash and sterilize your finger to ensure that the test area is clean.

Using a disposable blood collection needle, prick the edge of the finger to collect a whole blood sample of approximately 25 µL.

Using a matched capillary or transfer pipette, add the collected whole blood sample to the designated area on the test strip.

It may be necessary to add 1 to 2 drops of buffer solution to facilitate the reaction, as indicated in the instructions.

Wait for the reaction:

Whether it is a serum/plasma sample or a whole blood sample, wait for the designated amount of time (usually 15 minutes) to allow the test results to appear naturally.

Reading the results:


Observe the result indicator area on the test strip.

If a line similar in color to the Control line appears, and a line of the same color appears at the Test line, the result is usually "Positive. Positive.

If the color line appears only at the Control line, but not at the Test line, the result is usually Negative.

If there are no lines at the Control line, or if the result is unclear and difficult to determine, then the result may be Invalid and the test needs to be repeated.





Shenzhen Emeng Health Technology Co., LTD



+86 13760822077



kaylayi0612@foxmail.com



oxygenconcentratormedical.com

Room D, 7th Floor, Guang Long Building, 162 Ping Xin North Road, PingHu Street, LongGang District,
ShenZhen