

## 98% Accuracy Disposable H. Pylori Antigen Rapid Test Kit for Easy and Accurate Testing

Our Product Introduction

for more products please visit us on [oxygenconcentratormedical.com](http://oxygenconcentratormedical.com)

### Basic Information

- Place of Origin: China
- Brand Name: DVOT
- Certification: CE/ISO
- Model Number: H. Pylori
- Minimum Order Quantity: 360 test
- Price: 1.2USD
- Packaging Details: QTY: 360T/PER CTN Weight :11KGS CTN SIZE :46\*44\*44CM
- Delivery Time: 15-20 work days
- Payment Terms: D/P, T/T, Western Union, MoneyGram
- Supply Ability: 1000000pcs



### Product Specification

- Accuracy: 98%
- Certifications: CE ISO
- Test Time: 15mins
- Time Of Use: Disposable
- Keyword: Antigen Rapid Test
- Highlight: H. Pylori Antigen Rapid Test Kit, 98% Accuracy H. Pylori Antigen Test, Disposable H. Pylori Antigen Test



### More Images



## Product Description

### H. Pylori Antigen Rapid Test Kit for testing H.Pylori,accuracy 98% and easy to use

#### Product Description:

H. Pylori is a spiral-shaped, microaerobic bacterium that can survive on the gastric mucosa and cause infection. It is an important causative factor for gastritis, gastric ulcer and even gastric cancer. It is mainly transmitted through oral and fecal-oral transmission, sharing utensils, kissing, contaminated food, etc. In addition, it may also be transmitted through blood, mother-to-child transmission, etc., but it is not the main way. the symptoms of H. Pylori infection include indigestion, stomach pain, nausea and vomiting.

The H. Pylori Antigen Rapid Test Kit is a rapid test for the detection of *Helicobacter pylori* (H. pylori) antigen in humans. H. pylori is a bacterium that can survive in the highly acidic environment of the human stomach and is closely associated with a variety of gastric diseases, including gastritis, gastric ulcers, gastrointestinal tract infections and other gastric diseases. H. pylori is a bacterium that can survive in the highly acidic environment of the human stomach and is closely associated with a variety of gastric diseases, including gastritis, peptic ulcer and even gastric cancer. Therefore, accurate and rapid diagnosis of H. pylori infection is important for the prevention and treatment of related diseases.

#### Applications:

The easy-to-understand design allows users to complete the test without specialized training. Results are clear and intuitive and can be obtained in less than 8 minutes. A simple color change or indicator line determines whether you is infected with one of the H. Pylori viruses, making it easy to take action quickly.

In addition, it has a shelf life of 24 months and can be stored at without special storage conditions, making it convenient for home use and long-term storage. It ensures accurate and reliable testing and is commonly used in hospitals and clinics, homes and more. All reagents have been carefully screened and tested to ensure safety and harmlessness to people.

Product Name	Canine and Feline Giardia Antigen Test
Specimen	faeces
Time to Result	15mins
Shelf Life/Validity	24months
Accuracy	98%

#### Support and Services:

Our medical diagnostic test kits are designed for accuracy and ease of use. We offer technical support services to assist with any questions or issues that may arise during the testing process.

Our goal is to provide our customers with the highest level of support and service to ensure successful testing outcomes.



# H. Pylori Antigen Rapid Test Kit

## Why should I test for Helicobacter Pylori?

H. pylori infection can cause gastritis, peptic ulcer, gastric cancer, gastric lymphoma or iron deficiency anemia etc.

The potential of developing more severe illnesses from a H. pylori infection cannot be ruled out even if no symptom is observed. Discovering a H. pylori infection at an early stage and receiving treatment can lower the risk of potential stomach illnesses. It is important to eliminate H. pylori even if you recover from stomach illnesses as the H. pylori that remains in your stomach can lead to potential recurrence of illnesses.



## Features



**Painless stool sample  
collection**



**Results in 15 minutes**



**EU CE Certification**



**Self-test at home**



**For symptomatic &  
asymptomatic individuals**



**98% Accuracy**

## PERFORMANCE CHARACTERISTICS

DVOT kit	Commercial kit		Total
	Positive	Negative	
Positive	168	5	173
Negative	2	556	558
Total	170	561	731

Sensitivity of DVOT H. Pylori Antigen Rapid Test Kit: 98.82%

Specificity of DVOT H. Pylori Antigen Rapid Test Kit: 99.11%

## Test Procedure

**1**

Lift the toilet seat. Unfold the stool collection paper, tear the double-sided tapes on both ends of the stool collection paper, and stick it onto both sides of the toilet body.

**2**

Put down the toilet seat and carefully press the stool collection paper. Defecate onto the stool collection paper.

**3**

Remove the entire cap and use the spoon to take a spoonful sample (from 3 different positions off the stool, about 100mg) and add it into the stool collection bottle.

**4**

Close the entire bottle cap tightly and shake it in order to assure good sample dispersion.

**5**

Place the test device on a flat surface. Break off the top cap.

**6**

Hold the bottle vertically. Slowly squeeze 3 drops (80-100uL) of diluted sample into the sample well.

**7**

Leave for 15 minutes and read the result. Do not read the result after 20 minutes. Refer to the 'Interpretation of Results' section.

**15 Minutes Read Results**

**Interpretation of Results**

Positive

Negative

Invalid

## The packing details



CB	
QTY :	360 PCS
G.W :	12 KG
MEAS :	46X44X44CM



Shenzhen Emeng Health Technology Co., LTD

+86 13760822077

kaylayi0612@foxmail.com

oxygenconcentratormedical.com

Room D, 7th Floor, Guang Long Building, 162 Ping Xin North Road, PingHu Street, LongGang District, ShenZhen