

CE ISO Certified Disposable Test Kit for Rapid Detection of Cav Canine Adeno Virus in Pets

Our Product Introduction

for more products please visit us on oxygenconcentratormedical.com

Basic Information

- Place of Origin: China
- Brand Name: DVOT
- Certification: CE/ISO
- Model Number: CAV Ag
- Minimum Order Quantity: 360 test
- Price: 1.52USD
- Packaging Details: OPP Packing
- Delivery Time: 15-20 work days
- Payment Terms: D/P, T/T, Western Union, MoneyGram
- Supply Ability: 1000000pcs



Product Specification

- Accuracy: 99%
- Test Time: 8mins
- Certifications: CE ISO
- Time Of Use: Disposable
- Keyword: Antigen Rapid Test
- Highlight: **Pets Cav Disposable Test Kit ,
Canine Adeno Virus Disposable Test Kit ,
Rapid Detection Cav rapid Test Kit**



More Images



Product Description

Product Description:

CAV, known as Canine Adenovirus, or Canine Adenovirus, is a highly contagious virus that mainly affects the respiratory system, digestive system, and eyes of dogs. The virus is spread through direct or indirect contact, such as saliva, nasal secretions, urine and feces. When a dog is infected with CAV, it may show a diverse set of symptoms, such as coughing, diarrhea, and loss of appetite.

The Canine Adenovirus Antigen Test, a revolutionary rapid diagnostic tool in pet healthcare, is designed to accurately identify Canine Adenovirus (CAV) antigens in dogs. The test can be completed in a very short period of time, which not only ensures the accuracy of the test results, but also greatly shortens the diagnostic cycle, providing veterinarians with powerful diagnostic support for earlier identification and treatment of CAV-induced diseases.

Applications:

The easy-to-understand design allows users to complete the testing process without specialized training. Results are available in less than 8 minutes and are intuitively clear. A simple color change or indicator line will determine if your cat is infected with CAV, making it easy to take quick action.

In addition, the product has a shelf life of 24 months and can be stored at (2-30°C) without special storage conditions, making it convenient for home use and long-term storage. It ensures accurate and reliable testing and is commonly used in pet hospitals and clinics, pet owners' homes, pet boarding centers, etc. All reagents are rigorously screened and tested to ensure that they are safe and harmless to dogs and operators, giving you peace of mind when caring for your dog.

Product Name	Canine Adeno virus Antigen Test
Specimen	Dog Discharge
Time to Result	8mins
Shelf Life/Validity	24months
Storage	Between 2-30°C
Accuracy	99%

Support and Services:

Our medical diagnostic test kits are designed for accuracy and ease of use. We offer technical support services to assist with any questions or issues that may arise during the testing process.

Our goal is to provide our customers with the highest level of support and service to ensure successful testing outcomes.

Steps to use:

- 1, Sampling: Collecting samples from dog's secretions with swabs
- 2, Mixing: Insert the swab into the tube of assay buffer provided and mix thoroughly to obtain an effective sample leachate.
3. Place vertically and wait: Place the mixing test tube on a smooth surface and wait for 2 minutes as instructed. Do not shake or touch the test tube while waiting.
- 4, Drop into the test tube box: After capping the test tube, drop three drops of solution vertically into the appropriate position of the test tube box and wait for 8 minutes
- 5, Read the result: according to the color change of the test card and compare the result with the decision instructions to determine whether the dog is infected with the virus or not; if there is no specific band or abnormal color change, the test result is invalid.



PRODUCT CATALOGUE

KEEP YOUR PET HEALTHY



JUST RAPID TEST



Building solutions that facilitate living in harmony with nature

For the past 100 years, the world's population has increased by 1 billion every 13 years. As the population grows, ecosystem has been changed. The traditional buffer zones that separate humans from animals or from the pathogens that they harbor are notably reduced. And inevitably, the frequency and prevalence of these infectious diseases, defined as 'zoonosis', have increased.



The most effective way to minimize damages from zoonotic diseases is to detect them earlier.

We believe rapid tests (also known as lateral flow assays - LFAs - or immunochromatographic strip tests - ICTs), especially in the Point-of-Care Test (POCT), are crucial for managing zoonotic pandemics with the lowest cost.

Born in 2014, DVOT develops and produces high quality and highly efficient antigen-antibody reaction test kits that have shown consistent sensitivity, specificity, and accuracy.

And we will continue our effort to make faster and more accurate test kits by integrating advanced materials and technologies into our system.



Rapid diagnostic solution for

canine

Clinical signs	CDV	CIV	CPV	TOXO	CCV	RBV
Eye or nasal discharge	✓	✓				
Sneezing	✓	✓				
Coughing	✓	✓				
Poor appetite	✓	✓	✓	✓	✓	
Lethargic	✓	✓	✓	✓	✓	
Fever	✓	✓	✓	✓		
Oral ulcers						
Neurological signs	✓		✓	✓		
Abortions, stillbirths				✓		
Jaundice				✓		
Vomiting	✓		✓			
Diarrhoea	✓		✓		✓	
Hyperactivity, hydrophobia						✓

feline

Clinical signs	FCV	FHV	FPV	TOXO	FCOV	RBV
Eye or nasal discharge	✓	✓				
Sneezing	✓	✓				
Poor appetite	✓	✓	✓	✓		
Lethargic	✓	✓	✓	✓		
Fever	✓	✓	✓	✓		
Oral ulcers	✓					
Neurological signs			✓	✓		
Abortions, stillbirths				✓		
Jaundice				✓		
Vomiting			✓			
Diarrhoea			✓		✓	
Hyperactivity, hydrophobia						✓

Other available product list

Test Name	SPECIMEN	Pack size
Rabies Ab Rabies Antibody Test	whole blood/serum /plasma	1/5/25 test(s)/kit
RLN Pregnancy Relaxin Test	whole blood/serum /plasma	1/5/25 test(s)/kit
CPV Ab Canine Parvovirus Antibody Test	whole blood/serum /plasma	1/5/25 test(s)/kit
CDV Ab Canine Distemper Virus Antibody Test	whole blood/serum /plasma	1/5/25 test(s)/kit
CAV Ag Canine Infectious Hepatitis Antigen Test	ocular / nasal secretions	1/5/25 test(s)/kit
CAV I Ag Canine Adenovirus I Antigen Test	ocular / nasal secretions	1/5/25 test(s)/kit
CAV II Ag Canine Adenovirus II Antigen Test	ocular / nasal secretions	1/5/25 test(s)/kit
CPIV Ag Canine Parainfluenza Virus Antigen Test	ocular / nasal secretions	1/5/25 test(s)/kit
CIV Ag Canine Influenza Virus Antigen Test	ocular / nasal secretions	1/5/25 test(s)/kit
C.Bru Ab Canine Brucella Antibody Test	whole blood/serum /plasma	1/5/25 test(s)/kit
EHR Ab Ehrlichia Canis Antibody Test	whole blood/serum /plasma	1/5/25 test(s)/kit
FIV Ab Feline Immuno deficiency Virus Antibody Test	whole blood/serum /plasma	1/5/25 test(s)/kit

Test Name	SPECIMEN	Pack size
FeLV Ag Feline Leukemia Virus Antigen Test	whole blood/serum /plasma	1/5/25 test(s)/kit
LSH Ab Leishmania canis Antibody Test	whole blood/serum /plasma	1/5/25 test(s)/kit
BG Ab Babesia gibsoni Antibody Test	whole blood/serum /plasma	1/5/25 test(s)/kit
Lep Ab Leptospira Canis Antibody Test	whole blood/serum /plasma	1/5/25 test(s)/kit
Lyme Ab Lyme Disease Antibody Test	whole blood/serum /plasma	1/5/25 test(s)/kit



Shenzhen Emeng Health Technology Co., LTD



+86 13760822077



kaylayi0612@foxmail.com



oxygenconcentratormedical.com

Room D, 7th Floor, Guang Long Building, 162 Ping Xin North Road, PingHu Street, LongGang District, ShenZhen