

Digital Upper Arm/Wrist Electronic Blood Pressure Monitor Medical Automatic BP Machine

Our Product Introduction

for more products please visit us on oxygenconcentratormedical.com

Basic Information

- Place of Origin: China
- Brand Name: HANYUE
- Certification: CE/FCC/ Rohs /FDA
- Model Number: HY-01
- Minimum Order Quantity: 30pcs
- Price: 5USD ~ 15USD
- Packaging Details: OPP packing, 30PCS /Carton
- Delivery Time: 5-8 Working days
- Payment Terms: D/A, D/P, T/T, Western Union, MoneyGram
- Supply Ability: 100000pcs



Product Specification

- Type: Blood Pressure Monitor
- Display: Digital Liquid Crystal Display
- Product Weight: 470g
- Accessories Included: Warranty Card, Certificate, Wrist Band, Product Manual, Storage Bag
- Application: Medical Healthcare, Adult, Families And Adult, Hospital, Arm /obses
- Material: ABS, Nylon/PVC, Anodized Aluminum/zinc, Plastic, Nylon/ABS
- Highlight: **Upper Arm Blood Pressure Monitors, BP Blood Pressure Monitors, Medical Automatic BP Machine**

Product Description

Digital Upper Arm Blood Pressure Monitor Automatic BP Machine

An upper arm blood pressure monitor is a medical device that measures blood pressure using an inflatable cuff that is placed around the upper arm. The device works by inflating the cuff to temporarily stop the flow of blood through the brachial artery, and then slowly releasing the pressure to allow blood to flow again. The device then measures the pressure in the artery as the blood begins to flow through again.



Large LCD Display & Voice Function

Large LCD display can easily show you the readings and Read a your measurement at the same time. Friendly to the elders for home use.

One-button Operation

Measuring will start automatically once you pressing the START/STOP button. Bp Monitor Upper Arm Machine will provide you with accurate value within 1 minute.

Simple and convenient operation

Allows you to measure your heart health by pressing one button without connecting to a smartphone The Upper Arm Cuff itself can record and store up to 60 groups readings.

Comfortable Wide Range Cuff

20cm-42cm adjustable and soft wide range cuff is designed for a more comfortable measurement and better user experience.

Irregular Heartbeat Detector

Digital Blood Pressure Machine will automatically display and alert you (with your reading) to indicate the detection of an irregular heartbeat.

How can I measure my blood pressure to get accurate results?



1. Before turning on the device, please stay relaxed, then sit quietly and observe.
2. Place the cuff on the left arm above the elbow so that the edge is about 1 inch (2-3 cm) above the elbow joint and the tube is in the middle of the arm.
3. Raise the palm of your left hand and place your elbow on a table or flat surface. Sit upright and relax, do not move or talk, and breathe regularly.
4. Turn on the Blood Pressure Monitor, press the on/off button, stay relaxed and start the measurement.

The result will be displayed on the screen meanwhile voice broadcast.

*Why And When Do You Need To Choose An Upper Arm Blood Pressure Monitor ?

1, Old people whose vascular wall has poor flexibility. In this case, wrist blood pressure's result may be lower than the real. Because the wrist is the nerve endings, hypertensive patients may be complicated by arterial stenosis, blood clots and other vascular diseases, so the elderly use an Arm Blood Pressure Monitor that is more accurate than Wrist Blood Pressure Monitor!

2, Wrist blood pressure monitor is also accurate, but it is more suitable for young people, working people, travelling, gift. Because it is easy to carry.

What's in the box

1 x Electronic Blood Pressure Monitor

1 x Arm Band

1 x Storage Bag



Shenzhen Emeng Health Technology Co., LTD



+86 13760822077



kaylayi0612@foxmail.com



oxygenconcentratormedical.com

Room D, 7th Floor, Guang Long Building, 162 Ping Xin North Road, PingHu Street, LongGang District,
ShenZhen