



Automatic Blood Pressure Monitor Medical Automatic Upper Arm Blood Pressure Machine Arm Tensiometer

Our Product Introduction

for more products please visit us on oxygenconcentratormedical.com

Basic Information

- Place of Origin: China
- Brand Name: HANYUE
- Certification: CE/FCC/ Rohs /FDA
- Model Number: HY-01
- Minimum Order Quantity: 40pcs
- Price: 5USD ~ 15USD
- Packaging Details: OPP packing, 30PCS /Carton
- Delivery Time: 5-8 Working days
- Payment Terms: D/A, D/P, T/T, Western Union, MoneyGram
- Supply Ability: 100000pcs



Product Specification

- Type: Blood Pressure Monitor
- Measurable Arm Circumference: 22cm-32cm
- Power Supply: 4 AAA Batteries 1.5V
- Memory Function: 99sets
- Application: Medical Healthcare, Adult, Families And Adult, Hospital, Arm /obses
- Material: ABS, Nylon/PVC, Anodized Aluminum/zinc, Plastic, Nylon/ABS
- Highlight: medical blood pressure machine, medical blood pressure monitoring machine, adult blood pressure machine



More Images



Product Description

Medical Automatic Upper Arm Blood Pressure Machine Arm Tensiometer

Upper arm blood pressure monitors typically consist of a cuff that is wrapped around the upper arm, a pump that inflates the cuff, and a digital display that shows the blood pressure readings. Some devices also have additional features such as memory storage for multiple readings, irregular heartbeat detection, and the ability to indicate whether the blood pressure reading falls within normal range. These devices are commonly used in healthcare settings, as well as in home healthcare settings for patients who require regular monitoring of their blood pressure. They are also used by individuals who want to monitor their blood pressure at home for general health monitoring or to track the effectiveness of lifestyle changes or medication.

PRODUCT INFORMATION

INFRARED THERMOMETER

THE THERMOMETER IS A TYPE OF MAXIMUM THERMOMETER, ALSO KNOWN AS "MEDICAL THERMOMETER". GUANGDONG AND HONG KONG AND TAIWAN AREAS ARE CALLED HEAT PROBES. IT RECORDS THE HIGHEST TEMPERATURE EVER MEASURED BY THE THERMOMETER.



PRODUCT NAME: ELECTRONIC SPHYGMOMANOMETER

PRODUCT MODEL: X506

PULSE PARAMETERS: (40-160)

BLOOD PRESSURE (0-280)mmHg/
MEASUREMENT PARAMETERS: (0-37.33)kPa

POWER SUPPLY: 1.5V*2

WORKING CONDITIONS: -10°C-55°C

MEASUREMENT TIME: 8S

SIZE: 80X67X55MM

It is important to note that while upper arm blood pressure monitors can provide a useful tool for monitoring blood pressure, they are not a substitute for medical advice or treatment. If you have concerns about your blood pressure or other health issues, you should consult with a healthcare professional.

Large LCD Display & Voice Function

Large LCD display can easily show you the readings and read a your measurement at the same time. Friendly to the elders for home use.

One-button Operation

Measuring will start automatically once you pressing the START/STOP button. Bp Monitor Upper Arm Machine will provide you with accurate value within 1 minute.

Simple and convenient operation

Allows you to measure your heart health by pressing one button without connecting to a smartphone. The Upper Arm Cuff itself can record and store up to 60 groups readings.

Comfortable Wide Range Cuff

20cm-42cm adjustable and soft wide range cuff is designed for a more comfortable measurement and better user experience.

Irregular Heartbeat Detector

Digital Blood Pressure Machine will automatically display and alert you (with your reading) to indicate the detection of an irregular heartbeat.

How can I measure my blood pressure to get accurate results?



1. Before turning on the device, please stay relaxed, then sit quietly and observe.
2. Place the cuff on the left arm above the elbow so that the edge is about 1 inch (2-3 cm) above the elbow joint and the tube is in the middle of the arm.
3. Raise the palm of your left hand and place your elbow on a table or flat surface. Sit upright and relax, do not move or talk, and breathe regularly.
4. Turn on the Blood Pressure Monitor, press the on/off button, stay relaxed and start the measurement.

The result will be displayed on the screen meanwhile voice broadcast.

*Why And When Do You Need To Choose An Upper Arm Blood Pressure Monitor ?

1, Old people whose vascular wall has poor flexibility. In this case, wrist blood pressure's result may be lower than the real. Because the wrist is the nerve endings, hypertensive patients may be complicated by arterial stenosis, blood clots and other vascular diseases, so the elderly use an Arm Blood Pressure Monitor that is more accurate than Wrist Blood Pressure Monitor!

2, Wrist blood pressure monitor is also accurate, but it is more suitable for young people, working people, travelling, gift. Because it is easy to carry.

What's in the box

- 1 x Electronic Blood Pressure Monitor
- 1 x Arm Band
- 1 x Storage Bag



Shenzhen Emeng Health Technology Co., LTD



+86 13760822077



kaylayi0612@foxmail.com



oxygenconcentratormedical.com

Room D, 7th Floor, Guang Long Building, 162 Ping Xin North Road, PingHu Street, LongGang District, ShenZhen