

# Wireless Blood Pressure Monitors 6V For Patients Bluetooth WiFi Technology

## **Basic Information**

• Place of Origin:	China
Brand Name:	HANYUE
Certification:	CE/FCC/ Rohs /FDA
Model Number:	HY-01
Minimum Order Quantity:	40pcs

- Price: 5USE
- Packaging Details: OPF
- Delivery Time: 5-8 Working da
- Payment Terms:
- Supply Ability:

	HY-01
ntity:	40pcs
	5USD ~ 15USD
	OPP packing, 30PCS /Carton
	5-8 Working days
	D/A, D/P, T/T, Western Union, MoneyGram
	100000pcs



### **Product Specification**

- Prouduct Name:
- Accuracy:
- Working Voltage:
- Accessories Included:
- Application:
- Material:
- Highlight:

ded: Warranty Card, Certificate, Wrist Band,Product Manual, Storage Bag Medical Healthcare, Adult, Familes And Adult, Hospital, Arm /obses ABS, Nylon/PVC, Anodized Aluminum/zinc, Plastic, Nylon/ABS

**Blood Pressure Monitor** 

±3mmHg(±0.4kPa)

6V

Wireless Blood Pressure Monitors, Bluetooth Blood Pressure Monitors

#### Wireless Blood Pressure Monitors Offer Greater Convenience And Mobility For Patients

Wireless blood pressure monitors are becoming increasingly popular due to their convenience and mobility. Unlike traditional blood pressure monitor which are tethered to a machine and require a cord to connect to a cuff, wireless blood pressure monitors can be used anywhere without the need for corr or cables.

Wireless blood pressure monitors work by using Bluetooth or WiFi technology to connect the monitor to a smartphone, tablet, or other device. The device displays the blood pressure readings and may also provide additional features such as data tracking, reminders, and alerts. Some wireless blood pressure monitors may also come with their own app, which can provide additional features and functionality.

### How can I measure my blood pressure to get accurate results?



1.Before turning on the device, please stay relaxed, then sit quietly and observe.

2.Place the cuff on the left arm above the elbow so that the edge is about 1 inch (2-3 cm) above the elbow joint and the tube is in the middle of the arm.

3. Raise the palm of your left hand and place your elbow on a table or flat surface. Sit upright and relax, do not move or talk, and breathe regularly.

4. Turn on the Blood Pressure Monitor, press the on/off button, stay relaxed and start the measurement.

The result will be displayed on the screen meanwhile voice broadcast.

+86 13760822077

## \*Why And When Do You Need To Choose An Upper Arm Blood Pressure Monitor ?

1, Old people whose vascular wall has poor flexibility. In this case, wrist blood pressure's result may be lower than the real. Because the wrist is the nerve endings, hypertensive patients may be complicated by arterial stenosis, blood clots and other vascular diseases, so the elderly use an Arm Blood Pressure Monitor that is more accurate than Wrist Blood Pressure Monitor!

2, Wrist blood pressure monitor is also accurate, but it is more suitable for young people, working people, travelling, gift. Because it is easy to carry.

What's in the box

1 x Electronic Blood Pressure Monitor

1 x Arm Band

1 x Storage Bag

Shenzhen Emeng Health Technology Co., LTD

kaylayi0612@foxmail.com Covygenconcentratormedical.com

Room D, 7th Floor, Guang Long Building, 162 Ping Xin North Road, PingHu Street, LongGang District, ShenZhen