



Blood Pressure Wrist Monitor: Easy Operation & Clear LCD Display for Portable Digital Monitoring

Our Product Introduction

for more products please visit us on oxygenconcentratormedical.com

Basic Information

- Place of Origin: China
- Brand Name: HANYUE
- Certification: CE/FCC/ Rohs /FDA
- Model Number: BP-S10
- Minimum Order Quantity: 30pcs
- Price: 5USD ~ 15USD
- Packaging Details: OPP packing, 30PCS /Carton
- Delivery Time: 5-8 Working days
- Payment Terms: D/A, D/P, T/T, Western Union, MoneyGram
- Supply Ability: 100000pcs




Product Specification

- Type: Blood Pressure Monitor
- Working Voltage: 6V
- Shut Off The Current: <20uA
- Automatic Shutdown Time: 1 Mins
- Product Weight: 470g
- Box Volume: 12cmX10cmX15.4cm
- Accessories Included: Warranty Card, Certificate, Wrist Band, Product Manual, Storage Bag
- Highlight: **Blood Pressure Monitors Easy Operation, Digital Blood Pressure Monitors**



More Images



Product Description

Automatic Blood Pressure Monitor Upper Arm Portable BP Monitor

Large LCD Display & Voice Function

The blood pressure meter is equipped with a large, clear LCD display that provides easy-to-read measurements. The display typically shows systolic and diastolic blood pressure readings, as well as pulse rate. The large screen ensures that the readings are visible and accessible, even for individuals with visual impairments.

One-button Operation

The blood pressure meter is designed to automatically inflate and deflate the cuff, making the measurement process simple and convenient. The device typically includes a one-touch operation, allowing users to initiate the measurement with a single button press. Automatic operation ensures accurate and consistent readings.

Simple and convenient operation

Allows you to measure your heart health by pressing one button without connecting to a smartphone. The Upper Arm Cuff itself can record and store up to 60 groups readings.

Comfortable Wide Range Cuff

20cm-42cm adjustable and soft wide range cuff is designed for a more comfortable measurement and better user experience.

Irregular Heartbeat Detector

Digital Blood Pressure Machine will automatically display and alert you (with your reading) to indicate the detection of an irregular heartbeat.

How can I measure my blood pressure to get accurate results?



How can I measure my blood pressure to get accurate results?

1. Before turning on the device, please stay relaxed, then sit quietly and observe.
2. Place the cuff on the left arm above the elbow so that the edge is about 1 inch (2-3 cm) above the elbow joint and the tube is in the middle of the arm.
3. Raise the palm of your left hand and place your elbow on a table or flat surface. Sit upright and relax, do not move or talk, and breathe regularly.
4. Turn on the Blood Pressure Monitor, press the on/off button, stay relaxed and start the measurement.

The result will be displayed on the screen meanwhile voice broadcast.

*Why And When Do You Need To Choose An Upper Arm Blood Pressure Monitor ?

1, Old people whose vascular wall has poor flexibility. In this case, wrist blood pressure's result may be lower than the real. Because the wrist is the nerve endings, hypertensive patients may be complicated by arterial stenosis, blood clots and other vascular diseases, so the elderly use an Arm Blood Pressure Monitor that is more accurate than Wrist Blood Pressure Monitor!

2, Wrist blood pressure monitor is also accurate, but it is more suitable for young people, working people, travelling, gift. Because it is easy to carry.

What's in the box ?

- 1 x Electronic Blood Pressure Monitor
- 1 x Arm Band
- 1 x Storage Bag



OEM/ODM Multiple Custom

Logo / Brand



Package



Optional Cuff Size

Standard: 22-32cm

Kids Size: 17-22cm

LX Size: 22-42cm

LXX Size: 22-48cm



Shenzhen Emeng Health Technology Co., LTD

+86 13760822077

kaylayi0612@foxmail.com

oxygenconcentratormedical.com

Room D, 7th Floor, Guang Long Building, 162 Ping Xin North Road, PingHu Street, LongGang District,

