

Adults And Children Portable Asthma Nebulizer Quiet And Efficient Atomizer

Our Product Introduction

for more products please visit us on oxygenconcentratormedical.com

Basic Information

- Place of Origin: China
- Brand Name: HANYUE
- Certification: CE/FDA/ROHS/FCC
- Model Number: HY-08
- Minimum Order Quantity: 40pcs
- Price: 8~12USD
- Packaging Details: OPP Packing
- Delivery Time: 5-8 Working days
- Payment Terms: D/A, D/P, L/C, T/T, Western Union, MoneyGram
- Supply Ability: 10000pcs

三. Product picture display



Features:
Effective treatment, Low residual fluid, Low noise design, DC Power, Water quantity inspection, Food and materials



Product Specification

- Vibration Frequency: 120kHz+/-10%
- Power Supply Mode: Removable Battery
- Power Source: Electric, DC 3V, AC110-240V 50/60Hz
- Working Noise: ≤21dB
- Nebulization Rate: ≥0.2mL/min
- Safety Standard: GB/T18830-2009
- Size: 142*112*60 Mm
- Highlight: **Adults Portable Asthma Nebulizer, Children Portable Asthma Nebulizer**

Dyspnea
With this nebulizer it can help you to ease breathing difficulties



More Images



Product Description

Quiet And Efficient Atomizer Suitable For Adults And Children

Portable, Quiet, Efficient & Reliable

The APOWUS portable nebulizer is perfect for children and adults to use anywhere, anytime. This pocket-sized portable nebulizer is not only easy to carry, but also features a combination of micro-holes and ultrasonic mesh technology for quiet and effective relief of various breathing problems.



FOR TRAVELLING AND HOUSEHOLD USE



Adult mask



Child mask



Mouthpiece

Why are you choose Our nebulizer?

This portable nebulizer is not only small size easy to carry, but also features a combination of micro-holes and ultrasonic mesh technology for quiet and effective relief of various breathing problems.

MESH NEBULIZER HELP YOU HAVE BETTER HEALTH

Portable & Effective

- Quiet Working
- Modern Design
- Food-grade Material
- One-button Operation

Dimension Display.

Enjoy benefits from Our nebulizer

We are well aware that the health of the respiratory system is very important, especially for people with breathing problems. This is why our team is committed to researching better nebulizer. Our portable nebulizer has 2800 micro-holes, stable chip control and three adjustable modes to provide you with a soothing mist. Effectively release discomfort and help you breathe freely anytime, anywhere.

User-friendly design and a full range of accessories

The nebulizer has three light reminders which allows you to view the current operating mode. The one-button disassembly design and self-cleaning mode are convenient for daily cleaning. Accessories are provided for children and adults. 360-degree use enabled for any positions: lying or sitting down, especially for children when they feeling weak.

Operation Process



How to use a Handheld Portable Nebulizer ?

Prepare the Atomizer: Ensure that the atomizer is clean and free from any residue or liquid. If it has been used before, clean it thoroughly and allow it to dry.

Fill the Atomizer: Unscrew or remove the top of the atomizer to access the reservoir. Fill the reservoir with the liquid you want to spray. Be careful not to overfill it to avoid leakage when you screw the top back on.

Prime the Atomizer: Some handheld atomizers require priming before use. This involves pumping the spray mechanism a few times to create pressure inside the bottle. Check the instructions that came with your atomizer to see if priming is necessary, and follow the specific instructions provided.

Adjust the Nozzle: Most handheld atomizers have an adjustable nozzle that allows you to control the spray pattern. You can usually twist the nozzle to adjust between a fine mist and a more concentrated stream. Determine the desired spray pattern based on your needs.

Hold the Atomizer: Hold the atomizer in a comfortable grip, ensuring that your hand is not obstructing the nozzle or spray area.

Aim and Spray: Point the nozzle in the direction you want to spray. Depress the spray button or trigger to release the liquid. Move your hand in a sweeping motion or in short bursts, depending on the desired coverage and the type of liquid being sprayed.

Clean and Store: After using the atomizer, clean it thoroughly to prevent clogging or contamination. Rinse the reservoir with clean water, and if necessary, use a small brush to remove any residue. Allow the atomizer to dry completely before storing it in a cool, dry place.

Your trustworthy nebulizer

We committed to provide high efficient, fast, and easy-use nebulization for everyone. We serve the concept of becoming the guardian of everyone's health, and bring a healthy respiratory system to you and your family.



QTY: 40 pcs
N.W: 8.6 KG
G.W: 9.4 KG
SIZE: 572x297x262mm



N.W: 91.2 G
G.W: 214.5 G
SIZE: 142x112x60mm



Shenzhen Emeng Health Technology Co., LTD

+86 13760822077

kaylayi0612@foxmail.com

oxygenconcentratormedical.com

Room D, 7th Floor, Guang Long Building, 162 Ping Xin North Road, PingHu Street, LongGang District, ShenZhen