

## Portable Accurate Digital Sphygmomanometer With Voice Prompt Large Screen LCD Display's Digital BP Meter

Our Product Introduction

for more products please visit us on [oxygenconcentratormedical.com](http://oxygenconcentratormedical.com)

### Basic Information

- Place of Origin: China
- Brand Name: HANYUE
- Certification: CE/FCC/ Rohs /FDA
- Model Number: BP-S01
- Minimum Order Quantity: 30pcs
- Price: 5USD ~ 15USD
- Packaging Details: OPP packing, 30PCS /Carton
- Delivery Time: 5-8 Working days
- Payment Terms: D/A, D/P, T/T, Western Union, MoneyGram
- Supply Ability: 100000pcs




### Product Specification

- Display: Digital Liquid Crystal Display
- Automatic Shutdown Time: 1 Mins
- Accuracy:  $\pm 3\text{mmHg} (\pm 0.4\text{kPa})$
- Power Source: Battery \*4 Or Type-c
- Type: Blood Pressure Monitor
- Measurable Arm Circumference: 22cm-32cm
- Working Voltage: 6V
- Highlight: **Portable Digital Sphygmomanometer, Accurate Digital Sphygmomanometer, LCD Display Portable Sphygmomanometer**



### More Images



### Product Description

## Portable Large Screen LCD Display With Voice Prompts For Convenient And Fast Digital Blood Pressure Monitoring



### Large LCD Display & Voice Function

Large LCD display can easily show you the readings and Read a your measurement at the same time. Friendly to the elders for home use.

### One-button Operation

Measuring will start automatically once you pressing the START/STOP button. Bp Monitor Upper Arm Machine will provide you with accurate value within 1 minute.

### Simple and convenient operation

Allows you to measure your heart health by pressing one button without connecting to a smartphone The Upper Arm Cuff itself can record and store up to 60 groups readings.

### Comfortable Wide Range Cuff

20cm-42cm adjustable and soft wide range cuff is designed for a more comfortable measurement and better user experience.

### Irregular Heartbeat Detector

Digital Blood Pressure Machine will automatically display and alert you (with your reading) to indicate the detection of an irregular heartbeat.

## How can I measure my blood pressure to get accurate results?

Before turning on the device, please stay relaxed, then sit quietly and observe.

Place the cuff on the left arm above the elbow so that the edge is about 1 inch (2-3 cm) above the elbow joint and the tube is in the middle of the arm.

Raise the palm of your left hand and place your elbow on a table or flat surface. Sit upright and relax, do not move or talk, and breathe regularly.

Turn on the Blood Pressure Monitor, press the on/off button, stay relaxed and start the measurement.

The result will be displayed on the screen meanwhile voice broadcast.



### Why And When Do You Need To Choose An Upper Arm Blood Pressure Monitor ?

Old people whose vascular wall has poor flexibility. In this case, wrist blood pressure's result may be lower than the real. Because the wrist is the nerve endings, hypertensive patients may be complicated by arterial stenosis, blood clots and other vascular diseases, so the elderly use an Arm Blood Pressure Monitor that is more accurate than Wrist Blood Pressure Monitor!

Wrist blood pressure monitor is also accurate, but it is more suitable for young people, working people, travelling, gift. Because it is easy to carry.



### Why do you need a large LCD digital display voice function blood pressure monitor?

**Convenience:** An automatic blood pressure monitor allows you to measure your blood pressure at home or on the go, providing the convenience of monitoring your health without needing to visit a medical facility. The portable design makes it easy to carry the monitor with you and use it whenever needed.

**Accuracy:** Upper arm blood pressure monitors are generally considered more accurate than wrist or finger monitors. They tend to provide more reliable readings because the upper arm is closer to the heart and the major arteries, offering a better reflection of your true blood pressure.

**User-Friendly Interface:** The large LCD display on the monitor makes it easier to read and interpret the blood pressure readings. The clear and visible numbers ensure that you can quickly and accurately understand your blood pressure measurements, promoting better health monitoring.

**Voice Function:** The voice function further enhances the usability of the blood pressure monitor, especially for individuals with visual impairments or those who prefer auditory feedback. The monitor can audibly announce the blood pressure readings, making it accessible to a wider range of users.

**Tracking and Monitoring:** Many automatic blood pressure monitors come with memory storage capabilities, allowing you to track and monitor your blood pressure over time. This feature enables you to identify patterns, detect any changes, and share the data with your healthcare provider for better management of your blood pressure.

**Early Detection of Health Issues:** Regular monitoring of blood pressure using an automatic monitor can help in the early detection of hypertension (high blood pressure) or hypotension (low blood pressure). By identifying any abnormal readings, you can take appropriate measures, such as lifestyle changes seeking medical advice, to manage your blood pressure and reduce the risk of associated health issues.

**Improved Health Management:** With the ability to monitor your blood pressure at your convenience, you can be more proactive in managing your health. By tracking your readings regularly, you can make informed decisions about your lifestyle choices, including diet, exercise, stress management, and medication adherence, to maintain optimal blood pressure levels.

### What's in the box?

1 x Electronic Blood Pressure Monitor

1 x Arm Band

1 x Storage Bag

## Measure Your Blood Pressure At Home



Your Healthy Family Advisor





**Shenzhen Emeng Health Technology Co., LTD**

+86 13760822077

kaylayi0612@foxmail.com

oxygenconcentratormedical.com

Room D, 7th Floor, Guang Long Building, 162 Ping Xin North Road, PingHu Street, LongGang District,

