

#### Portable Blood Pressure Meter: Medical Grade Digital Blood Pressure Meter **Voice Prompt LCD Display Wrist BP Meter**

#### **Basic Information**

. Place of Origin: China • Brand Name: HANYUE

· Certification: CE/FCC/ Rohs /FDA

• Model Number: X506 Minimum Order Quantity: 30pcs

• Price: 5USD ~ 15USD

 Packaging Details: OPP packing, 30PCS /Carton

• Delivery Time: 5-8 Working days

• Payment Terms: D/A, D/P, T/T, Western Union, MoneyGram

• Supply Ability: 100000pcs



#### **Product Specification**

80 X 67 X 55 Mm • Box Volume:

• Display: Digital Liquid Crystal Display

• Operating Temperature: -10°C To 55°C Accuracy: ±3mmHg(±0.4kPa) **Blood Pressure Monitor** . Type:

• Power Supply: 2\* 1.5V

• Highlight: Medical Digital Blood Pressure Meter,

Wrist Digital Blood Pressure Meter,

**Medical Wrist BP Meter** 



#### More Images



#### Wrist Type Voice Prompt LCD Display Medical Grade Digital Blood Pressure Meter

## The new experience is small and convenient

# A sphygmomanometer that fits into a pocket

The machine body is small and light, with a thoughtful portable box, regardless of business trip Take blood pressure anytime and anywhere when traveling, working, etc



#### Large LCD Display & Voice Function

Large LCD display can easily show you the readings and Read a your measurement at the same time. Friendly to the elders for home use.

#### **One-button Operation**

Measuring will start automatically once you pressing the START/STOP button. Bp Monitor Upper Arm Machine will provide you with accurate value within 1 minute.

#### Simple and convenient operation

Allows you to measure your heart health by pressing one button without connecting to a smartphone The Upper Arm Cuff itself can record and store up to 60 groups readings.

#### **Comfortable Wide Range Cuff**

20cm-42cm adjustable and soft wide range cuff is designed for a more comfortable measurement and better user experience.

#### **Irregular Heartbeat Detector**

Digital Blood Pressure Machine will automatically display and alert you (with your reading) to indicate the detection of an irregular heartbeat.

#### How can I measure my blood pressure to get accurate results?

Positioning: Place the cuff of the sphygmomanometer on your wrist. The cuffs should fit snugly, but not too tightly Activating the device: Turn on the sphygmomanometer by pressing the power button. The LCD display will typically show the current d and time or indicate that it is ready for measurement.

Start measuring: Most wrist digital blood pressure monitors have auto-inflation. Press the start button and the cuff will begin to inflate. When the cuffs are tightened, you may feel a slight pressure on your wrists.

Blood pressure measurement: Once the cuff is inflated, the device will begin measuring your blood pressure. It uses an oscilloscope method to detect the pulsation of blood vessels. As the device measures your blood, the cuff will gradually deflate.

Show results: After the measurement is complete, the sphygmomanometer will show your systolic blood pressure (higher number) and diastolic blood pressure (lower number) on the LCD display. Some devices may also show your pulse rate.

Record results: If needed, you can record measurements manually or use any of the device's built-in memory or data storage features track blood pressure readings

## **MEASURE POSTURE CORRECTLY**

PLEASE MEDITATE FOR 5 MINUTES BEFORE THE TEST TO RELAX YOUR BODY, WRISTS AND FINGERS



1. THE BANDAGE SHALL BE CORRECTLY WOUND AND NOT TOO TIGHT



**POSTURE** 



2. CHOOSE A RELAXED 3. THE BANDAGE IS AT THE SAME HEIGHT AS THE HEART

### THREE-STEP QUICK **MEASUREMENT**

1.5 seconds to start



Wrap the sphygmomanometer tothe left wrist







button to start measuring blood pressure

#### Why And When Do You Need To Choose An Upper Arm Blood Pressure Monitor?

Old people whose vascular wall has poor flexibility. In this case, wrist blood pressure's result may be lower than the real. Because the wrist is the nerve endings hypertensive patients may be complicated by arterial stenosis, blood clots and other vascular diseases, so the elderly use an Arm Blood Pressure Monitor that is mon accurate than Wrist Blood Pressure Monitor!

Wrist blood pressure monitor is also accurate, but it is more suitable for young people, working people, travelling, gift. Because it is easy to carry.

Why do I need your wrist mounted electronic digital display blood pressure monitor?

Portability and Convenience: Wrist blood pressure monitors are compact and lightweight, making them highly portable. They are ideal individuals who need to monitor their blood pressure on the go or while traveling. Their small size also makes them easy to store and carry around.

Ease of Use: Wrist monitors are generally easy to use and require minimal setup. They typically have a simple interface with clear digit displays, making it easy to read and interpret the blood pressure readings. Many models also feature large buttons and intuitive contro enhancing user-friendliness.

Comfortable and Non-Invasive: Wrist monitors are less invasive than traditional upper-arm monitors, as they don't require a cuff to be wrapped around the upper arm. Instead, they use a cuff that wraps around the wrist, which can be more comfortable for some individue especially those with larger arms or mobility issues.

Quick and Convenient Measurements: Wrist monitors often provide quick and convenient measurements. They usually have automate inflation and deflation features, which means you can obtain blood pressure readings with a push of a button. This can save time and effort compared to manual inflation and deflation methods.

Memory and Data Storage: Many wrist blood pressure monitors come with built-in memory and data storage capabilities. They can sto multiple readings over time, allowing users to track their blood pressure trends and share the data with their healthcare providers for better monitoring and management of their condition.

Hypertension Indicators and Alerts: Some wrist monitors have additional features such as hypertension indicators and irregular heartbedetection. These indicators can alert users to potential high blood pressure or irregular heart rhythms.

#### What's in the box?

1 x Electronic Blood Pressure Monitor

SIZE: 80X67X55MM


1 x Storage Bag

# PRODUCT INFORMATION INFRARED THERMOMETER PRODUCT NAME: ELECTRONIC SPHYGMOMANOMETER PRODUCT MODEL: X506 PULSE PARAMETERS: (40-160) **BLOOD PRESSURE** (0-280)mmHg/ MEASUREMENT PARAMETERS: (0-37.33)kPa POWER SUPPLY: 1.5V\*2 WORKING CONDITIONS: -10°C-55°C **MEASUREMENT TIME:8S**









# One touch measurement voice guidance Mom and Dad's Pocket Doctor

From power on to the analysis of measurement results, one click operation is required to complete the whole process





**VOICE GUIDANCE** 



**USB CHARGING** 



REMINDER OF TIMED PRESSURE MEASUREMENT



**HD BACKLIGHT** 



MEMORY VALUE OF 2 \* 99 GROUPS



**COMPACT** AND CONVENIENT



+86 13760822077

kaylayi0612@foxmail.com

e oxygenconcentratormedical.com

Room D, 7th Floor, Guang Long Building, 162 Ping Xin North Road, PingHu Street, LongGang District, ShenZhen