

Hbot Hard Type Hyperbaric Chamber 1.3 Ata For Medical/Therapeutic Applications And Benefits

Our Product Introduction

For more products please visit us on oxygenconcentratormedical.com

Basic Information

- Place of Origin: SHENZHEN EMENG HEALTH
- Brand Name: O2EMENG
- Certification: CE/FDA/FC/ROHS/UKCA/SGS
- Model Number: EMY03
- Minimum Order Quantity: 1PC
- Price: \$23499
- Packaging Details: Wooden packing
- Delivery Time: 15-20 work days
- Payment Terms: D/P, T/T, Western Union, MoneyGram
- Supply Ability: 100pcs




Product Specification

- Temperature Range: Cooling Adjustable
- Weight: 900kgs
- Maximum Working: 1500w
- Maximum Pressure: 1.2-1.3ATA
- Sound Level: Greater Than 43dB
- Chamber Oxygen Concentration: 28% 32%
- Highlight: medical hyperbaric oxygen chamber, medical oxygen hyperbaric chamber, therapeutic hyperbaric oxygen chamber



More Images



Product Description

O2EMENG hbot hard type hyperbaric chamber 1.3 ata for Medical/Therapeutic Applications and Benefits

Product Description:

Designed for maximum safety and comfort, the Hyperbaric Chamber is equipped with advanced features that ensure a constant pressure is maintained throughout the treatment. This feature is essential to ensure that patients receive the full benefits of hyperbaric oxygen therapy.

Voltage: 220V, 50Hz9, Rated power: 1500W, Oxygen concentration in combustion chamber: 28% to 32%

One of the most significant advantages of the hyperbaric chamber is its portability. The hyperbaric chamber can be easily transported and installed anywhere, making it the best choice for patients who require regular hyperbaric oxygen therapy. Its portability also makes it ideal for hospitals and clinics that need to provide hyperbaric oxygen therapy to a large number of patients.

All in all, the hyperbaric chamber is the perfect solution for anyone who needs to receive hyperbaric oxygen therapy safely and effectively. With its impressive size, advanced features, and portability, the hyperbaric chamber is sure to provide patients with the best possible treatment experience. What are you waiting for? Invest in a hyperbaric chamber today and experience the benefits of hyperbaric oxygen therapy for yourself!

Applications:

The O2EMENG Hyperbaric Oxygen Chamber is an ideal product for individuals who suffer from a variety of medical conditions, including chronic pain, fatigue, and stress. The hyperbaric chamber can also be used to treat a range of other medical conditions, including respiratory disorders, stroke, and traumatic brain injury.

The portable hyperbaric chamber is perfect for individuals who want to enjoy the benefits of oxygen therapy from the comfort of their own home. The chamber is easy to set up and can be used in any room of the house. It is also easy to transport, which makes it ideal for individuals who want to take it with them when they travel.

The O2EMENG Hyperbaric Oxygen Chamber is also suitable for use in a range of healthcare settings, including hospitals, clinics, and rehabilitation centers. The chamber is designed to be easy to use.

Meanwhile this oxygen chamber has the following functions:

1.Enhance the speed of recovery:

Hyperbaric oxygen chambers can help patients recover quickly from fatigue, illness, etc. O2EMENG Hyperbaric Oxygen Chambers may accelerate the speed of recovery by stimulating aerobic metabolic processes in cells by increasing the level of oxygen in the blood and the distance over which the blood oxygen diffuses.

2.Immune Enhancement:

Hyperbaric Oxygen Therapy can help the elderly and the physically weak to increase blood oxygen levels and strengthen the body. In conclusion, the O2EMENG Hyperbaric Oxygen Chamber Model EMY03 is a versatile and high-quality product that is suitable for a wide range of applications. Whether you are looking to improve your health and wellbeing or are seeking a product to use in a healthcare setting, the O2EMENG Hyperbaric Oxygen Chamber is an excellent choice.

Customization:

Brand Name: O2EMENG

Model Number: EMY03

Place of Origin: Shenzhen Emeng Health Technology Co., LT-D(O2EMENG)located in Shenzhen

Rated Voltage: 220V

Material:Reclining seated with hardware material

Size: 220*150*190CM

Sound Level ≥ 43dB

Support and Services:

The Hyperbaric Oxygen Chamber product comes with comprehensive technical support and services to ensure that customers have a seamless experience with our product. Our technical support team is available to assist with any queries or issues that may arise during installation, setup, or usage of the product. We also provide regular maintenance and repair services to keep the product in optimal

condition and extend its lifespan. Our services include:

1.Demand customization: according to your use environment and expectations, recommend the optimal core products and life aesthetics composition scheme.

2.Exclusive design: the advanced technology team and the top design team jointly, from the size and color material matching, the actual program customization landing.

3. Craftsmanship Manufacturing: The core components and other materials of hyperbaric chamber series products are manufactured in aerospace material factories.

4. Lifetime maintenance: provide perfect professional after-sales service guarantee

We are committed to providing the highest level of support and services to our customers to ensure their satisfaction with our Hyperbaric Oxygen Chamber product.

Reasons to choose us:

1. Professionalism and authority:

As a manufacturer or brand of hyperbaric chamber, O2EMENG has rich industry experience and technology accumulation, and can provide professional and reliable hyperbaric oxygen therapy equipment.

Its products may have undergone stringent quality control and certification and comply with relevant industry standards to ensure safety and effectiveness in use.

2. Technological advancement:

O2EMENG hyperbaric chambers may adopt the latest technology and design concepts, such as air pressurization or pure oxygen pressurization technology, as well as advanced air supply and exhaust (oxygen) systems, air conditioning systems and control systems.

These technologies can ensure that the human body can safely and effectively inhale a high concentration of oxygen in a high-pressure environment for therapeutic purposes.

3. User feedback and reputation:

When choosing hyperbaric chambers, user feedback and reputation are also important reference factors, and O2EMENG hyperbaric chambers may be widely praised for their good treatment effect and user experience.



CHAMBER CAPSULE DESCRIPTION



1. Automatic digital display real-time temperature, humidity, oxygen purity and pressure.
2. Rotary timer system, use of time can be pressed in advance Multiple noise reduction technology.
3. Excellent Silence performance, the noise is $\leq 63\text{dB}$. Offer you a peaceful experience environment.
4. Double pressure relief device Spiral step less super smooth rotary valve. Spiral ergonomic hand-held design. Pressure relief is smoother and more convenient. Internal and external two-way integrated polycarbonate valve. One valve, two-way operation.
5. Automatic pressurize relief valves
Two automatic high precision constant pressure valve. Ensure the product is safe and effective by Aluminum alloy steel, high hardness, wear resistance.
6. Communicating system.
Intercom telephone is convenient for users to talk with operator outside.
7. Machine start-up protection system.
The chamber only can be turned on when the chamber door is completely closed, avoiding the danger of the door being left unlocked.
8. Automatic pressurize and depressurize system.
Based on the principle of "Cylinder (engine)" , Pressurize automatically after the door is completely closed.
9. Made by Aluminum alloy steel; With environmental spray paint;
Strong weight bearing capacity, Durable and wear-resistant; non-deformable; long lifespan.
10. In case of unexpected power failure, the door can be opened by emergency safety from inside the chamber.



SHIPPING AND DELIVERY



FAQ

1. How does hyperbarics work?

Hyperbarics is a technology in which the air pressure in the environment is increased. When a person's body is placed in a more pressure environment, it absorbs more oxygen molecules per volume of compressed air.

The body normally transports oxygen via the hemoglobin of the red blood cells. By increasing the air pressure, oxygen is then driven into the body's fluids, allowing a super-saturation of tissues and organs with oxygen. The increased pressure infuses the body with oxygen, even reaching injuries with damaged circulation.

2. What are the health benefits of hyperbaric oxygen chambers?

Studies have shown that hyperbaric oxygen chambers can treat hypoxic diseases, relieve fatigue, delay aging, prevent cardiovascular disease, improve sleep quality, accelerated wound repair and enhance immunity. We have customers serviced shows hyperbaric oxygen therapy makes very positive results.

3. How long is the treatment?

Usually, we recommend our clients receive the HBOT therapy 5-6 sessions per week, 60-90mins per session. Using frequency or duration can be adjusted according to individual circumstances.

4. Is HBOT safe?

Hyperbaric oxygen therapy is prescribed by a physician and performed under medical supervision. Although there are minor risks like all medical treatments, overall hyperbaric oxygen therapy is extremely safe.

5. Is there anyone who can not use our chamber?

Pregnant ladies, patients who catch a cold or Flu, or those who have rhinitis or otitis media are not recommended to use the chamber.

6. How do you feel inside a hyperbaric oxygen chamber?

When lying in the oxygen chamber, you will feel comfortable and at ease. Also, You may feel a little squeezing of the ears because of the pressure. It is quite normal and can be relieved by chewing or swallowing.

7. After bought the chamber, do I need to by any other things?

No, the price of our chamber includes the chamber capsule, all machines, namely the air compressor, oxygen concentrator, air cooler, and accessories like frame, tubes, etc, so you don't have to buy anything else.

8. Is it complex to install the chamber or change the filter?

Of course not, we will provide you with a video to show you how to install the chamber and change the filter.

We also support online services to help you with all the issues, so don't worry.

9. What is the general time warranty expected for ZHIKANG oxygen chambers?

We will provide one-year warranty and lifelong after-sales services.

If any quality problem/fault in material/design is under correct operation within one year,

If it's easy to fix, we will send new components for free and teach you how to repair them.


(Even after the one-year warranty period, you can also contact us freely and we will be glad to help you.)



O2EMENG

After Sale Service

WhatsApp: +86 13760822077
E-mail: kaylayi0612@foxmail.com
Hotline: +86 13760822077
Official website: o2oxygen.com



Shenzhen Emeng Health Technology Co., LTD



+86 13760822077



kaylayi0612@foxmail.com



oxygenconcentratormedical.com

Room D, 7th Floor, Guang Long Building, 162 Ping Xin North Road, PingHu Street, LongGang District,
ShenZhen