


Personalized Hyperbaric Oxygen Chamber Treatment For Rehabilitation Therapy

Our Product Introduction

For products please visit us on oxygenconcentratormedical.com


Basic Information

- Place of Origin: SHENZHEN EMENG HEALTH
- Brand Name: O2EMENG
- Certification: CE/FDA/FC/ROHS/UKCA/SGS
- Model Number: EMS04
- Minimum Order Quantity: 1PC
- Price: \$3750
- Packaging Details: wooden packing
- Delivery Time: 15-20 work days
- Payment Terms: D/P, T/T, Western Union, MoneyGram
- Supply Ability: 100PCS



Product Specification

- Properties: Rehabilitation Therapy Supplies
- Function: Hyperbaric Oxygen Therapy
- Pressure Reduction Time: 5-10 Mins
- Usage: Medical/Therapeutic
- Noise Level: Less Than 63 DB
- Model: EMS04
- Oxygen Concentration: $\geq 90\%$ (can Be Set)
- Working Pressure: 1.3-1.5ATA (can Be Set)
- Highlight: Personalized Hyperbaric Oxygen Chamber, Rehabilitation Hyperbaric Oxygen Chamber, Oxygen Hyperbaric Chamber Treatment



More Images



Product Description

O2EMENG Personalized Hyperbaric Oxygen Chamber for Rehabilitation Therapy

Product Description:

The working pressure of this Hyperbaric Oxygen Chamber ranges from 1.3-1.5ATA, which can be set according to the patient's needs. This range of pressure is ideal for hyperbaric oxygen therapy and ensures that patients receive the most effective treatment possible. The model of this chamber is EMS04, which is a model that has been designed with the latest technology and research in mind. This model is ideal for rehabilitation therapy supplies and is perfect for medical professionals who want to provide the best possible care for their patients.

The Hyperbaric Oxygen Chamber is used for hyperbaric oxygen therapy, which is a medical treatment that involves breathing pure oxygen while inside a pressurized chamber. This treatment is used for a wide range of conditions, including decompression sickness, infections, wounds that won't heal, and more. Hyperbaric oxygen therapy is also used to promote healing and recovery after surgery or injury. The Hyperbaric Oxygen Chamber is designed to deliver this therapy in a safe and effective manner.

The total weight of this Hyperbaric Oxygen Chamber is 83kgs. This makes it easy to move and transport as needed. Despite its lightweight design, the chamber is sturdy and durable, ensuring that it can withstand frequent use and provide reliable results for patients. The Hyperbaric Oxygen Chamber is an excellent investment for medical professionals who want to provide their patients with the best possible care.

Features:

Product Name: Hyperbaric Oxygen Chamber

Function: Hyperbaric Oxygen Therapy

Usage: Medical/Therapeutic

Properties: Rehabilitation Therapy Supplies

Working pressure: 1.3-1.5ATA (can Be Set)

Pressure reduction time: 5-10 Mins

Keywords: Hyperbaric oxygen treatment chamber, Hyperbaric oxygen treatment chamber, High-pressure oxygen chamber

Technical Parameters:

Function	Hyperbaric Oxygen Therapy
Usage	Medical/Therapeutic
Machine CTN size	127*55*105cm
Pressure reduction time	5-10 Mins
Noise Level	Less Than 63 DB
Model	EMS04
Working pressure	1.3-1.5ATA (can Be Set)
Oxygen concentration	≥90%(can Be Set)
Total Weight	83kgs
Properties	Rehabilitation Therapy Supplies

Applications:

1. Rehabilitation Therapy: The O2EMENG EMS04 is a valuable tool for rehabilitation therapists who work with patients recovering from a variety of injuries and illnesses. Hyperbaric therapy is known to accelerate the healing process and reduce inflammation, making it an ideal treatment for many conditions. The EMS04 model can be set to a working pressure of 1.3-1.5ATA, which is ideal for most rehabilitation applications.

2. Sports Medicine: Athletes who suffer from injuries or chronic pain can benefit from hyperbaric therapy in an oxygen-rich environment. The O2EMENG EMS04 Hyperbaric Oxygen Chamber provides the ideal conditions for this type of treatment, allowing athletes to recover faster and get back to their training routines more quickly.

3. Wellness Centers: Many wellness centers and spas offer hyperbaric oxygen therapy as a way to promote relaxation and stress relief. The EMS04 model is easy to use and can be set up quickly, making it an ideal option for these types of facilities.

4. Hospitals and Clinics: Hyperbaric therapy is often used in hospitals and clinics to treat a variety of conditions, including

decompression sickness, carbon monoxide poisoning, and diabetic wounds. The O2EMENG EMS04 is a high-quality hyperbaric chamber of oxygen that can be used in these settings to provide effective treatment to patients.

5. Home Use: The O2EMENG EMS04 Hyperbaric Oxygen Chamber is also suitable for home use. With its compact size and easy-to-use controls, it can be set up in a spare room or garage, providing a convenient way for individuals to receive hyperbaric therapy on their own schedule.

No matter the occasion or scenario, the O2EMENG EMS04 Hyperbaric Oxygen Chamber is a reliable and effective solution for hyperbaric therapy. With its impressive features and high-quality construction, this hyperbaric chamber of oxygen is sure to meet the needs of a wide range of users. The machine CTN size is 127*55*105cm, and pressure reduction time is 5-10 Mins.

Support and Services:

1.Demand customization: according to your use environment and expectations, recommend the optimal core products and life aesthetics composition scheme.

2.Exclusive design: the advanced technology team and the top design team jointly, from the size and color material matching, the actual program customization landing.

3. Craftsmanship Manufacturing: The core components and other materials of hyperbaric chamber series products are manufactured in aerospace material factories.

4. Lifetime maintenance: provide perfect professional after-sales service guarantee

Our Hyperbaric Oxygen Chamber product comes with comprehensive technical support and services to ensure optimal performance and customer satisfaction. Our team of experienced technicians and engineers are available to assist with installation, maintenance, and troubleshooting of the product.

In addition, we offer training sessions for customers to learn how to properly operate and maintain the Hyperbaric Oxygen Chamber. This includes training on safety protocols and best practices for usage.

Our technical support and services also include regular maintenance and repair services to keep the product functioning at its best. We offer both on-site and remote support options to better serve our customers.

If you have any questions or concerns about our technical support and services for the Hyperbaric Oxygen Chamber, please do not hesitate to contact us.

O2EMENG

Hyperbaric Oxygen Chamber

HYPERBARIC OXYGEN CHAMBER


- IMPROVE BRAIN BLOOD CIRCULATION
- SPEED UP METABOLISM
- HIGH PRESSURE OXYGEN INHALATION



WHY ARE MORE AND MORE PEOPLE CHOOSING HYPERBARIC OXYGEN CHAMBERS?



In some parts of the body, the capillaries are so narrow that blood cells cannot pass through them, bound oxygen cannot reach them, and they can only rely on dissolved oxygen. The oxygen generator has little change in the dissolved oxygen in the blood, which is also the reason why the ordinary oxygen absorption effect in the plateau area is general.



SOFT TYPE HYPERBARIC OXYGEN CHAMBER



Innovative technology,quality assurance



Portable, Small Footprint
Fast pressurization



Made of multi-layer TPU composite
environmentally friendly material



Fresh Air System
Achieve cabin air circulation effect





CONCLUSION



It is more feasible to increase the dissolved oxygen in the blood.
The hyperbaric oxygen chamber can maintain the oxygen partial pressure within a safe and effective range by flexibly adjusting the oxygen concentration or pressure, and gently and effectively improve the oxygen delivery to tissues and organs by increasing dissolved oxygen.




INTEGRATION OF SUPERCHARGING AND OXYGENATION



Overturn the complexity of traditional household



New hyperbaric
oxygen equipment



Oxygen equipment
+
air power system
+
Air Cooler

Intelligent all in one machine, no need to
operate three devices,
easy one-button operation saves time and effort.

FIVE FUNCTIONS OF THE ALL-IN-ONE MACHINE



01

Oxygen Production

02

Pressurization

03

Dehumidification

04

Refrigeration

05

Negative Ion Sterilization



One-button Start

Oxygen Flow Meter

304 Stainless Steel Body

Medical Grade Filter

Medical Grade
Humidification Bottle

360° Universal Wheel

PRODUCT DETAIL



Fresh Air System

The all-in-one machine continuously pressurizes to ensure air circulation in the chamber.



≥95%±3% Medical Oxygen Purity

In the chamber, the oxygen concentration will rise to 93%.



Emergency Pressure Relief Valve

An emergency valve lets you get out of the chamber quickly in an urgent situation.



Two Zippers

Two zippers let people enter the chamber easily and provide better sealing.



Two Side Pressure Gauge

Two side pressure gauges make it easier and safer to monitor the pressure both from inside and outside.



Easy operating by one people without assistant


PRODUCT PARAMETER

Product Model	EMS04 Soft Type Hyperbaric Oxygen Chamber
All-in-one Machine Weight	85KG
All-in-one Machine Dimensions	40*80*87cm (Length, Width and Height)
Noise	≤63dB
Voltage	AC110V/220V
Rated Power	1500W
Oxygen Production Flow	10L/Min
Oxygen Purity	≥95%±3%
Chamber Capsule Color	White/Grey (Accept Customization)
Total Chamber Weight	13KG
Chamber Size (Length, Width and Height)	210*80*80cm
Disposition	Chamber,All-in-one Machine,Dual Refrigeration System,Sleeping Pad,Bracket, Medical Grade Filter,Negative Ion Disinfection System,Air Dehumidification

APPLICATION SCENARIO



MEDICAL INSTITUTION



OXYGEN BAR CLUB



ELDERLY CARE CENTER



FITNESS INDUSTRY



BEAUTY INDUSTRY



HOTEL



PLATEAU AREA



INDIVIDUAL FAMILY



WORK PLACE

SHIPPING AND DELIVERY



FAQ

1. How does hyperbarics work?

Hyperbarics is a technology in which the air pressure in the environment is increased. When a person's body is placed in a more pressure environment, it absorbs more oxygen molecules per volume of compressed air.

The body normally transports oxygen via the hemoglobin of the red blood cells. By increasing the air pressure, oxygen is then driven into the body's fluids, allowing a super-saturation of tissues and organs with oxygen. The increased pressure infuses the body with oxygen, even reaching injuries with damaged circulation.

2. What are the health benefits of hyperbaric oxygen chambers?

Studies have shown that hyperbaric oxygen chambers can treat hypoxic diseases, relieve fatigue, delay aging, prevent cardiovascular disease, improve sleep quality, accelerated wound repair and enhance immunity. We have customers serviced shows hyperbaric oxygen therapy makes very positive results.

3. How long is the treatment?

Usually, we recommend our clients receive the HBOT therapy 5-6 sessions per week, 60-90mins per session. Using frequency or duration can be adjusted according to individual circumstances.

4. Is HBOT safe?

Hyperbaric oxygen therapy is prescribed by a physician and performed under medical supervision. Although there are minor risks like all medical treatments, overall hyperbaric oxygen therapy is extremely safe.

5. Is there anyone who can not use our chamber?

Pregnant ladies, patients who catch a cold or Flu, or those who have rhinitis or otitis media are not recommended to use the chamber.

6. How do you feel inside a hyperbaric oxygen chamber?

When lying in the oxygen chamber, you will feel comfortable and at ease. Also, You may feel a little squeezing of the ears because of the pressure. It is quite normal and can be relieved by chewing or swallowing.

7. After bought the chamber, do I need to by any other things?

No, the price of our chamber includes the chamber capsule, all machines, namely the air compressor, oxygen concentrator, air cooler, and accessories like frame, tubes, etc, so you don't have to buy anything else.

8. Is it complex to install the chamber or change the filter?

Of course not, we will provide you with a video to show you how to install the chamber and change the filter.

We also support online services to help you with all the issues, so don't worry.

9. What is the general time warranty expected for ZHIKANG oxygen chambers?

We will provide one-year warranty and lifelong after-sales services.

If any quality problem/fault in material/design is under correct operation within one year,

If it's easy to fix, we will send new components for free and teach you how to repair them.

(Even after the one-year warranty period, you can also contact us freely and we will be glad to help you.)



Precision Zipper

Using the YKK brand in the zipper industry, the cabin is hand sewn and heat sealed.

TPU Material

Nanocomposite fibers can maintain good stability when used for a long time or in high-intensity environments, and have a long service life, which can last more than 5 years.



Multi Directional Visual Window

There are observation windows all around to observe the situation of personnel inside the cabin at any time. The cabin is full of light. It can eliminate the tension of the cabin personnel.

Oil Free Compressor

Oil-free, low noise, upgraded heat dissipation (air duct + efficient fan). Stable performance, strong power, rapid oxygen production.




**SOFT TYPE
HYPERBARIC OXYGEN CHAMBER
ACCESSORIES**



2 Automatic Exhaust Valves 2 with Base



1 Manual Exhaust Valve 1 with Base



3 Air Inlets




1 Plug (spare port)



Internal Pressure Gauges



External Pressure Gauges



1 Sponge Silencer



Air Tube and Oxygen Tube



Oxygen Adaptor

SOFT TYPE CHAMBER PRODUCT DISPLAY



EMS04



Lying Soft Type
Hyperbaric Oxygen Chamber

EMS03



Sitting Soft Type
Hyperbaric Oxygen Chamber

EMS01



Sitting Soft Type
Hyperbaric Oxygen Chamber

EMS06



Spherical Soft Type
Hyperbaric Oxygen Chamber

OUR CERTIFICATE





Shenzhen Emeng Health Technology Co., LTD(O2EMENG) located in Shenzhen, is a professional manufacturer specializing in the independent development, production and sales of Hyperbaric Oxygen Chambers.



We are confident in making customized designs to meet our customers' ODM and OEM requirements.



O2EMENG specialize in all kinds of Hyperbaric Oxygen Chambers ranging from soft and hard lying type, soft sitting type, to veterinary hard type chambers, which are widely used in homecare, clinics, health centers, spas, sports industry and high altitude areas, etc..

OUR FACTORY

QC DEPARTMENT

Including incoming inspection, semi-finished products inspection, finished product inspection, and pre-delivery inspection, your satisfaction will be O2 EMENG honor



O2EMENG

After Sale Service

WhatsApp: +86 13760822077

E-mail: kaylayi0612@foxmail.com

Hotline: +86 13760822077

Official website: emoxygen.com





Shenzhen Emeng Health Technology Co., LTD



+86 13760822077



kaylayi0612@foxmail.com



oxygenconcentratormedical.com

Room D, 7th Floor, Guang Long Building, 162 Ping Xin North Road, PingHu Street, LongGang District,
ShenZhen