

Hyperbaric Pressure Chamber For Medical Treatment 1.5 - 2.0ATA

Basic Information

- Place of Origin: SHENZHEN EMENG HEALTH
 Brand Name: O2EMENG
 Certification: CE/FDA/FC/ROHS/UKCA/SGS
 Model Number: EMY01
 Minimum Order 1PC Quantity:
- Price: \$24000
 Packaging Details: wooden packing
 Delivery Time: 15-20 work days
 Payment Terms: D/P, T/T, Western Union, MoneyGram

100PCS

Supply Ability:



Product Specification

- Oxygen Concentration: ≥90%(can Be Set)
- Working Pressure: 1.5-2.0ATA adjustable
- Cabin Szie: 185*88*152CM
- Model:

• Highlight:

- Cabin Weight: 1000KGS
- Usage: Medical/Therapeutic
- Noise Level: Less Than 63 DB
- Properties: Rehabilitation Therapy Supplies

EMY01

medical hyperbaric pressure chamber, medical high pressure hyperbaric chamber, therapeutic hyperbaric pressure chamber



More Images







Product Description

O2EMENG Pressure Hyperbaric pressure chamber for Medical Treatment 1.5-2.0ATA

Product Description:

The Hyperbaric Oxygen Chamber is built with the highest quality materials, ensuring durability and reliability. With a cabin size of 185*88*152CM, this Hyperbaric pressure chamber can accommodate people of different sizes, making it an excellent choice for individual and group use.

The model of this Hyperbaric chamber is EMY01, and it is classified under Rehabilitation Therapy Supplies. It is a must-have for individuals who want to enhance their physical health and mental wellbeing.

The Hyperbaric Oxygen Chamber provides an oxygen concentration of \geq 90%, which can be easily set according to the user's requirements. This feature makes it suitable for individuals with respiratory issues, as it helps to improve oxygen intake and reduce breathing difficulties.

The working pressure of this Hyperbaric Oxygen Chamber is adjustable between 1.5-2.0ATA. This feature allows users to customize the pressure according to their unique needs, making it an ideal option for individuals who require different levels of pressure for their treatment.

Hyperbaric oxygen treatment is a popular therapy used by athletes, patients with chronic illnesses, and individuals seeking to improve their overall health. This Hyperbaric chamber of oxygen provides a safe and comfortable environment for individuals to undergo this treatment.

In conclusion, the Hyperbaric Oxygen Chamber is an excellent investment for individuals who want to enhance their physical health and mental wellbeing. It is built with quality materials, adjustable working pressure, and a comfortable cabin size. This Hyperbaric chamber of oxygen provides a safe and effective way to undergo Hyperbaric oxygen treatment, making it a must-have for anyone looking to improve their health.

Features:

Product Name: Hyperbaric Oxygen Chamber Model: EMY01 Cabin Size: 185*88*152CM Oxygen Concentration: ≥90%(can Be Set) Cabin Weight: 1000KGS Usage: Medical/Therapeutic Features: Oxygen therapy chamber Hyperbaric oxygen therapy unit HBOT chamber

Technical Parameters:

Usage	Medical/Therapeutic
Properties	Rehabilitation Therapy Supplies
Cabin Size	185*88*152CM
Model	EMY01
Function	Hyperbaric Oxygen Therapy
Oxygen Concentration	≥90% (can be set)
Working Pressure	1.5-2.0ATA (adjustable)
Cabin Weight	1000KGS
Noise Level	Less Than 63 DB

Applications:

The O2EMENG Hyperbaric Oxygen Chamber is a durable and reliable product that is perfect for use in hospitals, clinics, and other medical facilities. It is a safe and effective way to provide hyperbaric oxygen therapy to individuals who need it. This product is designed to be easy to use and is suitable for use by medical professionals.

The O2EMENG Hyperbaric Oxygen Chamber has a cabin weight of 1000KGS and a cabin size of 185*88*152CM. The noise level of this product is less than 63 DB which makes it a quiet machine that can be used in any environment without causing disturbance. This product is ideal for use in a wide range of scenarios including: Rehabilitation centers

Sports medicine clinics Wound care centers Hospitals Private clinics Home use The O2EMENG Hyperbaric Oxygen Chamber is a versatile product that can be used to help treat a range of medical conditions. Hyperbaric oxygen therapy units are used to treat a variety of conditions including: Diabetic wounds Carbon monoxide poisoning Decompression sickness

Gas gangrene Burns

Necrotizing infections

The O2EMENG Hyperbaric Oxygen Chamber is a reliable and effective product that can help individuals recover faster from certain medical conditions. If you are looking for a high-quality hyperbaric oxygen therapy unit, then the O2EMENG Hyperbaric Oxygen Chamber model number EMY01 is definitely worth considering.

Support and Services:

Our Hyperbaric Oxygen Chamber product is designed to provide therapeutic benefits to patients by increasing the amount of oxygen delivered to their tissues. Our technical support team is available 24/7 to assist with any issues related to the operation or maintenance of the chamber. We also offer a range of services to ensure that our customers can get the most out of their investment, including installation, training, and ongoing support.

Advantage:

HIGH QUALITY MATERIALS & DESIGN: The O2EMENG is manufactured from 304 stainless steel and high strength, highly permeable explosion-proof materials. The door offers high sealing and stability, providing a safe and reliable environment for use. These features not only enhance the durability and safety of the product, but also ensure pressure stabilization and fast response in high-pressure environments.

Intelligent Design and Voice Control System: O2EMENG is equipped with an intelligent voice control system that provides a convenient user experience. In addition, a USB-A charging port is provided to facilitate charging and connecting other devices in the high-pressure chamber.

Data Monitoring and Regulation: The O2EMENG is equipped with a data monitoring function that monitors the room temperature, humidity, pressure and oxygen concentration in real time to ensure the user's safety and comfort in a high voltage environment.

Low Noise and Energy Efficiency: The O2EMENG operates at a sound level of 43 decibels or more, ensuring a quiet and comfortable indoor environment. In addition, the product's power rating of 1500 watts provides high energy efficiency for long periods of use.

These advantages make the O2EMENG a full-featured, safe, reliable, comfortable and quiet high-pressure room product suitable for a variety of application scenarios, from medical to scientific research.

OUR CERTIFICATE

...



SHIPPING AND DELIVERY





How does hyperbarics work?

Hyperbarics is a technology in which the air pressure in the environment is increased. When a person's body is placed in a more pressure environment, it absorbs more oxygen molecules per volume of compressed air.

The body normally transports oxygen via the hemoglobin of the red blood cells. By increasing the air pressure, oxygen is then driven into the body's fluids, allowing a super-saturation of tissues and organs with oxygen. The increased pressure infuses the body with oxygen, even reaching injuries with damaged circulation.

2. What are the health benefits of hyperbaric oxygen chambers?

Studies have shown that hyperbaric oxygen chambers can treat hypoxic dieases, relieve fatigue, delay aging, prevent cardiovascular disease, improve sleep quality, accelerated wound repair and enhance immunity. We have customers serviced shows hyperbaric oxygen therapy makes very positive results.

3. How long is the treatment?

Usually, we recommend our clients receive the HBOT therapy 5-6 sessions per week, 60-90mins per session. Using frequency or duration can be adjusted according to individual circumstances.

4. Is HBOT safe?

Hyperbaric oxygen therapy is prescribed by a physician and performed under medical supervision. Although there are minor risks like all medical treatments, overall hyperbaric oxygen therapy is extremely safe.

5. Is there anyone who can not use our chamber?

Pregnant ladies, patients who catch a cold or Flu, or those who have rhinitis or otitis media are not recommended to use the chamber.

6. How do you feel inside a hyperbaric oxygen chamber?

When lying in the oxygen chamber, you will feel comfortable and at ease. Also, You may feel a little squeezing of the ears because of the pressure. It is quite normal and can be relieved by chewing or swallowing.

7 After bought the chamber, do I need to by any other things?

No, the price of our chamber includes the chamber capsule, all machines, namely the air compressor, oxygen concentrator, air cooler, and accessories like frame, tubes, etc, so you don't have to buy anything else.

8. Is it complex to install the chamber or change the filter?

Of course not, we will provide you with a video to show you how to install the chamber and change the filter.

We also support online services to help you with all the issues, so don't worry.

9. What is the general time warranty expected for ZHIKANG oxygen chambers?

We will provide one-year warranty and lifelong after-sales services. If any quality problem/fault in material/design is under correct operation within one year, If it's easy to fix, we will send new components for free and teach you how to repair them. (Even after the one-year warranty period, you can also contact us freely and we will be glad to help you.)

CHAMBER CAPSULE DESCRIPTION

...

1. Automatic digital display real-time temperature, humidity, oxygen purity and pressure.

- 2 Rotary timer system, use of time can be pressed in advance Multiple noise reduction technology.
- 3. Excellent Silence performance, the noise is <63dB. Offer you a peaceful experience environment.

4. Double pressure relief device Spiral step less super smooth rotary valve. Spiral ergonomic handheld design. Pressure relief is smoother and more convenient.Internal and external two-way integrated polycarbonate valve.One valve, two-way operation.

5. Automatic pressurize relief valves

Two automatic high precision constant pressure valve. Ensure the product is safe and effective by Aluminum alloy steel, high hardness, wear resistance.

6. Communicating system.

Intercom telephone is convenient for users to talk with operator outside.

7 Machine start-up protection system.

The chamber only can be turned on when the chamber door is completely closed, avoiding the danger of the door being left unlocked.

8. Automatic pressurize and depressurize system.

Based on the principle of "Cylinder (engine)", Pressurize automatically after the door is completely closed.

9. Made by Aluminum alloy steel; With environmental spray paint;

Strong weight bearing capacity, Durable and wear-resistant; non-deformable; long lifespan.

10. In case of unexpected power failure, the door can be opened by emergency safety from inside the chamber.



WHY ARE MORE AND MORE PEOPLE CHOOSING **HYPERBARIC OXYGEN CHAMBERS?**

0

n some parts of the body, the capillaries are so narrow that blood cells cannot pass through mem, bound oxygen cannot reach them, and they can only rely on dissolved oxygen. The oxygen generator has little change in the dissolved oxygen in the blood, which is also the asson why the ordinary oxygen absorption effect in the plateau area is general.

G Shenzhen Emeng Health Technology Co., LTD

+86 13760822077

kaylayi0612@foxmail.com

02EMENG

© oxygenconcentratormedical.com

Room D, 7th Floor, Guang Long Building, 162 Ping Xin North Road, PingHu Street, LongGang District, ShenZhen