



Shenzhen Emeng Health Technology Co., LTD

oxygenconcentratormedical.com

EMS04 Noise Level Less Than 63 DB Hyperbaric Oxygen Therapy Chamber for Rehabilitation

Our Product Introduction

For more products please visit us on oxygenconcentratormedical.com

Basic Information

- Place of Origin: SHENZHEN EMENG HEALTH
- Brand Name: O2EMENG
- Certification: CE/FDA/FC/ROHS/UKCA/SGS
- Model Number: EMS04
- Minimum Order Quantity: 1PC
- Price: \$3750
- Packaging Details: wooden packing
- Delivery Time: 15-20 work days
- Payment Terms: D/P, T/T, Western Union, MoneyGram
- Supply Ability: 100PCS



Product Specification

- Noise Level: Less Than 63 DB
- Total Weight: 83kgs
- Properties: Rehabilitation Therapy Supplies
- Model: EMS04
- Oxygen Concentration: $\geq 90\%$ (can Be Set)
- Function: Hyperbaric Oxygen Therapy
- Working Pressure: 1.3-1.5ATA (can Be Set)
- Machine CTN Size: 127*55*105cm
- Highlight: 1.5ata hyperbaric oxygen chamber,
1.5ata hyperbaric oxygen therapy chamber



More Images



Product Description

EMS04 Noise Level Less Than 63 DB Hyperbaric Oxygen Therapy Chamber for Rehabilitation

Product Description:

One of the key features of the EMS04 model is its low noise level. This chamber operates at less than 63 dB, which means that you can use it comfortably without any disturbance. Whether you're using it at home or in a professional setting, this pressurized oxygen chamber won't cause any noise-related issues.

Of course, the main function of this chamber is to provide hyperbaric oxygen therapy. This type of therapy involves breathing in pure oxygen while inside a pressurized chamber. The high pressure forces more oxygen into your body's tissues, which can help with a variety of health conditions.

The EMS04 model is designed to deliver hyperbaric oxygen therapy in a safe and effective way. It is made from high-quality materials and is built to withstand the pressures involved in this type of therapy. Whether you're using it for medical or wellness purposes, you can trust that this pressurized oxygen chamber will deliver the results you need.

Despite its sturdy construction, the EMS04 model is also relatively lightweight. At 83 kgs, it is easy to move around and set up wherever you need it. This makes it an excellent choice for anyone who needs a portable hyperbaric chamber of oxygen.

Finally, the working pressure of this chamber can be set between 1.3-1.5ATA. This means that you can customize your hyperbaric oxygen therapy to suit your individual needs. Whether you need a higher or lower pressure, this pressurized oxygen chamber can be adjusted accordingly.

Technical Parameters:

Product Attribute	Value
Oxygen Concentration	≥90% (can be set)
Noise Level	Less than 63 DB
Function	Hyperbaric Oxygen Therapy
Total Weight	83kgs
Machine CTN size	127*55*105cm
Model	EMS04
Working pressure	1.3-1.5ATA (can be set)
Properties	Rehabilitation Therapy Supplies
Usage	Medical/Therapeutic

Applications:

Manufactured in SHENZHEN EMENG HEALTH, this product is made with the finest materials to ensure durability and longevity. The O2EMENG EMS04 has a working pressure of 1.3-1.5ATA, which can be set to your desired level. It has a total weight of 83kgs, making it easy to move around and position in your facility.

One of the standout features of this product is its noise level, which is less than 63 DB. This means that it operates quietly, ensuring that patients can relax and recover without any disturbance. Additionally, the oxygen concentration of the O2EMENG EMS04 is ≥90%, which can also be set to your desired level.

The O2EMENG EMS04 Hyperbaric Oxygen Chamber is ideal for a wide range of medical and therapeutic applications. It is perfect for patients with injuries, chronic pain, and other health conditions that require oxygen therapy. Additionally, it can be used to speed up the healing process after surgery or other medical procedures.

Some of the common scenarios where this product can be used include post-stroke rehabilitation, wound healing, sports injury recovery, and chronic fatigue syndrome treatment. It can also be used to treat conditions such as autism, cerebral palsy, and multiple sclerosis.

Overall, the O2EMENG EMS04 Hyperbaric Oxygen Chamber is a high-pressure oxygen chamber that is designed to provide top-quality oxygen therapy to patients. It is a reliable and effective product that is perfect for medical and therapeutic applications.

Advantage:

Highly Adjustable Oxygen Concentration: The O2EMENG can be adjusted to ensure that the oxygen concentration reaches or exceeds 90%, adapting to the oxygen needs of different users and ensuring adequate oxygen supply, especially in high-pressure environments.

Rapid Decompression Function: The design allows O2EMENG to quickly reduce pressure within 5-10 minutes, effectively reducing the discomfort caused by high-pressure environments, and improving the user's adaptability and comfort.

Intelligent Dual Screen Control System: The O2EMENG is equipped with a dual screen intelligent control system, which provides a more convenient operating interface and function selection, allowing users to easily control and monitor various parameters during use.

High-quality materials and design: Imported from Japan, the O2EMENG is made of YKK medical-grade zippers and medical-grade TPU materials to ensure the durability and safety of the product. It is also equipped with metal bracket and customizable shell color to enhance the practicality and aesthetics of the product.

These advantages make the O2EMENG a safe and reliable choice in high-pressure environments, suitable for users who enter high-pressure environments for long periods of time or frequently, providing a comfortable and convenient experience.

Support and Services:

The Hyperbaric Oxygen Chamber product comes with comprehensive technical support and services to ensure that you get the most out of your investment. Our team of experienced technicians and engineers are available to provide support and assistance whenever you need it.

Our technical support services include:

1.Demand customization: according to your use environment and expectations, recommend the optimal core products and life aesthetics composition scheme.

2.Exclusive design: the advanced technology team and the top design team jointly, from the size and color material matching, the actual program customization landing.

3. Craftsmanship Manufacturing: The core components and other materials of hyperbaric chamber series products are manufactured in aerospace material factories.

4. Lifetime maintenance: provide perfect professional after-sales service guarantee

At Hyperbaric Oxygen Chamber, we are committed to providing the highest level of support and service to our customers. Contact us today to learn more about our technical support and services.

O2EMENG
Hyperbaric Oxygen Chamber

HYPERBARIC OXYGEN CHAMBER


- IMPROVE BRAIN BLOOD CIRCULATION
- SPEED UP METABOLISM
- HIGH PRESSURE OXYGEN INHALATION



WHY ARE MORE AND MORE PEOPLE CHOOSING HYPERBARIC OXYGEN CHAMBERS?



In some parts of the body, the capillaries are so narrow that blood cells cannot pass through them, bound oxygen cannot reach them, and they can only rely on dissolved oxygen. The oxygen generator has little change in the dissolved oxygen in the blood, which is also the reason why the ordinary oxygen absorption effect in the plateau area is general.



SOFT TYPE HYPERBARIC OXYGEN CHAMBER



Innovative technology,quality assurance



Portable, Small Footprint
Fast pressurization



Made of multi-layer TPU composite
environmentally friendly material



Fresh Air System
Achieve cabin air circulation effect





CONCLUSION



It is more feasible to increase the dissolved oxygen in the blood.
The hyperbaric oxygen chamber can maintain the oxygen partial pressure within a safe and effective range by flexibly adjusting the oxygen concentration or pressure, and gently and effectively improve the oxygen delivery to tissues and organs by increasing dissolved oxygen.




INTEGRATION OF SUPERCHARGING AND OXYGENATION



Overturn the complexity of traditional household



New hyperbaric
oxygen equipment



Oxygen equipment
+
air power system
+
Air Cooler

Intelligent all in one machine, no need to
operate three devices,
easy one-button operation saves time and effort.

FIVE FUNCTIONS OF THE ALL-IN-ONE MACHINE



01

Oxygen Production

02

Pressurization

03

Dehumidification

04

Refrigeration

05

Negative Ion Sterilization



One-button Start

Oxygen Flow Meter

304 Stainless Steel Body

Medical Grade Filter

Medical Grade
Humidification Bottle

360° Universal Wheel

PRODUCT DETAIL



Fresh Air System

The all-in-one machine continuously pressurizes to ensure air circulation in the chamber.



≥95%±3% Medical Oxygen Purity

In the chamber, the oxygen concentration will rise to 93%.



Emergency Pressure Relief Valve

An emergency valve lets you get out of the chamber quickly in an urgent situation.



Two Zippers

Two zippers let people enter the chamber easily and provide better sealing.



Two Side Pressure Gauge

Two side pressure gauges make it easier and safer to monitor the pressure both from inside and outside.



Easy operating by one person without assistant


PRODUCT PARAMETER

Product Model	EMS04 Soft Type Hyperbaric Oxygen Chamber
All-in-one Machine Weight	85KG
All-in-one Machine Dimensions	40*80*87cm (Length, Width and Height)
Noise	≤63dB
Voltage	AC110V/220V
Rated Power	1500W
Oxygen Production Flow	10L/Min
Oxygen Purity	≥95%±3%
Chamber Capsule Color	White/Grey (Accept Customization)
Total Chamber Weight	13KG
Chamber Size (Length, Width and Height)	210*80*80cm
Disposition	Chamber,All-in-one Machine,Dual Refrigeration System,Sleeping Pad,Bracket, Medical Grade Filter,Negative Ion Disinfection System,Air Dehumidification

APPLICATION SCENARIO



MEDICAL INSTITUTION



OXYGEN BAR CLUB



ELDERLY CARE CENTER



FITNESS INDUSTRY



BEAUTY INDUSTRY



HOTEL



PLATEAU AREA



INDIVIDUAL FAMILY



WORK PLACE

SHIPPING AND DELIVERY



FAQ

1. How does hyperbarics work?

Hyperbarics is a technology in which the air pressure in the environment is increased. When a person's body is placed in a more pressure environment, it absorbs more oxygen molecules per volume of compressed air.

The body normally transports oxygen via the hemoglobin of the red blood cells. By increasing the air pressure, oxygen is then driven into the body's fluids, allowing a super-saturation of tissues and organs with oxygen. The increased pressure infuses the body with oxygen, even reaching injuries with damaged circulation.

2. What are the health benefits of hyperbaric oxygen chambers?

Studies have shown that hyperbaric oxygen chambers can treat hypoxic diseases, relieve fatigue, delay aging, prevent cardiovascular disease, improve sleep quality, accelerated wound repair and enhance immunity. We have customers serviced shows hyperbaric oxygen therapy makes very positive results.

3. How long is the treatment?

Usually, we recommend our clients receive the HBOT therapy 5-6 sessions per week, 60-90mins per session. Using frequency or duration can be adjusted according to individual circumstances.

4. Is HBOT safe?

Hyperbaric oxygen therapy is prescribed by a physician and performed under medical supervision. Although there are minor risks like all medical treatments, overall hyperbaric oxygen therapy is extremely safe.

5. Is there anyone who can not use our chamber?

Pregnant ladies, patients who catch a cold or Flu, or those who have rhinitis or otitis media are not recommended to use the chamber.

6. How do you feel inside a hyperbaric oxygen chamber?

When lying in the oxygen chamber, you will feel comfortable and at ease. Also, You may feel a little squeezing of the ears because of the pressure. It is quite normal and can be relieved by chewing or swallowing.

7. After bought the chamber, do I need to by any other things?

No, the price of our chamber includes the chamber capsule, all machines, namely the air compressor, oxygen concentrator, air cooler, and accessories like frame, tubes, etc, so you don't have to buy anything else.

8. Is it complex to install the chamber or change the filter?

Of course not, we will provide you with a video to show you how to install the chamber and change the filter.

We also support online services to help you with all the issues, so don't worry.

9. What is the general time warranty expected for ZHIKANG oxygen chambers?

We will provide one-year warranty and lifelong after-sales services.

If any quality problem/fault in material/design is under correct operation within one year,

If it's easy to fix, we will send new components for free and teach you how to repair them.

(Even after the one-year warranty period, you can also contact us freely and we will be glad to help you.)



Precision Zipper

Using the YKK brand in the zipper industry, the cabin is hand sewn and heat sealed.

TPU Material

Nanocomposite fibers can maintain good stability when used for a long time or in high-intensity environments, and have a long service life, which can last more than 5 years.



Multi Directional Visual Window

There are observation windows all around to observe the situation of personnel inside the cabin at any time. The cabin is full of light. It can eliminate the tension of the cabin personnel.

Oil Free Compressor

Oil-free, low noise, upgraded heat dissipation (air duct + efficient fan). Stable performance, strong power, rapid oxygen production.




**SOFT TYPE
HYPERBARIC OXYGEN CHAMBER
ACCESSORIES**



2 Automatic Exhaust Valves 2 with Base



1 Manual Exhaust Valve 1 with Base



3 Air Inlets




1 Plug (spare port)



Internal Pressure Gauges



External Pressure Gauges



1 Sponge Silencer



Air Tube and Oxygen Tube



Oxygen Adaptor

SOFT TYPE CHAMBER PRODUCT DISPLAY



EMS04



Lying Soft Type
Hyperbaric Oxygen Chamber

EMS03



Sitting Soft Type
Hyperbaric Oxygen Chamber

EMS01



Sitting Soft Type
Hyperbaric Oxygen Chamber

EMS06



Spherical Soft Type
Hyperbaric Oxygen Chamber

OUR CERTIFICATE



Shenzhen Emeng Health Technology Co., LTD

+86 13760822077

kaylayi0612@foxmail.com

oxygenconcentratormedical.com

Room D, 7th Floor, Guang Long Building, 162 Ping Xin North Road, PingHu Street, LongGang District,

