

DVOT (8 IN ONE) Antigen Combo Test Kit for testing SARS-CoV-2 & RSV & MP & ADV & Flu A/B & PIV1/3&PIV

Our Product Introduction

for more products please visit us on oxygenconcentratormedical.com

Basic Information

- Place of Origin: china
- Brand Name: DVOT
- Certification: CE/ISO
- Model Number: 8 in one
- Minimum Order Quantity: 720Test
- Price: 1.8 ~ 2 USD Per test
- Packaging Details: 720T/PER CTN Weight :25KGS CTN SIZE :53*50*43CM
- Delivery Time: 15-20 working days
- Payment Terms: D/P, T/T, Western Union, MoneyGram
- Supply Ability: 1000000pcs



Product Specification

- Product Name: Multiplex Respiratory Antigen Test Kit
- Time Of Use: One Time Off Use Only
- Reaction Time: 10-15 Mins
- Shelf Life: 24 Months



More Images



Product Description

DVOT (8 IN ONE) Antigen Combo Test Kit for testing SARS-CoV-2 & RSV & MP & ADV & Flu A/B & PIV1/3&PIV2

Advantages:

During the first few days of a Covid-19 illness, symptoms tend to be nonspecific and may include cough, fever and sore throat. These symptoms can be indistinguishable from other viral infections, such as influenza, respiratory syncytial virus (RSV), or viruses associated with the common cold. Because of this, combo diagnostic testing is needed and important, not only to determine the cause of the illness, but also to determine whether treatment is appropriate. DVOT® SARS-CoV-2 & RSV & MP & ADV & Flu A/B/PIV 1/3 & PIV 2 Antigen Combo Test Kit Antigen Combo Test Kit is an easy and rapid way to identify which infections cause the symptoms. It can be used anywhere and anytime by yourself.





Key benefits:

1. Reliable test performance
2. Fast answers in just 15 minutes
3. Supports differential detection of SARS-CoV-2
4. Simple, 4-step process
5. Fully integrated assay, requiring only one application area and one sample
6. Optimized combinations to identify antigens that cause symptoms of infection
7. Can be used anywhere and anytime by yourself

Steps to use:

Step 1: Insert a sterile swab into the patient's nostril. Rotate the swab 5 times over the nasopharyngeal surface and then leave the swab in place for 15 seconds to absorb secretions. Remove the swab from the nasal cavity.

Step 2: Immerse the swab in the sample container, making sure the swab is in full contact with the fluid, but do not apply too much force to avoid splashing. Cover and gently shake the sample container.

Step 3: Place the test instrument on a flat surface and place 4 drops of the immersion sample into the sample hole of the test instrument at a 90° angle.

Step 4: Read test results within 15 to 30 minutes. Do not touch or move the test apparatus until results can be read.

QUICK USER GUIDE



😊 DO USE

- 😊 As an aid in the diagnosis of current COVID-19 / Influenza A/B / ADV/MP&RSV/PIV1&3 / PIV2 test.
- 😊 If you have COVID-19 / Influenza A/B / ADV / MP&RSV/ PIV1&3 / PIV2 specific symptoms within the first 7 days of symptom onset.

😞 DO NOT USE WHEN:

- 😞 Baby under 2 years old.
- 😞 If you are prone to nose bleeds.
- 😞 If you have had a facial or head injury/surgery in the last 6 months.

03/07/2024
LOT: 4A24070301
02/07/2026



6 973083 141189

Infected or not? Test right now!



New
Omicron
EG.5 .CH.1.1-.JN.1
Detected

Multiplex Respiratory Antigen Test Kit

SARS-CoV-2 & RSV & MP & ADV & Flu A/B & PIV1/3 & PIV2

8 In 1



Anterior **Nasal** Swab
Anteriorer Nasenabstrich

SARS-CoV-2 **94.55%**
Sensitivity

Influenza A **90.48%**
Sensitivity

Influenza B **91.67%**
Sensitivity

RSV **93.62%**
Sensitivity



MP **90.82%**
Sensitivity

ADV **91.18%**
Sensitivity

PIV1/3 **92.23%**
Sensitivity

PIV2 **90.20%**
Sensitivity



About the Antigen Combo Test Kit

During the first few days of a Covid-19 illness, symptoms tend to be nonspecific and may include cough, fever and sore throat. These symptoms can be indistinguishable from other viral infections, such as influenza, respiratory syncytial virus (RSV), or viruses associated with the common cold. Because of this, combo diagnostic testing is needed and important, not only to determine the cause of the illness, but also to determine whether treatment is appropriate.

DVOT® SARS-CoV-2 & RSV & MP & ADV & Flu A/B & PIV1/3 & PIV2 Antigen Combo Test Kit is an easy and rapid way to identify which infections cause the symptoms. It can be used anywhere and anytime by yourself.



Get result at home



Easy self-testing



Time-saving

Anytime, Anywhere, Test by yourself !



Hospital



Test Site



Airport



Station



Hotel



Corporation



Mass Screening

HAVE YOU EVER HAD SUCH CONCERNS?

Worried about inconvenient to test

- No hospitals or test site available nearby.
- Physical discomfort, queue for a long time to test.

Worried about can not get result immediately

- Long time waiting for a PCR test result.
- Want to know the result immediately.

Worried about the family

- Anyone in my family been infected?
- How to test elder and infants?

Worried about infection

- Infection risk in crowded hospitals or test site.
- If there are close contacts around you?



Let us
SOLVE
your troubles!


It can detect mutated virus

High accuracy detection $\geq 96\%$

- ✓ Test result is available in 8 minutes.
- ✓ Self-testing is convenient at home, or anywhere else.
- ✓ Risk control for easy life and working.
- ✓ Shelf-life of 1.5 years, you can store for a long time.


Easy detection will make you feel ease, people around you also feel comfortable.

Simple steps, easy-handling for anybody!





Result Reading

For Cov & Rsv



For MP & ADV





Shipping Details



CB	
QTY :	360 PCS
G.W :	14 KG
MEAS :	46X44X44CM



Shenzhen Emeng Health Technology Co., LTD

+86 13760822077

kaylayi0612@foxmail.com

oxygenconcentratormedical.com

Room D, 7th Floor, Guang Long Building, 162 Ping Xin North Road, PingHu Street, LongGang District, ShenZhen