

**Shenzhen Emeng Health Technology Co., LTD** oxygenconcentratormedical.com

# **DVOT 8 IN ONE Rapid Antigen Test Kit 98% Accuracy CE and ISO Certified** that use for testing SARS-CoV-2 & RSV & MP etc.

### **Basic Information**

• Place of Origin: china DVOT • Brand Name: · Certification: CE/ISO • Model Number: 8 in one • Minimum Order 720 test

Quantity:

• Price: 1.8 ~ 2 USD Per test

 Packaging Details: 720T/PER CTN Weight :25KGS CTN SIZE

:53\*50\*43CM

• Delivery Time: 15-20 work days

• Payment Terms: D/P, T/T, Western Union, MoneyGram

• Supply Ability: 1000000pcs



# **Product Specification**

• Test Time: 15 Minutes

• Result Time: 15 Minutes To Get Rapid Test Kit • Type: · Keyword: Antigen Rapid Test

• Time Of Use: Disposable · Accuracy: 98% CE ISO · Certifications:

• Highlight: ce clungene antigen rapid test, ce medical diagnostic test kits,

disposable clungene antigen rapid test



# More Images









### **Product Description**

# DVOT (8 IN ONE)Rapid Antigen Test Kit 98% Accuracy CE and ISO Certified that use for testing SARS-CoV-2 & RSV & MP etc.

#### **Product Description:**

Our Antigen Rapid Test Kit is designed to detect specific proteins on the surface of the virus. This test is ideal for individuals who are experiencing symptoms of COVID-19. With a testing time of only 15 minutes, this test kit is perfect for organizations and businesses that require fast and accurate results. The Antigen Rapid Test Kit is a straightforward and easy-to-use diagnostic tool that requires no specialized equipment or training.

In addition to our Antigen Rapid Test Kit, we also offer an Antibody Rapid Test Kit. This test is designed to detect antibodies that are produced by the immune system in response to a COVID-19 infection. The Antibody Rapid Test Kit is ideal for individuals who may have been exposed to COVID-19 but are not currently experiencing any symptoms. This test is also quick and easy to use, with results available in only 15 minutes.

Our Medical Diagnostic Test Kits are designed to meet the needs of businesses, organizations, and individuals who require fast, accurate, and reliable COVID-19 testing. Our Antigen Rapid Test Kit and Antibody Rapid Test Kit are both simple and easy to use, with no specialized equipment or training required. With the certification of CE and ISO, you can trust the quality and accuracy of our products.

#### **Technical Parameters:**

Infectious Disease Test Kit Type:	Rapid Test Kit
Diagnostic Test Kits Certification:	CE, ISO
Antigen Home Test Kit Keyword:	Antigen Rapid Test
Time of Use:	Disposable
Test Time:	15 Minutes
Accuracy:	98%
Result Time:	15 Minutes to Get

#### **Applications:**

These Antibody Rapid Test Kits are easy to use and provide results in just 15 minutes, making them ideal for use in a variety of settings. They are manufactured in China, ensuring high quality and reliability.

DOVT 8 in one Medical Diagnostic Test Kits are suitable for use in a range of occasions and scenarios, including: **Healthcare facilities:** These kits are an essential tool for healthcare professionals working in hospitals, clinics, and other medical settings. They can be used to quickly identify and isolate infected patients, preventing the spread of disease.

**Travel:** Many countries require travelers to provide proof of a negative COVID-19 test before entering. These kits are a convenient and reliable way to get tested before traveling.

**Schools and universities:** With the ongoing pandemic, it's important for schools and universities to have a way to quickly test students and staff for COVID-19. These kits can be used for regular testing to help prevent outbreaks. **Workplaces:** Many employers are implementing testing programs to help keep their employees safe. These kits are a great option for on-site testing.

DOVT 8 in one Medical Diagnostic Test Kits have an accuracy rate of 98%, and are certified by CE and ISO. They are an essential tool for anyone looking to quickly and accurately test for infectious diseases.

#### **Support and Services:**

Our Medical Diagnostic Test Kits come with comprehensive technical support and services to ensure optimal performance and accuracy. We are committed to ensuring that our customers have the necessary support and resources to effectively utilize our Medical Diagnostic Test Kits.

#### Steps to use:


Step 1: Insert a sterile swab into the patient's nostril. Rotate the swab 5 times over the nasopharyngeal surface and then leave the swab in place for 15 seconds to absorb secretions. Remove the swab from the nasal cavity.

**Step 2:** Immerse the swab in the sample container, making sure the swab is in full contact with the fluid, but do not apply too much force to avoid splashing. Cover and gently shake the sample container.

Step 3: Place the test instrument on a flat surface and place 4 drops of the immersion sample into the sample hole of the test instrument at a 90° angle.


Step 4: Read test results within 15 to 30 minutes. Do not touch or move the test apparatus until results can be read.







СВ	
QTY:	360 PCS
G.W:	14 KG
MEAS:	46X44X44CM





+86 13760822077

kaylayi0612@foxmail.com

Oxygenconcentratormedical.com

Room D, 7th Floor, Guang Long Building, 162 Ping Xin North Road, PingHu Street, LongGang District, ShenZhen