

ISO Certified Disposable Medical Diagnostic Test Kits for Rapid Testing Malaria Pf/Pan

Our Product Introduction

For products please visit us on oxygenconcentratormedical.com

Basic Information

- Place of Origin: china
- Brand Name: DVOT
- Certification: CE/ISO
- Model Number: Malaria Pf/Pan Rapid Test Kit
- Minimum Order Quantity: 360 test
- Price: 1-1.2USD
- Packaging Details: 360T/PER CTN Weight :11KGS CTN SIZE :46*44*44CM
- Delivery Time: 15-20 work days
- Payment Terms: D/P, T/T, Western Union, MoneyGram
- Supply Ability: 1000000pcs





Product Specification

- Test Time: 15 Minutes
- Result Time: 15 Minutes To Get
- Type: Rapid Test Kit
- Time Of Use: Disposable
- Certifications: CE ISO
- Accuracy: 98%
- Keyword: Antigen Rapid Test
- Highlight: disposable clungene antigen rapid test, disposable medical diagnostic test kits, ce clungene antigen rapid test



More Images



Product Description

ISO Certified Disposable Medical Diagnostic Test Kits for Rapid Testing Malaria Pf/Pan

Product Description:

Our Antigen Rapid Test Kits have been certified with both CE and ISO certifications, ensuring that they meet the highest quality standards. Plus, our Antigen Swab Test Kit is disposable, making it easy to use and dispose of safely.

When it comes to accuracy, our Antigen Rapid Test Kits are incredibly reliable, with an accuracy rate of 98%. This makes them ideal for use in a variety of settings, including healthcare facilities, schools, workplaces, and more.

Whether you're looking to test yourself or someone else, our Antigen Swab Test Kit is a convenient and effective way to get fast results. So why wait? Order your Antigen Rapid Test Kits today and start testing with confidence!

Technical Parameters:

Type	Antigen Rapid Self Test Kit
Accuracy	98%
Certifications	CE, ISO
Result Time	15 Minutes
Time Of Use	Disposable
Test Time	15 Minutes
Keyword	Antigen Rapid Test

Applications:

The DOVT Malaria Pf/Pan Rapid Test Kit is ideal for use in various medical settings, including hospitals, clinics, and laboratories. It is also suitable for use in remote and resource-limited areas where access to sophisticated diagnostic equipment is limited. The kit is easy to use and does not require any special training, making it suitable for use by healthcare professionals and non-professionals alike.

The DOVT Malaria Pf/Pan Rapid Test Kit is CE and ISO certified, which ensures that it meets the highest quality standards. The kit is designed to provide accurate and reliable results, with a sensitivity of 98%. This makes it an excellent option for diagnosing Malaria infections.

In addition to medical settings, the DOVT Malaria Pf/Pan Rapid Test Kit is also suitable for home testing. It is an antigen rapid test kit that provides quick and accurate results, making it an ideal option for people who want to check their Malaria status in the comfort of their own homes. The kit is easy to use and comes with detailed instructions, making it suitable for use by non-professionals.

The DOVT Malaria Pf/Pan Rapid Test Kit is suitable for use in various scenarios, including:

Medical settings such as hospitals, clinics, and laboratories

Remote and resource-limited areas

Home testing for individuals who want to check their Malaria status

Overall, the DOVT Malaria Pf/Pan Rapid Test Kit is a reliable and accurate diagnostic tool for detecting Malaria infections. Its ease of use, quick result time, and high accuracy make it an excellent option for use in various medical settings and home testing scenarios.

Support and Services:

Our medical diagnostic test kits are designed to provide accurate and reliable results for a variety of medical conditions. We offer technical support and services to ensure that our customers are able to use our products effectively and efficiently.

Steps to use the product:

Prepare the tools:

- Make sure you have matching capillary pipettes (capillary pipette) and clean blood collection tools (e.g., blood

collection needles, alcohol pads, etc.).

Collect a blood sample:

- Clean the tip of your finger (or other blood collection site) with an alcohol pad and allow the alcohol to evaporate before collecting blood.
- Collect an appropriate amount of blood sample (usually 25 microliters) by gently pricking the tip of the finger with a blood collection needle.
- Aspirate the blood sample using a matched capillary pipette.

Perform the test:

- Add 2 drops of the just drawn blood to the test card
- Then add the buffer solution

Read Results:

- At the end of the waiting time, observe the results area on the test card or container.

How to read test results:

Observe the results area:

- Look for the results interpretation area or color comparison chart according to the product description.
- Observe the results interpretation area on the test card or container for color changes, lines appearing, or other recognizable indications.

Compare Results:


- A "Positive" result indicates the presence of the substance or condition to be tested for in the blood sample and may require consultation with a physician on a case-by-case basis.
- If the result is "Negative", the substance or condition to be tested for is not detected in the blood sample.
- If the result is "Invalid", the test may need to be repeated due to improper handling, sample contamination or invalid test card/reagent.

Caution:

- Ensure that the time requirements in the product instructions are followed when reading the results.
- If there is any doubt or uncertainty about the results, it is advisable to consult a physician or relevant healthcare provider for professional advice.
- Blood collection kits and test cards/containers should be disposed of properly after use to avoid contamination of the environment.







Shenzhen Emeng Health Technology Co., LTD

+86 13760822077

kaylayi0612@foxmail.com

oxygenconcentratormedical.com

Room D, 7th Floor, Guang Long Building, 162 Ping Xin North Road, PingHu Street, LongGang District, ShenZhen