



Disposable Dengue Ns1 Ag & IgG/IgM Rapid Antigen At Home Test Kit with 15-Minute Result

Our Product Introduction

For more products please visit us on oxygenconcentratormedical.com

Basic Information

- Place of Origin: china
- Brand Name: DVOT
- Certification: CE/ISO
- Model Number: Dengue Ns1 Ag & IgG/IgM Combo Test Kit
- Minimum Order Quantity: 360 test
- Price: 0.9-1 USD
- Packaging Details: 360T/PER CTN Weight :11KGS CTN SIZE :46*44*44CM
- Delivery Time: 15-20 work days
- Payment Terms: D/P, T/T, Western Union, MoneyGram
- Supply Ability: 1000000pcs



Product Specification

- Test Time: 15 Minutes
- Accuracy: 98%
- Keyword: Antigen Rapid Test
- Certifications: CE ISO
- Type: Rapid Test Kit
- Result Time: 15 Minutes To Get
- Time Of Use: Disposable
- Highlight: disposable clungene antigen rapid test, disposable medical diagnostic test kits, ce clungene antigen rapid test



More Images



Product Description

Disposable Dengue Ns1 Ag & IgG/IgM Rapid Antigen At Home Test Kit with 15-Minute Result

Product Description:

Dengue fever is an acute insect-borne infectious disease caused by dengue virus, which is mainly transmitted by *Aedes aegypti* and *Aedes albopictus* mosquitoes. The incubation period of dengue fever is usually 1-14 days, mostly 5-9 days. Patients usually have an acute onset of illness, with the first symptom being sudden onset of high fever, which may be accompanied by chills, and a temperature of up to 40°C within 24 hours. It usually lasts for 3~7 days. In some patients, the temperature may drop to normal after 3-5 days of fever, but rise again after 1 day, a phenomenon known as bimodal or saddle fever. If you are not sure whether you are infected with dengue virus, you can use this antigen test.

The Antigen Rapid Test Kit is a disposable product that can be used quickly and easily without the need for special equipment or a laboratory setting. Test results are obtained in as little as 15 minutes, making it a very efficient testing option. This At Home Antigen Test Kit is very convenient for patients who need to be tested at home, such as those who are elderly or have underlying health conditions that make it difficult for them to go to a testing clinic.

The Antigen Rapid Test Kit is an excellent alternative to the Antibody Test Kit. Antibody Test Kits are designed to detect the presence of antibodies in a person's body, indicating that they have been infected with a particular disease. However, Antibody Test Kits can only detect the virus after the body has had time to produce antibodies, which can take several days or even weeks. The Antigen Rapid Test Kit, on the other hand, can detect the virus early on, making it an ideal tool for early detection and containment of infectious diseases.

Overall, the Antigen Rapid Test Kit is a reliable and accurate product that can be used in a variety of settings. Whether you are a medical professional seeking to detect infectious diseases early on or a patient seeking an At Home Antigen Test Kit, the Antigen Rapid Test Kit is an excellent choice.

Technical Parameters:

Product Attribute	Description
Keyword	Covid 19 Test Kits, Influenza Test Kits, Antigen Rapid Testing Kit
Type	Rapid Test Kit
Result Time	15 Minutes To Get
Time Of Use	Disposable
Certifications	CE,ISO
Accuracy	98%
Test Time	15 Minutes

Applications:

The DVOT Dengue Ns1 Ag & IgG/IgM Combo Test Kit can be used in a variety of medical applications. It is particularly useful in diagnosing various diseases, such as dengue fever. Additionally, this test kit is ideal for use in detecting infections caused by bacteria, viruses, and fungi.

With the ongoing Dengue Ns1 Ag & IgG/IgM pandemic, the DOVT Dengue Ns1 Ag & IgG/IgM Combo Test Kit has become a crucial tool in the fight against the virus. It can be used to detect Dengue Ns1 Ag & IgG/IgM infections in patients, providing fast results in just 15 minutes. This Antigen Home Test Kit is a valuable asset for healthcare providers, as well as individuals seeking to monitor their own health.

The DVOT Dengue Ns1 Ag & IgG/IgM Combo Test Kit is designed for one-time use, making it a disposable and cost-effective solution for medical diagnosis. The quick result time and ease of use make it an ideal choice for medical professionals in a variety of settings, including hospitals, clinics, and laboratories.

Support and Services:

Our Medical Diagnostic Test Kits are designed to provide accurate and reliable results for a variety of medical conditions. Our team is available to assist with any questions or concerns you may have regarding our products.

Steps to use:

Preparation Stage:

Ensure that your hands are clean and that you have the required blood collection kit (e.g., blood collection needle, alcohol pads, etc.). Remove the test card from the box and check that the package is intact and that it is within the expiration date. Follow the instructions on the package to confirm the type of sample and sample size required for the test card.

Collect the blood sample:

Clean the blood collection site (usually the tip of the finger) with an alcohol pad and allow the alcohol to evaporate. Use a blood collection needle to gently prick the tip of the finger and collect the appropriate amount of blood sample. Take care to avoid contaminating the sample.

Perform the test:

Step 1: Add a drop of blood sample to the designated area of the test card.
Step 2: Wait a certain amount of time (usually 15 minutes) according to the instructions so that the test card can respond adequately.
Step 3: After completing the wait, observe the results on the test card.

Read the results:

Observe the indicated areas or lines on the test card and judge the results according to the explanations provided in the instructions. Take care to distinguish between the different results and their significance for Dengue Negative, Dengue Positive and Dengue NS1 Test.



Dvot

DENGUE

Dengue NS1Ag & IgG/IgM Combo Test kit

Dengue fever, also known as breakbone fever, is a mosquito-borne tropical disease caused by the dengue virus. Symptoms include fever, headache, muscle and joint pains, and a characteristic skin rash that is similar to measles.

FEATURE

Accuracy, visual result, simple and safe to use

SPECIMEN

Whole Blood/Serum/Plasma

REACTION TIME

15mins

STORAGE TEMPERATURE

4-30°C

SHELF LIFE

24 months




TEST PROCEDURE


• For serum/plasma sample testing



Dengue NS1 Ag




• For whole blood sample testing



15 Minutes
Read Results

Dengue IgG/IgM



EXCELLENT CLINICAL COMPARISON RESULT

Items	Sensitivity	Specificity
Dengue NS1 Antigen Test	95.8%	96.1%
Dengue IgG/IgM Antibody Test	94.3%	98.3%

The packing details



CB	
QTY :	360 PCS
G.W :	12 KG
MEAS :	46X44X44CM



Shenzhen Emeng Health Technology Co., LTD

+86 13760822077

kaylayi0612@foxmail.com

oxygenconcentratormedical.com

Room D, 7th Floor, Guang Long Building, 162 Ping Xin North Road, PingHu Street, LongGang District, ShenZhen