

Shenzhen Emeng Health Technology Co., LTD oxygenconcentratormedical.com

98% Accuracy Dengue Ns1 Ag & IgG/IgM Antigen Rapid Test Kit with CE and ISO Certifications

Basic Information

Place of Origin: chinaBrand Name: DVOTCertification: CE/ISO

Model Number: Dengue Ns1 Ag & lgG/lgM Combo Test Kit

 Minimum Order Quantity: 360 test

• Price: 0.9-1 USD

• Packaging Details: 360T/PER CTN Weight :11KGS CTN SIZE

:46*44*44CM

• Delivery Time: 15-20 work days

• Payment Terms: D/P, T/T, Western Union, MoneyGram

• Supply Ability: 1000000pcs



Product Specification

• Type: Rapid Test Kit

• Accuracy: 98%

Keyword: Antigen Rapid Test

• Time Of Use: Disposable

• Result Time: 15 Minutes To Get

Test Time: 15 MinutesCertifications: CE ISO

• Highlight: ce clungene antigen rapid test,

ce medical diagnostic test kits,

disposable clungene antigen rapid test



98% Accuracy Dengue Ns1 Ag & IgG/IgM Antigen Rapid Test Kit with CE and ISO Certifications

Product Description:

Dengue fever is an acute insect-borne infectious disease caused by dengue virus, which is mainly transmitted by Aedes aegypti and Aedes albopictus mosquitoes. The incubation period of dengue fever is usually 1-14 days, mostly 5-9 days. Patients usually have an acute onset of illness, with the first symptom being sudden onset of high fever, which may be accompanied by chills, and a temperature of up to 40°C within 24 hours. It usually lasts for 3~7 days. In some patients, the temperature may drop to normal after 3-5 days of fever, but rise again after 1 day, a phenomenon known as bimodal or saddle fever. If you are not sure whether you are infected with dengue virus, you can use this antigen test.

Our Infectious Disease Rapid Test Kits have been certified by CE and ISO, ensuring that they meet the highest quality and safety standards. The Antigen Rapid Test is a disposable test, making it convenient and hygienic for both patients and healthcare professionals.

Our Antigen Rapid Test is suitable for use in a variety of settings, including hospitals, clinics, and laboratories. It can also be used in non-medical settings such as workplaces, schools, and airports to help prevent the spread of infectious diseases.

Overall, our Antigen Rapid Test is a reliable and efficient way to detect the presence of infectious agents in the body. With its quick test time of just 15 minutes, CE and ISO certifications, and convenient disposable design, it is an excellent choice for anyone in need of an accurate and easy-to-use diagnostic test.

Technical Parameters:

Product Attribute	Value	
Product Name	Diagnostic Antigen Test Kit	
Keyword	Antigen Home Test Kit, Antigen Rapid Test Kit	
Туре	Rapid Test Kit	
Time of Use	Disposable	
Result Time	15 Minutes	
Accuracy	98%	
Certifications	CE, ISO	

Applications:

The DVOT Dengue Ns1 Ag & IgG/IgM Combo Test Kit is a Rapid Test Kit that can be used on patients suspected of having Dengue fever. It is also suitable for use in hospitals, clinics, and laboratories. This test kit can detect both early and late-stage infections, which makes it a valuable tool in the diagnosis of Dengue virus.

One of the key features of this Antibody Rapid Test Kit is that it is disposable. This means that it can be used only once, which eliminates the risk of contamination and ensures accurate results. The DVOT Dengue Ns1 Ag & IgG/IgM Combo Test Kit is also easy to use, and it does not require any special training or equipment.

The DVOT Dengue Ns1 Ag & IgG/IgM Combo Test Kit is suitable for a wide range of product application occasions and scenarios. For example, it can be used in hospitals and clinics as a diagnostic tool for patients suspected of having Dengue fever. It can also be used by travelers who have been to areas where Dengue virus is prevalent and want to check if they have been infected.

Furthermore, the DVOT Dengue Ns1 Ag & IgG/IgM Combo Test Kit can be used in research settings to identify the prevalence of Dengue virus in a particular population. This Antibody Diagnostic Kit can also be used by public health officials to monitor the spread of Dengue virus and to implement appropriate control measures.

In summary, the DVOT Dengue Ns1 Ag & IgG/IgM Combo Test Kit is a valuable tool in the diagnosis and monitoring of Dengue virus. It is an accurate and reliable Antibody Rapid Test Kit that can deliver results in just 15 minutes. This Home Testing Kit can be used in a wide range of product application occasions and scenarios, making it a versatile and valuable tool in the fight against Dengue fever.

Support and Services:

Our Medical Diagnostic Test Kits are designed to provide accurate and reliable results for a variety of medical conditions. Our team is available to assist with any questions or concerns you may have regarding our products.

Steps to use:

Preparation Stage:

Ensure that your hands are clean and that you have the required blood collection kit (e.g., blood collection needle, alcohol pads, etc.). Remove the test card from the box and check that the package is intact and that it is within the expiration date. Follow the instructions on the package to confirm the type of sample and sample size required for the test card.

Collect the blood sample:

Clean the blood collection site (usually the tip of the finger) with an alcohol pad and allow the alcohol to evaporate.

Use a blood collection needle to gently prick the tip of the finger and collect the appropriate amount of blood sample. Take care to avoid contaminating the sample.

Perform the test:


Step 1: Add a drop of blood sample to the designated area of the test card.

Step 2: Wait a certain amount of time (usually 15 minutes) according to the instructions so that the test card can respond adequately.

Step 3: After completing the wait, observe the results on the test card.

Read the results:

Observe the indicated areas or lines on the test card and judge the results according to the explanations provided in the instructions. Take care to distinguish between the different results and their significance for Dengue Negative, Dengue Positive and Dengue NS1 Test.



FEATURE

Accuracy, visual result, simple and safe to use

SPECIMEN

Whole Blood/Serum/Plasma

REACTION TIME

15mins

STORAGE TEMPERATURE

4~30°C


SHELF LIFE

24 months




TEST PROCEDURE






· For whole blood sample testing







EXCELLENT CLINICAL COMPARISON RESULT

Items	Sensitivity	Specificity
Dengue NS1 Antigen Test	95.8%	96.1%
Dengue IgG/IgM Antibody Test	94.3%	98.3%

The packing details



СВ				
QTY:	360 PCS			
G.W:	12 KG			
MEAS:	46X44X44CM			

Shenzhen Emeng Health Technology Co., LTD



+86 13760822077



kaylayi0612@foxmail.com



Oxygenconcentratormedical.com

Room D, 7th Floor, Guang Long Building, 162 Ping Xin North Road, PingHu Street, LongGang District, ShenZhen