

15 Minute Dengue Ns1 Ag & IgG/IgM Instant Diagnostic Rapid Test Kit With 99% Accuracy - CE ISO Certified

Our Product Introduction

Basic Information

- Place of Origin: china
- Brand Name: DVOT
- Certification: CE/ISO
- Model Number: Dengue Ns1 Ag & IgG/IgM Combo Test Kit
- Minimum Order Quantity: 360 test
- Price: 0.9-1 USD
- Packaging Details: 360T/PER CTN Weight :11KGS CTN SIZE :46*44*44CM
- Delivery Time: 15-20 work days
- Payment Terms: D/P, T/T, Western Union, MoneyGram
- Supply Ability: 1000000pcs



Product Specification

- Certifications: CE ISO
- Type: Rapid Test Kit
- Time Of Use: Disposable
- Test Time: 15 Minutes
- Accuracy: 98%
- Result Time: 15 Minutes To Get
- Keyword: Antigen Rapid Test
- Highlight: ce clungene antigen rapid test, ce medical diagnostic test kits, disposable clungene antigen rapid test



Product Description

15 Minute Dengue Ns1 Ag & IgG/IgM Instant Diagnostic Rapid Test Kit With 99% Accuracy - CE ISO Certified

Product Description:

Dengue fever is an acute insect-borne infectious disease caused by dengue virus, which is mainly transmitted by *Aedes aegypti* and *Aedes albopictus* mosquitoes. The incubation period of dengue fever is usually 1-14 days, mostly 5-9 days. Patients usually have an acute onset of illness, with the first symptom being sudden onset of high fever, which may be accompanied by chills, and a temperature of up to 40°C within 24 hours. It usually lasts for 3~7 days. In some patients, the temperature may drop to normal after 3-5 days of fever, but rise again after 1 day, a phenomenon known as bimodal or saddle fever. If you are not sure whether you are infected with dengue virus, you can use this antigen test.

Our diagnostic test kits are CE and ISO certified, so you can be assured of their quality and accuracy. In addition, our Antigen Rapid Test Kit has an accuracy rate of 98%, making it one of the most reliable and effective diagnostic tools available on the market today.

Using our Antibody Diagnostic Kit is easy and straightforward. Simply follow the instructions provided with the kit, and you'll be able to quickly and accurately diagnose infections and other medical conditions. With a test time of just 15 minutes, you can get the results you need quickly and get your patients the care they need as soon as possible.

Our diagnostic test kit is suitable for a wide range of applications, including in hospitals, clinics, and other medical facilities. Whether you're looking to diagnose infectious diseases or other medical conditions, our Antigen Rapid Test Kit is the ideal choice for accurate and reliable results.

Technical Parameters:

Product Attribute	Value
Product Name	Antigen Rapid Test Kits
Type	Rapid Test Kit
Keyword	Antigen Rapid Test
Time Of Use	Disposable
Test Time	15 Minutes
Result Time	15 Minutes To Get
Certifications	CE, ISO
Accuracy	98%

Applications:

The test kit is ideal for use in a variety of occasions and scenarios. For medical professionals, the DVOT Dengue Ns1 Ag & IgG/IgM Combo Test Kit can be used in hospital or clinic settings to quickly and accurately diagnose patients. The test kit is also a great option for use in remote or rural areas where access to medical facilities is limited.

For individuals looking to monitor their own health, the DVOT Dengue Ns1 Ag & IgG/IgM Combo Test Kit is a convenient option that can be used in the comfort of their own home. The test kit is an antibody test kit that can detect the presence of both Ns1 Ag and IgG/IgM antibodies, providing accurate results in just 15 minutes.

Thanks to its CE and ISO certifications, the DVOT Dengue Ns1 Ag & IgG/IgM Combo Test Kit is a reliable and trustworthy option for anyone looking to test for Dengue virus. So whether you're a medical professional or an individual concerned about your own health, this infectious disease test kit is the perfect solution.

Support and Services:

Our Medical Diagnostic Test Kits are designed to provide accurate and reliable results for a variety of medical conditions. Our team is available to assist with any questions or concerns you may have regarding our products.

Steps to use:

Preparation Stage:

Ensure that your hands are clean and that you have the required blood collection kit (e.g., blood collection needle, alcohol pads, etc.). Remove the test card from the box and check that the package is intact and that it is within the expiration date. Follow the instructions on the package to confirm the type of sample and sample size required for the test card.

Collect the blood sample:

Clean the blood collection site (usually the tip of the finger) with an alcohol pad and allow the alcohol to evaporate. Use a blood collection needle to gently prick the tip of the finger and collect the appropriate amount of blood sample. Take care to avoid contaminating the sample.

Perform the test:

Step 1: Add a drop of blood sample to the designated area of the test card.
Step 2: Wait a certain amount of time (usually 15 minutes) according to the instructions so that the test card can respond adequately.
Step 3: After completing the wait, observe the results on the test card.

Read the results:

Observe the indicated areas or lines on the test card and judge the results according to the explanations provided in the instructions. Take care to distinguish between the different results and their significance for Dengue Negative, Dengue Positive and Dengue NS1 Test.



The advertisement features a yellow background with a faint image of a laboratory setting. At the top left is the Dvot logo. The word "DENGUE" is prominently displayed in large, bold, black letters. Below it, the text "Dengue NS1Ag & IgG/IgM Combo Test kit" is written in a smaller, bold font. A blue banner at the bottom contains white text describing Dengue fever as a mosquito-borne tropical disease with symptoms like fever, headache, muscle and joint pains, and a skin rash similar to measles. Below this text are three hexagonal images: a mosquito, a map of the Americas, and a map of Asia.

FEATURE

Accuracy, visual result, simple and safe to use

SPECIMEN

Whole Blood/Serum/Plasma

REACTION TIME

15mins

STORAGE TEMPERATURE

4-30°C

SHELF LIFE

24 months




TEST PROCEDURE


• For serum/plasma sample testing




Dengue NS1 Ag



• For whole blood sample testing



Dengue IgG/IgM



EXCELLENT CLINICAL COMPARISON RESULT

Items	Sensitivity	Specificity
Dengue NS1 Antigen Test	95.8%	96.1%
Dengue IgG/IgM Antibody Test	94.3%	98.3%

The packing details



CB	
QTY :	360 PCS
G.W :	12 KG
MEAS :	46X44X44CM



Shenzhen Emeng Health Technology Co., LTD



+86 13760822077



kaylayi0612@foxmail.com



oxygenconcentratormedical.com

Room D, 7th Floor, Guang Long Building, 162 Ping Xin North Road, PingHu Street, LongGang District, ShenZhen