



EMY06 Hyperbaric Chamber for Comprehensive Care in 1.2-1.3ATA Oxygen Therapy

Our Product Introduction

For more products please visit us on oxygenconcentratormedical.com

Basic Information

- Place of Origin: china
- Brand Name: O2EMENG
- Certification: CE/FDA/FC/ROHS/UKCA/SGS
- Model Number: EMY06
- Minimum Order Quantity: 1PC
- Price: \$25000
- Packaging Details: wooden packing
- Delivery Time: 15-20 work days
- Payment Terms: D/P, T/T, Western Union, MoneyGram
- Supply Ability: 100PCS



Product Specification

- Model: EMY06
- Cabin Weight: 1100KGS
- Cabin Size: 200*170*190CM
- Function: Hyperbaric Oxygen Therapy
- Working Pressure: 1.2-1.3ATA
- Usage: Medical/Therapeutic
- Highlight: **1.3ata hyperbaric oxygen chamber, 1.3ata oxygen hyperbaric chamber, medical hyperbaric oxygen chamber**



More Images



Product Description

EMY06 Hyperbaric Chamber for Comprehensive Care in 1.2-1.3ATA Oxygen Therapy

Product Description:

Our Hyperbaric Oxygen Chamber is made from high-quality materials, ensuring that it is durable and long-lasting. With a cabin weight of 1100KGS, it is built to withstand repeated use, making it perfect for clinics and hospitals that require frequent treatments.

The EMY06 Hyperbaric Oxygen Chamber is perfect for individuals who are seeking a natural way to recover from injuries and illnesses. HBOT has been proven to be effective in treating a variety of conditions, such as traumatic brain injury, autism, and even stroke. With our product, you can now enjoy the benefits of HBOT in the comfort of your own home, or as part of your clinic or hospital's treatment program.

The EMY06 Hyperbaric Oxygen Chamber is easy to use, and you can expect to see results within a few sessions. The therapy works by increasing the amount of oxygen in your body, which helps to promote healing and repair damaged tissues. It also helps to reduce inflammation, which is a common cause of pain and discomfort.

Our product is designed with your safety in mind, and we have implemented several safety features to ensure that you have a safe and comfortable experience. The chamber is equipped with a ventilation system that helps to regulate the temperature and humidity levels inside the cabin, ensuring that you are always comfortable during your therapy sessions.

Technical Parameters:

Usage:	Medical/Therapeutic
Model:	EMY06
Cabin Weight:	1100KGS
Cabin size:	200*170*190CM
Working pressure:	1.2-1.3ATA
Function:	Hyperbaric Oxygen Therapy

Advantages:

The O2EMMENG EMY06 model hyperbaric chamber has several significant differences and advantages over traditional hyperbaric chambers:

- Advanced Technology and Control Systems:** The EMY06 model utilizes the latest in hyperbaric technology with a precise oxygen concentration and pressure control system. Traditional hyperbaric chambers may use more basic controls that do not provide the same precise regulation and monitoring capabilities.
- New Hyperbaric Oxygen Equipment:** Intelligent all-in-one machine, unlike traditional equipment, does not require the operation of three machines, Oxygen Equipment + Aerodynamic Equipment + Air Cooler, realizing easy one-button operation and saving time and effort.
- Comfort and User Experience:** The EMY06 model is designed with a more spacious and comfortable interior space, taking into account the patient's comfort and convenience during treatment.
- Durability and Reliability:** The EMY06 model undergoes rigorous quality control and testing to ensure stable operation over long periods of time without frequent maintenance.
- Technological innovations and safety:** The EMY06 model integrates new technological innovations that enhance safety and therapeutic efficacy. It allows for more precise control of oxygen delivery to support medical treatments more effectively.

Support and Services:

- ✓ **Demand customization:** according to your use environment and expectations, recommend the optimal core products and life aesthetics composition scheme.
- ✓ **Exclusive design:** the advanced technology team and the top design team jointly, from the size and color material

matching, the actual program customization landing.

✓**Craftsmanship Manufacturing:** The core components and other materials of hyperbaric chamber series products are manufactured in aerospace material factories.

✓**Lifetime maintenance:** provide perfect professional after-sales service guarantee



CHAMBER CAPSULE DESCRIPTION



1. Automatic digital display real-time temperature, humidity, oxygen purity and pressure.
2. Rotary timer system, use of time can be pressed in advance Multiple noise reduction technology.
3. Excellent Silence performance, the noise is ≤ 63 dB. Offer you a peaceful experience environment.
4. Double pressure relief device Spiral step less super smooth rotary valve. Spiral ergonomic hand-held design. Pressure relief is smoother and more convenient. Internal and external two-way integrated polycarbonate valve. One valve, two-way operation.
5. Automatic pressurize relief valves
Two automatic high precision constant pressure valve. Ensure the product is safe and effective by Aluminum alloy steel, high hardness, wear resistance.
6. Communicating system.
Intercom telephone is convenient for users to talk with operator outside.
7. Machine start-up protection system.
The chamber only can be turned on when the chamber door is completely closed, avoiding the danger of the door being left unlocked.
8. Automatic pressurize and depressurize system.
Based on the principle of "Cylinder (engine)" , Pressurize automatically after the door is completely closed.
9. Made by Aluminum alloy steel; With environmental spray paint;
Strong weight bearing capacity, Durable and wear-resistant; non-deformable; long lifespan.
10. In case of unexpected power failure, the door can be opened by emergency safety from inside the chamber.



HARD TYPE CHAMBER PRODUCT DISPLAY



EMY06



(Two people) Hard Type
Hyperbaric Oxygen Chamber

EMY03



(Two people) Hard Type
Hyperbaric Oxygen Chamber

EMY05



(One person) Hard Type
Hyperbaric Oxygen Chamber

EMY01



(One person) Hard Type
Hyperbaric Oxygen Chamber

CONCLUSION



It is more feasible to increase the dissolved oxygen in the blood. The hyperbaric oxygen chamber can maintain the oxygen partial pressure within a safe and effective range by flexibly adjusting the oxygen concentration or pressure, and gently and effectively improve the oxygen delivery to tissues and organs by increasing dissolved oxygen.



APPLICATION SCENARIO




MEDICAL INSTITUTION



OXYGEN BAR CLUB



ELDERLY CARE CENTER



FITNESS INDUSTRY



BEAUTY INDUSTRY



HOTEL



PLATEAU AREA



INDIVIDUAL FAMILY



WORK PLACE



Shenzhen Emeng Health Technology Co., LTD

+86 13760822077

kaylayi0612@foxmail.com

oxygenconcentratormedical.com

Room D, 7th Floor, Guang Long Building, 162 Ping Xin North Road, PingHu Street, LongGang District, ShenZhen