



Shenzhen Emeng Health Technology Co., LTD

oxygenconcentratormedical.com


EMY06 Hyperbaric Chamber Control 1.3ATA Double Hard Seat Type for Health Care

Our Product Introduction

For products please visit us on oxygenconcentratormedical.com

Basic Information

- Place of Origin: china
- Brand Name: O2EMENG
- Certification: CE/FDA/FC/ROHS/UKCA/SGS
- Model Number: EMY06
- Minimum Order Quantity: 1PC
- Price: \$25000
- Packaging Details: wooden packing
- Delivery Time: 15-20 work days
- Payment Terms: D/P, T/T, Western Union, MoneyGram
- Supply Ability: 100PCS



Product Specification

- Working Pressure: 1.2-1.3ATA
- Cabin Size: 200*170*190CM
- Usage: Medical/Therapeutic
- Function: Hyperbaric Oxygen Therapy
- Cabin Weight: 1100KGS
- Model: EMY06
- Highlight: Health Care Hyperbaric Chamber Control , Double Hard Seat Hyperbaric Chamber Control , 1.3ATA Hyperbaric Chamber Control

Product Description

Medical 1.2-1.3ATA Hyperbaric Oxygen Chamber EMY06,1-3 person use

Product Description:

The Hyperbaric Oxygen Chamber is equipped with a Compact Portable Oxygen Concentrator, which ensures that patients receive an optimal level of oxygen therapy while inside the chamber. This concentrator is designed to deliver a high flow of oxygen to patients, ensuring that they receive the right amount of oxygen to improve their health. The oxygen concentrator is easy to operate and is designed to deliver a steady supply of oxygen to patients throughout the therapy session.

In addition to the oxygen concentrator, the Hyperbaric Oxygen Chamber is also equipped with a High Flow Oxygen Concentrator. This concentrator is designed to deliver a high flow of oxygen to patients, ensuring that they receive the right amount of oxygen to improve their health. The concentrator is designed to operate quietly, ensuring that patients can relax and enjoy their therapy session without any distractions.

The Hyperbaric Oxygen Chamber is a heavy-duty device that is built to last. It has a cabin weight of 1100KGS, making it one of the durable and sturdy HBOT chambers available in the market. It is designed to offer superior performance and is perfect for patients who need a high level of oxygen therapy to improve their health.

Technical Parameters:

Model	EMY06
Cabin Weight	1100KGS
Usage	Medical/Therapeutic
Working Pressure	1.2-1.3ATA
Cabin Size	200*170*190CM
Function	Hyperbaric Oxygen Therapy

What are they used to treat?

1. Asthma, bronchitis, pneumonia, acute carbon monoxide poisoning, post-stroke, cerebral hemorrhage, arteriosclerosis, diabetes, heart disease.
2. Treat Alzheimer's disease, stroke, muscle paralysis, neurological disorders, autism, schizophrenia, headache and migraine.
3. Treat insomnia, improve poor blood circulation and prevent Alzheimer's disease.
4. Therapeutic effect on diabetic retinitis and keratitis.
5. Promote swelling reabsorption stimulation, cell regeneration after micro-integration surgery, promote muscle and bone wound healing.
6. Treatment of traumatic infections, acute soft tissue injuries, ulcers, burns, wounds, fractures, etc.
7. Acne, herpes simplex, herpes zoster, cicatricial alopecia, dermatologic treatment, etc.
8. auxiliary therapeutic effect on infertility, and obvious auxiliary therapeutic effect on postpartum recovery and postpartum restoration.
9. It has obvious therapeutic effect on the rehabilitation of children with congenital cranial nerve defects and congenital brain development problems.
10. Accelerates skin cell repair and regeneration, improves skin cell activity, promotes lipolytic enzyme activity, and reduces body weight.

Advantages:

The O2EMMENG EMY06 model hyperbaric chamber has several significant differences and advantages over traditional hyperbaric chambers:

1. **Advanced Technology and Control Systems:** The EMY06 model utilizes the latest in hyperbaric technology with a precise oxygen concentration and pressure control system. Traditional hyperbaric chambers may use more basic controls that do not provide the same precise regulation and monitoring capabilities.

2. New Hyperbaric Oxygen Equipment: Intelligent all-in-one machine, unlike traditional equipment, does not require the operation of three machines, Oxygen Equipment + Aerodynamic Equipment + Air Cooler, realizing easy one-button operation and saving time and effort.

3. Comfort and User Experience: The EMY06 model is designed with a more spacious and comfortable interior space, taking into account the patient's comfort and convenience during treatment.

4. Durability and Reliability: The EMY06 model undergoes rigorous quality control and testing to ensure stable operation over long periods of time without frequent maintenance.

5. Technological innovations and safety: The EMY06 model integrates new technological innovations that enhance safety and therapeutic efficacy. It allows for more precise control of oxygen delivery to support medical treatments more effectively.

Support and Services:

✓**Demand customization:** according to your use environment and expectations, recommend the optimal core products and life aesthetics composition scheme.

✓**Exclusive design:** the advanced technology team and the top design team jointly, from the size and color material matching, the actual program customization landing.

✓**Craftsmanship Manufacturing:** The core components and other materials of hyperbaric chamber series products are manufactured in aerospace material factories.

✓**Lifetime maintenance:** provide perfect professional after-sales service guarantee



CHAMBER CAPSULE DESCRIPTION



1. Automatic digital display real-time temperature, humidity, oxygen purity and pressure.
2. Rotary timer system, use of time can be pressed in advance Multiple noise reduction technology.
3. Excellent Silence performance, the noise is $\leq 63\text{dB}$. Offer you a peaceful experience environment.
4. Double pressure relief device Spiral step less super smooth rotary valve. Spiral ergonomic hand-held design. Pressure relief is smoother and more convenient. Internal and external two-way integrated polycarbonate valve. One valve, two-way operation.
5. Automatic pressurize relief valves
Two automatic high precision constant pressure valve. Ensure the product is safe and effective by Aluminum alloy steel, high hardness, wear resistance.
6. Communicating system.
Intercom telephone is convenient for users to talk with operator outside.
7. Machine start-up protection system.
The chamber only can be turned on when the chamber door is completely closed, avoiding the danger of the door being left unlocked.
8. Automatic pressurize and depressurize system.
Based on the principle of "Cylinder (engine)" , Pressurize automatically after the door is completely closed.
9. Made by Aluminum alloy steel; With environmental spray paint;
Strong weight bearing capacity, Durable and wear-resistant; non-deformable; long lifespan.
10. In case of unexpected power failure, the door can be opened by emergency safety from inside the chamber.



HARD TYPE CHAMBER PRODUCT DISPLAY

● ● ●

EMY06



(Two people) Hard Type
Hyperbaric Oxygen Chamber

EMY03



(Two people) Hard Type
Hyperbaric Oxygen Chamber

EMY05



(One person) Hard Type
Hyperbaric Oxygen Chamber

EMY01



(One person) Hard Type
Hyperbaric Oxygen Chamber

CONCLUSION



It is more feasible to increase the dissolved oxygen in the blood.
The hyperbaric oxygen chamber can maintain the oxygen partial pressure within a safe and effective range by flexibly adjusting the oxygen concentration or pressure, and gently and effectively improve the oxygen delivery to tissues and organs by increasing dissolved oxygen.



APPLICATION SCENARIO



MEDICAL INSTITUTION



OXYGEN BAR CLUB



ELDERLY CARE CENTER



FITNESS INDUSTRY



BEAUTY INDUSTRY



HOTEL



PLATEAU AREA



INDIVIDUAL FAMILY



WORK PLACE



Shenzhen Emeng Health Technology Co., LTD



+86 13760822077



kaylayi0612@foxmail.com



oxygenconcentratormedical.com

Room D, 7th Floor, Guang Long Building, 162 Ping Xin North Road, PingHu Street, LongGang District, ShenZhen