



## Hyperbaric Oxygen Therapy Medical Oxygen Concentrator with Less Than 63 DB Noise Level White

### Our Product Introduction

For products please visit us on [oxygenconcentratormedical.com](http://oxygenconcentratormedical.com)

#### Basic Information

- Place of Origin: China
- Brand Name: Tyleneolmed
- Certification: CE,ISO13485,SGS ,FCC,RoHs
- Model Number: ND-ZY10SL
- Minimum Order Quantity: Negotiable
- Price: Negotiable
- Packaging Details: OPP packing ,1pcs per Carton. Carton size:540\*370\*695mm,Gross Weight:21kgs.
- Delivery Time: 7-10 Working days
- Payment Terms: T/T, Western Union
- Supply Ability: 1000pcs



#### Product Specification

- Product Name: O2EMENG Medical Therapy Machine  
1.3ATA Hyperbaric Oxygen Chamber For Healing
- Chamber Colour: White /Black
- Cabin Weight: 1000KGS
- Noise Level: Less Than 63 DB
- Chamber Oxygen Concentration: 28% 32%
- Flow Rate: 10Liter/Min
- Product Category: Medical Oxygen Concentrator
- Model: EMY01
- Highlight: **Hyperbaric oxygen therapy concentrator ,  
Medical oxygen concentrator with warranty,  
Quiet oxygen concentrator under 63 DB**

## Product Description

### Product Description:

A Therapeutic Oxygen Generator is an essential medical device that ensures a continuous and reliable supply of oxygen for patients who require oxygen therapy. The Medical Oxygen Concentrator is a state-of-the-art machine designed to provide high-quality oxygen for individuals with respiratory conditions. This Clinic Oxygen Producer is equipped with advanced features to deliver the necessary oxygen levels efficiently and effectively.

With an input power of 700VA, this Medical Oxygen Concentrator is capable of operating efficiently while consuming low power. The device has a flow rate of 10 liters per minute, ensuring a steady and consistent supply of oxygen to meet the patient's needs. The Chamber Oxygen Concentration of this oxygen concentrator ranges from 28% to 32%, making it suitable for a wide range of therapeutic applications.

One of the key features of this Medical Oxygen Generator is its low noise level, measuring less than 63 dB during operation. This ensures a quiet and peaceful environment for patients while the device is in use. Additionally, the Power Consumption of the oxygen concentrator is 1500w, making it energy-efficient and cost-effective to operate.

The Medical Oxygen Concentrator is designed to be user-friendly, with intuitive controls and a clear display that allows for easy monitoring of oxygen levels. Its compact and portable design makes it ideal for use in various healthcare settings, including hospitals, clinics, and homecare environments.

Patients who rely on oxygen therapy can benefit greatly from the reliable performance and advanced features of this Therapeutic Oxygen Generator. The device is engineered to deliver a constant and steady flow of oxygen, helping patients manage their respiratory conditions more effectively and improve their overall quality of life.

Whether used in a clinical setting or at home, the Medical Oxygen Concentrator provides a dependable source of oxygen that patients can trust. Its high chamber oxygen concentration ensures that patients receive the necessary oxygen levels to support their respiratory needs, promoting better health outcomes and enhanced well-being.

In conclusion, the Medical Oxygen Concentrator is a top-of-the-line device that offers exceptional performance and reliability for individuals in need of oxygen therapy. Its advanced features, including a high flow rate, chamber oxygen concentration, low noise level, and energy-efficient operation, make it a preferred choice for healthcare professionals and patients alike. Invest in this Clinic Oxygen Producer to ensure a consistent and effective supply of therapeutic oxygen for improved patient care and well-being.

### Features:

Product Name: Medical Oxygen Concentrator

Cabin Weight: 1000KGS

Function: Hyperbaric Oxygen Therapy

Flow Rate: 10Liter/Min

Noise Level: Less Than 63 DB

Rating: 1500W

Keywords: Oxygen Concentration Machine, Therapeutic Oxygen Generator, Clinic Oxygen Producer

### Technical Parameters:

Flow Rate	10 Liter/Min
Color	White, Black
Rating	1500W
Chamber Colour	White / Black
Function	Hyperbaric Oxygen Therapy
Chamber Oxygen Concentration	28% 32%
Cabin Weight	1000KGS
Model	EMY01
Usage	Medical/Therapeutic
Power Consumption	1500W

### Applications:

Originating from China, the Tylenolmed ND-ZY10SL is suitable for a wide range of application occasions and scenarios. It is ideal for use in hospitals, clinics, nursing homes, and even in-home healthcare settings. The therapeutic oxygen generated by this device is crucial for patients with respiratory conditions such as COPD, asthma, and pneumonia.

The minimum order quantity and price for the Tylenolmed ND-ZY10SL Medical Oxygen Concentrator are negotiable, making it accessible for various healthcare facilities and individuals. Payment terms include T/T and Western Union, providing flexibility for purchasing.

With a supply ability of 1000pcs and a delivery time of 7-10 working days, the Tylenolmed ND-ZY10SL ensures a steady and timely availability of this essential healthcare oxygen separator. The packaging details include OPP packing with 1 piece per carton. The carton size is 540\*370\*695mm, with a gross weight of 21kgs.

The chamber of the Tylenolmed ND-ZY10SL comes in white or black, offering options to suit different preferences and environments. The model, EMY01, features a chamber oxygen concentration of 28% to 32%, meeting the therapeutic oxygen needs of patients.

Operating at voltages of 110V/50Hz and 220V/50Hz, the O2EMENG Medical Therapy Machine 1.3ATA Hyperbaric Oxygen Chamber, powered by the Tylenolmed ND-ZY10SL, is a reliable and effective solution for delivering therapeutic oxygen. It is a must-have product for healthcare facilities and individuals looking to enhance patient care and well-being.

### Support and Services:

Our Product Technical Support and Services for the Medical Oxygen Concentrator include:

- Troubleshooting assistance for any issues related to the oxygen concentrator
- Guidance on setup and proper usage of the device
- Product maintenance tips and recommendations
- Information on cleaning and care of the oxygen concentrator
- Assistance with product warranties and service agreements

### Packing and Shipping:

Product Packaging for Medical Oxygen Concentrator:

The medical oxygen concentrator will be securely packaged in a durable cardboard box with foam inserts to prevent damage during transit. The packaging will also include clear instructions for setup and use.

Shipping Information:

Once your order is processed, the medical oxygen concentrator will be carefully packed and shipped using a reputable courier service. You will receive a tracking number to monitor the shipment status and ensure a timely delivery.

### FAQ:

**Q: What is the brand name of this Medical Oxygen Concentrator?**

A: The brand name is Tylenolmed.

**Q: What is the model number of this Medical Oxygen Concentrator?**

A: The model number is ND-ZY10SL.

**Q: What certifications does this Medical Oxygen Concentrator have?**

A: It is certified with CE, ISO13485, SGS, FCC, and RoHs.

**Q: Where is this Medical Oxygen Concentrator manufactured?**

A: It is manufactured in China.

**Q: What are the payment terms accepted for purchasing this Medical Oxygen Concentrator?**

A: The accepted payment terms are T/T and Western Union.

O2EMENG

# HYPERBARIC OXYGEN CHAMBER

- IMPROVE BRAIN BLOOD CIRCULATION
- SPEED UP METABOLISM
- HIGH PRESSURE OXYGEN INHALATION



## CHAMBER CAPSULE DESCRIPTION



1. Automatic digital display real-time temperature, humidity, oxygen purity and pressure.
2. Rotary timer system, use of time can be pressed in advance Multiple noise reduction technology.
3. Excellent Silence performance, the noise is  $\leq 63\text{dB}$ . Offer you a peaceful experience environment.
4. Double pressure relief device Spiral step less super smooth rotary valve. Spiral ergonomic hand-held design. Pressure relief is smoother and more convenient. Internal and external two-way integrated polycarbonate valve. One valve, two-way operation.
5. Automatic pressurize relief valves  
Two automatic high precision constant pressure valve. Ensure the product is safe and effective by Aluminum alloy steel, high hardness, wear resistance.
6. Communicating system.  
Intercom telephone is convenient for users to talk with operator outside.
7. Machine start-up protection system.  
The chamber only can be turned on when the chamber door is completely closed, avoiding the danger of the door being left unlocked.
8. Automatic pressurize and depressurize system.  
Based on the principle of "Cylinder (engine)" , Pressurize automatically after the door is completely closed.
9. Made by Aluminum alloy steel; With environmental spray paint;  
Strong weight bearing capacity, Durable and wear-resistant; non-deformable; long lifespan.
10. In case of unexpected power failure, the door can be opened by emergency safety from inside the chamber.





## WHY ARE MORE AND MORE PEOPLE CHOOSING HYPERBARIC OXYGEN CHAMBERS?

In some parts of the body, the capillaries are so narrow that blood cells cannot pass through them, bound oxygen cannot reach them, and they can only rely on dissolved oxygen. The oxygen generator has little change in the dissolved oxygen in the blood, which is also the reason why the ordinary oxygen absorption effect in the plateau area is general.



## CONCLUSION



It is more feasible to increase the dissolved oxygen in the blood.  
The hyperbaric oxygen chamber can maintain the oxygen partial pressure within a safe and effective range by flexibly adjusting the oxygen concentration or pressure, and gently and effectively improve the oxygen delivery to tissues and organs by increasing dissolved oxygen.



# HARD TYPE HYPERBARIC OXYGEN CHAMBER



Innovative technology, quality assurance



The cabin door is made of aviation PC material

Multiple safety technologies  
One-button pressure relief protection system during power outage

AI intelligent voice control system



## OUR CERTIFICATE



## SHIPPING AND DELIVERY



## FAQ

### 1. How does hyperbarics work?

Hyperbarics is a technology in which the air pressure in the environment is increased. When a person's body is placed in a more pressure environment, it absorbs more oxygen molecules per volume of compressed air.

The body normally transports oxygen via the hemoglobin of the red blood cells. By increasing the air pressure, oxygen is then driven into the body's fluids, allowing a super-saturation of tissues and organs with oxygen. The increased pressure infuses the body with oxygen, even reaching injuries with damaged circulation.

### 2. What are the health benefits of hyperbaric oxygen chambers?

Studies have shown that hyperbaric oxygen chambers can treat hypoxic diseases, relieve fatigue, delay aging, prevent cardiovascular disease, improve sleep quality, accelerated wound repair and enhance immunity. We have customers serviced shows hyperbaric oxygen therapy makes very positive results.

### 3. How long is the treatment?

Usually, we recommend our clients receive the HBOT therapy 5-6 sessions per week, 60-90mins per session. Using frequency or duration can be adjusted according to individual circumstances.

### 4. Is HBOT safe?

Hyperbaric oxygen therapy is prescribed by a physician and performed under medical supervision. Although there are minor risks like all medical treatments, overall hyperbaric oxygen therapy is extremely safe.

### 5. Is there anyone who can not use our chamber?

Pregnant ladies, patients who catch a cold or Flu, or those who have rhinitis or otitis media are not recommended to use the chamber.

### 6. How do you feel inside a hyperbaric oxygen chamber?

When lying in the oxygen chamber, you will feel comfortable and at ease. Also, You may feel a little squeezing of the ears because of the pressure. It is quite normal and can be relieved by chewing or swallowing.

### 7. After bought the chamber, do I need to by any other things?

No, the price of our chamber includes the chamber capsule, all machines, namely the air compressor, oxygen concentrator, air cooler, and accessories like frame, tubes, etc, so you don't have to buy anything else.

### 8. Is it complex to install the chamber or change the filter?

Of course not, we will provide you with a video to show you how to install the chamber and change the filter.

We also support online services to help you with all the issues, so don't worry.

### 9. What is the general time warranty expected for ZHIKANG oxygen chambers?

We will provide one-year warranty and lifelong after-sales services.

If any quality problem/fault in material/design is under correct operation within one year,

If it's easy to fix, we will send new components for free and teach you how to repair them.

(Even after the one-year warranty period, you can also contact us freely and we will be glad to help you.)



**Shenzhen Emeng Health Technology Co., LTD**



+86 13760822077



kaylayi0612@foxmail.com



oxygenconcentratormedical.com

Room D, 7th Floor, Guang Long Building, 162 Ping Xin North Road, PingHu Street, LongGang District, ShenZhen