

Blood Pressure Digital Monitor With LED display easy handling

Our Product Introduction

for more products please visit us on oxygenconcentratormedical.com

Basic Information

- Place of Origin: China
- Brand Name: HANYUE
- Certification: CE/FCC/ Rohs /FDA
- Model Number: HY-10
- Minimum Order Quantity: 60pcs
- Price: 5USD ~ 15USD
- Packaging Details: OPP packing, 30PCS /Carton
- Delivery Time: 5-8 Working days
- Payment Terms: D/A, D/P, T/T, Western Union, MoneyGram
- Supply Ability: 100000pcs



Product Specification

- Product Name: Blood Pressure Monitor
- Display: Digital Liquid Crystal Display
- Power Supply: 4 AAA Batteries 1.5V
- Memory Function: 99sets
- Accessories Included: Warranty Card, Certificate, Wrist Band, Product Manual, Storage Bag
- Application: Medical Healthcare, Adult, Families And Adult, Hospital, Arm /obses
- Highlight: **digital blood pressure machine, digital blood pressure monitors, adult blood pressure machine**

Product Description

How Accurate are Wrist Blood Pressure Monitors?

Wrist blood pressure monitors are a convenient option for people who want to monitor their blood pressure at home or on-the-go. However, the accuracy of wrist blood pressure monitors can be affected by several factors.

One of the main factors that can affect the accuracy of wrist blood pressure monitors is the positioning of the monitor on the wrist. The monitor should be worn at heart level, which can be difficult to achieve with a wrist monitor. If the monitor is not positioned correctly, it can provide inaccurate readings.

Another factor that can affect the accuracy of wrist blood pressure monitors is user error. It's important to follow the manufacturer's instructions carefully and to make sure that the cuff is wrapped snugly around the wrist. Additionally, wrist blood pressure monitors are sensitive to movement and may provide inaccurate readings if the user moves or talks during the measurement.

Research studies have shown that wrist blood pressure monitors can provide accurate readings when used correctly. However, it is important to note that the accuracy of wrist blood pressure monitors can vary between different models and brands. It's important to choose a high-quality monitor from a reputable manufacturer and to use it correctly to ensure accurate readings.

If you're concerned about the accuracy of your wrist blood pressure monitor, you may want to consider having your blood pressure checked by a healthcare professional to compare the readings. Additionally, it's important to follow up with your healthcare provider regularly to monitor your blood pressure and to discuss any concerns or questions you may have.



*Why And When Do You Need To Choose An Upper Arm Blood Pressure Monitor ?

1, For the Elder people whose vascular wall has poor flexibility. In this case, wrist blood pressure's result may be lower than the real. Because the wrist is the nerve endings, hypertensive patients may be complicated by arterial stenosis, blood clots and other vascular diseases, so the elderly use an Arm Blood Pressure Monitor that is more accurate than Wrist Blood Pressure Monitor!

2, Wrist blood pressure monitor is also accurate, but it is more suitable for young people, working people, travelling, gift. Because it is easy to carry.

What's in the box

- 1 x Electronic Blood Pressure Monitor
- 1 x Arm Band
- 1 x Storage Bag



Shenzhen Emeng Health Technology Co., LTD



+86 13760822077



kaylayi0612@foxmail.com



oxygenconcentratormedical.com

Room D, 7th Floor, Guang Long Building, 162 Ping Xin North Road, PingHu Street, LongGang District, ShenZhen