

Wireless Medical Pulse Oximeter For Patient Care Measure PR SpO2

Basic Information

- Place of Origin:
- Brand Name:
- Certification: CE/CFDA/RoHs/FCC/ISO13485/FDA

China

ACCARE

FS-10E

- Model Number:
- Minimum Order Quantity: 200PCS
- Price: 8~12USD
- Packaging Details:

🔊 accurate



Product Specification

- Type:
- Warranty:
- Safety Standard:
- Material:
- Instrument Classification:
 - Measure PR, S
- Function:Highlight:
- Class II Measure PR, SpO2 Wireless Medical Pulse Oximeter, Patient Care Medical Pulse Oximeter,

SpO2 Wireless Pulse Oximeter

Fingertip Pulse Oximeter

EN ISO13485.2016

12months

Plastic

OPP packing , 200pcs Per Carton

Our Product Introduction

Medical Pulse Oximeter Industry Adopts Wireless Technology For Enhanced Monitoring And Patient Care

Medical breast pump donation programs have been helping mothers in need gain access to breast pumps for some time now. These programs provide an essential service, as many mothers may not be able to afford a breast pump, or may not have access to one due to their geographical location or other circumstances.

Recently, medical breast pump donation programs have expanded their reach to support mothers worldwide. These programs are being implemented in low-income countries, where access to breast pumps is limited, and where the benefits of breast milk are particularly important for infant health.

International organizations and non-profits are leading the charge in expanding medical breast pump donation programs to support mothers worldwide. These programs typically involve the donation of new or gently used breast pumps, which are then distributed to mothers in need through local healthcare providers or community organizations.

In addition to providing breast pumps, these programs often include education and support for nursing mothers, to help them maximize the benefits of breastfeeding and ensure safe and effective milk expression.

The expansion of medical breast pump donation programs to support mothers worldwide is a positive development, as it helps to promote breastfeeding and improve infant health outcomes in low-income countries. By providing access to breast pumps, these programs can help reduce the risk of infections and malnutrition among infants, and improve long-term health outcomes for both mothers and babies.

Operating Guide:

Batteries must be inserted in the proper direction and then replaced by the cover.

The clip must be open as visualized below. The finger must be placed in a proper position.

The patient's finger must be placed aligned with the clip's rubber cushions in the right position and then clip the finger.

The switch on the front panel should be pressed.

During the process, the patient should be kept stationary and avoid the finger from shaking. Information is presented directly on the screen display.



Automatic Shut-down

Package Contents:

1*Pulse Oximeter

1*On	e hanging rope
1*Use	or Manual
1*Sto	rage bag
	l home care products (not including accessories) you purchase from Homecarewholesale.com, get a worry-free 1-year warranty iendly 24 hours customer service. Enjoy your shopping freely at Homecarewholesale!
	Shenzhen Emeng Health Technology Co., LTD
	🔇 +86 13760822077 🔄 kaylayi0612@foxmail.com 📀 oxygenconcentratormedical.com
	Room D, 7th Floor, Guang Long Building, 162 Ping Xin North Road, PingHu Street, LongGang District, ShenZhen