

## Babies Fetal Heart Monitor Machine Safe Heart Monitor For Pregnant Women

### Our Product Introduction

for more products please visit us on [oxygenconcentratormedical.com](http://oxygenconcentratormedical.com)

#### Basic Information

- Place of Origin: China
- Brand Name: HANYUE
- Certification: CE/FDA/ISO13485/FCC/RoHs
- Model Number: JSL-T501
- Minimum Order Quantity: 60pcs
- Price: 7USD
- Packaging Details: Each Boxes OPP packing /12.6 KGS / Per carton: 53\*29\*37cm / 60PCS Per Carton.
- Delivery Time: 3-5 Working days
- Payment Terms: T/T, Western Union, MoneyGram
- Supply Ability: 100000PCS



#### Product Specification

- Product Name: Doppler Fetal Monitor
- Type: Doppler Ultrasound Equipment
- Color: White , Pink,Blue
- Function: Baby Heart Rate Monitor
- Working Frequency: 2.0 MHz  $\pm$ 10%
- Audio Output Power: Er:  $\leq$ 1.0 W
- Highlight: **Babies Fetal Heart Monitor Machine,  
Pregnant Fetal Heart Monitor Machine,  
Safe Heart Monitor For Pregnant Women**

## Product Description

### Fetal Heart Detector Is Safe And Reliable Way To Monitor Babies Health During Pregnanc

Fetal heart detector technology, also known as fetal doppler or fetal monitor, is a safe and reliable way to monitor the health of a developing fetus during pregnancy. It uses ultrasound technology to detect and monitor the fetal heartbeat, providing valuable information about the health and well-being of the fetus. This technology has become an essential tool in prenatal care, allowing healthcare providers to identify and address potential problems early in pregnancy, ensuring the best possible outcomes for mother and baby.

Fetal heart detector technology is a safe and non-invasive way to monitor the health of a developing fetus during pregnancy. It uses ultrasound technology to detect and monitor the fetal heartbeat, providing valuable information about the health and well-being of the fetus. Fetal heart detector technology is typically used during routine prenatal appointments to check the fetal heart rate and rhythm. This information can help healthcare providers identify potential problems, such as fetal distress or growth issues, and intervene early if necessary to provide appropriate care.

The device consists of a handheld probe that is placed on the mother's abdomen, and it emits high-frequency sound waves that bounce off the fetal heart and create a sound that can be heard through the device's speakers. The fetal heart rate can be monitored throughout pregnancy, and changes in the heart rate can indicate potential problems such as fetal distress or growth issues. If the fetal heart rate is abnormal, healthcare providers may intervene with appropriate treatment or delivery.

Fetal heart detector technology is considered safe for both the mother and the fetus. It does not use any radiation and has no known harmful effects. The procedure is usually done by a trained sonographer or obstetrician and takes about 20-30 minutes.

In addition to monitoring the fetal heart rate, fetal heart detector technology can also be used to listen to the fetal heartbeat for the purpose of bonding with the baby. Many expectant parents find it comforting to hear their baby's heartbeat and feel more connected to their pregnancy. Overall, fetal heart detector technology is a safe and reliable way to monitor the health of a developing fetus during pregnancy. It provides valuable information that can help ensure a healthy pregnancy and delivery, making it an essential tool in prenatal care.

### Why might I need fetal heart monitoring?

Fetal heart rate monitoring is especially helpful if you have a high-risk pregnancy. Your pregnancy is high risk if you have diabetes or high blood pressure. It is also high risk if your baby is not developing or growing as it should.

Fetal heart rate monitoring may be used to check how preterm labor medicines are affecting your baby. These are medicines are used to help keep labor from starting too early.

Fetal heart rate monitoring may be used in other tests, including:

**Nonstress test.** This measures the fetal heart rate as your baby moves.

**Contraction stress test.** This measures fetal heart rate along with uterine contractions. Contractions are started with medicine or other methods.

**A biophysical profile (BPP).** This test combines a nonstress test with ultrasound.

Things that may affect the fetal heart rate during labor:

Uterine contractions

Pain medicines or anesthesia given to you during labor

Tests done during labor

Pushing during the second stage of labor

Your healthcare provider may have other reasons to use fetal heart rate monitoring.

## Nine weeks

Fetal heart sounds  
can be measured as  
early as 9 weeks

9

So you can hear  
Darling' s voice  
sooner



Shenzhen Emeng Health Technology Co., LTD

+86 13760822077

kaylayi0612@foxmail.com

oxygenconcentratormedical.com

Room D, 7th Floor, Guang Long Building, 162 Ping Xin North Road, PingHu Street, LongGang District,  
ShenZhen