



Continuous Medical Oxygen Concentrator High Concentration 10L Oxygen Generator

Our Product Introduction

for more products please visit us on oxygenconcentratormedical.com

Basic Information

- Place of Origin: China
- Brand Name: SS
- Certification: CE, ISO13485, SGS, FCC, RoHs
- Model Number: Y-X01W
- Minimum Order Quantity: Negotiable
- Price: Negotiable
- Packaging Details: OPP packing, 1pcs per Carton. Carton size: 540*370*695mm, Gross Weight: 21kgs.
- Delivery Time: 7-10 Working days
- Payment Terms: T/T, Western Union
- Supply Ability: 1000pcs



Product Specification

- Color: Black
- Oxygen Concentration: 93% \pm 3% 1L/min
- Product Name: Oxygen Concentrator For Medical & Home Use
- Oxygen Flow: 10L/min Adjustable
- Safety Standard: GB/T 18830-2009
- Power Source: Electric
- Display Control: LCD Screen
- Highlight: Continuous Medical Oxygen Concentrator, Medical Oxygen Concentrator 10L, High Concentration 10L Oxygen Generator



More Images



Product Description

Efficient and high concentration 10L continuous oxygen generator suitable for the elderly

A medical oxygen concentrator is a medical device that concentrates and purifies oxygen from the surrounding air, delivering it to individuals who require supplemental oxygen as part of their medical treatment. It's designed to be a safe and effective solution for oxygen therapy at home.

Our Product Introduction

Y-X01W

10L FLOW RATE ADJUSTABLE

Oxygen concentration 93% ± 3%

Oxygen concentration

Atomization

Smaller size
Lighter weight

Continuous and stable oxygen supply for 48 HOURS

How does an oxygen generator generate oxygen?

Oil-free pure copper compressor

Independent R&D, multiple leading technologies
Built-in overheat protection function; low noise and strong power

Independent R&D, High adsorption molecular sieve

Super adsorption capacity, slow decay

Super adsorption capacity

Pure oxygen outlet

Stable

Long use life

8 level filtering system

4 level filtering for air inlet and denitrogenation + 4 level filtering for oxygen outlet
Make sure the oxygen out let is pure and clean

Air intake: The POC draws in ambient air from the surroundings using a built-in compressor. This air contains approximately 21% oxygen, along with other gases like nitrogen, carbon dioxide, and trace elements.

Filtration: The incoming air passes through a series of filters to remove impurities, dust, and other particulate matter. These filters ensure that the air entering the concentrator is clean and free from contaminants that could affect the user's health.

Compression: The filtered air is then compressed using a compressor. The compressor increases the pressure of the air, allowing it to be more efficiently processed in subsequent stages.

Sieve bed adsorption: The compressed air is directed into a molecular sieve bed, which is typically filled with a material called zeolite. Zeolite has the ability to selectively adsorb nitrogen from the air while allowing oxygen to pass through. As a result, the nitrogen is trapped within the sieve bed, and the oxygen is concentrated.

Oxygen collection: The concentrated oxygen is collected and directed into a reservoir or storage chamber. This reservoir acts as a buffer, ensuring a continuous and stable supply of oxygen even when the user inhales rapidly or the demand fluctuates.

Oxygen delivery: The concentrated oxygen is delivered to the user through a nasal cannula or a mask. The user can breathe in the enriched oxygen, which helps to increase the oxygen levels in their bloodstream.

Waste gas release: The nitrogen and other waste gases that were adsorbed by the sieve bed during the adsorption process are released back into the environment. This allows the concentrator to continue functioning and producing concentrated oxygen.

Some of the benefits of using a medical oxygen concentrator for oxygen therapy at home include:

Convenience: Medical oxygen concentrators are designed to be used at home and provide a continuous supply of oxygen,

eliminating the need for regular trips to a healthcare facility.

Efficiency: Medical oxygen concentrators are an efficient way to provide oxygen therapy since they do not require the use of bulky and heavy oxygen tanks.

Cost-effective: Medical oxygen concentrators can be more cost-effective than other forms of oxygen therapy since they do not require frequent refilling or replacement of oxygen tanks.

Improved quality of life: Supplemental oxygen therapy can improve quality of life for individuals with respiratory conditions by reducing symptoms such as shortness of breath, fatigue, and difficulty sleeping.

Safety: Medical oxygen concentrators are designed to be safe and reliable, with built-in safety features such as alarms and automatic shut-off in case of a malfunction.

It's important to note that a medical oxygen concentrator should only be used under the guidance of a healthcare professional. They can help determine if oxygen therapy is appropriate for an individual's specific medical condition and provide guidance on how to properly use the device.

Overall, a medical oxygen concentrator can provide a safe, convenient, and effective solution for individuals who require supplemental oxygen therapy at home.

Home-use Health Care And Medical Treatment



HYPOXIC PEOPLE

Quickly increase blood oxygen saturation, increase oxygen content in the body, and help tp do djuvant therapy.



NORMAL HEALTH-CARE

When pregnant women, office workers or students are dizzy and fatigued, they can lessen their symptoms by inhaling oxygen. And it also helps to maintain their energy.



DOUBLE OXYGEN INHALATION

Medical level oxygen rate which concentration and flow can support two people to inhale oxygen at the same time.

Product Parameter

...



Product name MEDRIS 10L Oxygen Concentrator	Product model Y-X01W	Flow 1-10L/min adjustable flow rate
Oxygen concentration 93%±3% (10L/min)	Oxygen pressure range 86kPa-106kPa	Method Pressure swing adsorption (PSA)
Atomization rate ≥0.2ml/min (Only attached within Y-509W)	Operating noise ≤42dB(A weighting)	Operation mode Continuous operation
Power supply AC 220V/50Hz 110V/60Hz	Dimensions 375*340*650mm (length, width and height)	Net weight About 27.5KG


HIGH OXYGEN
CONCENTRATION


1-10L/MIN ADJUSTABLE
FLOW RATE


OXYGEN CONCENTRATION
REAL-TIME MONITORING


SMART ALERT


HD LARGE SCREEN


ATOMIZATION
FUNCTION


6-LEVEL LOW
NOISE TECHNOLOGY


MOLECULAR SIEVE


TWO-CYLINDER
PUMP COMPRESSOR

Wire-less Remote Control

Portable and steady in your hand



Packaging

01



02



03



 **Shenzhen Emeng
Health Technology Co.,
LTD**

+86 13760822077

Room D, 7th
Floor, Guang
Long Building,
162 Ping Xin
North Road,
PingHu Street,
LongGang
District,
ShenZhen

kaylayi0612@foxmail.com

oxyge