


Therapeutic Hyperbaric Pressure Chamber 220*150*190CM Cabin Size for Recovery

Our Product Introduction

ore products please visit us on oxygenconcentratormedical.com

Basic Information

- Place of Origin: SHENZHEN EMENG HEALTH
- Brand Name: O2EMENG
- Certification: CE/FDA/FC/ROHS/UKCA/SGS
- Model Number: EMY03
- Minimum Order Quantity: 1PC
- Price: \$23499
- Packaging Details: wooden packing
- Delivery Time: 15-20 work days
- Payment Terms: D/P, T/T, Western Union, MoneyGram
- Supply Ability: 100PCS




Product Specification

- Cabin Weight: 900KGS
- Properties: Rehabilitation Therapy Supplies
- Function: Hyperbaric Oxygen Therapy
- Chamber Oxygen Concentration: 28% 32%
- Usage: Medical/Therapeutic
- Cabin Size: 220*150*190CM
- Noise Level: Less Than 63 DB
- Rated Power: 1500W
- Highlight: therapeutic hyperbaric pressure chamber, therapeutic high pressure hyperbaric chamber, 1500w hyperbaric pressure chamber



More Images



Product Description

O2EMENG Therapeutic Hyperbaric Pressure Chamber 220*150*190CM Cabin Size for Recovery

Product Description:

The Hyperbaric Oxygen Chamber model EMY03 has a rated power of 1500W, which makes it a powerful tool for medical treatment. This hyperbaric chamber of oxygen is designed to allow patients to breathe in pure oxygen while inside a pressurized chamber. The oxygen-rich environment stimulates the body's natural healing processes, promoting the growth of new tissue and reducing inflammation.

The Hyperbaric Oxygen Chamber functions as a decompression chamber for oxygen therapy. As the pressure inside the chamber is slowly increased, the body's tissues absorb more oxygen. This promotes faster healing and recovery from a variety of medical conditions. The Hyperbaric Oxygen Chamber is commonly used for the following medical conditions:

Chronic wounds and ulcers
Burns and skin grafts
Carbon monoxide poisoning
Multiple sclerosis
Autism
Stroke
Sports injuries
Post-surgery recovery

The Hyperbaric Oxygen Chamber is a safe and non-invasive therapy that can enhance the healing process for a variety of medical conditions. This product is easy to use and can be operated by medical professionals or trained therapists. The Hyperbaric Oxygen Chamber is a highly recommended product for those looking for an effective and efficient therapy for their patients.

Features:

Product Name: Hyperbaric Oxygen Chamber

Properties: Rehabilitation Therapy Supplies

Function: Hyperbaric Oxygen Therapy

Cabin Weight: 900KGS

Noise Level: Less Than 63 DB

Chamber Oxygen Concentration: 28% 32%

Keywords: Hyperbaric pressure chamber, Hyperbaric oxygen therapy unit, Pressurized oxygen chamber

Technical Parameters:

Model	EMY03
Chamber Oxygen Concentration	28% 32%
Properties	Rehabilitation Therapy Supplies
Usage	Medical/Therapeutic
Cabin Weight	900KGS
Noise Level	Less Than 63 DB
Rated power	1500W
Cabin size	220*150*190CM
Function	Hyperbaric Oxygen Therapy

Applications:

The O2EMENG EMY03 Hyperbaric Oxygen Chamber is designed to offer a comfortable and spacious cabin size of 220*150*190CM, making it suitable for a wide range of patients. The cabin is also equipped with a range of advanced features, including a decompression chamber that helps to reduce the risk of decompression sickness during oxygen therapy.

The O2EMENG EMY03 Hyperbaric Oxygen Chamber is a powerful and reliable solution that can help you to achieve the best possible results from your oxygen therapy. With a rated power of 1500W, this high-pressure oxygen chamber is capable of delivering consistent and reliable results. It is also designed to be durable and long-lasting, with a cabin weight of 900KGS.

So whether you are looking for a high-pressure oxygen chamber for medical or therapeutic purposes, the O2EMENG EMY03

Hyperbaric Oxygen Chamber is the perfect choice. With its advanced features and reliable performance, it is sure to meet all your needs and exceed your expectations.

Support and Services:

The Hyperbaric Oxygen Chamber product comes with comprehensive technical support and services to ensure that customers have a seamless experience with our product. Our technical support team is available to assist with any queries or issues that may arise during installation, setup, or usage of the product. We also provide regular maintenance and repair services to keep the product in optimal

condition and extend its lifespan. Our services include:

1.Demand customization: according to your use environment and expectations, recommend the optimal core products and life aesthetics composition scheme.

2.Exclusive design: the advanced technology team and the top design team jointly, from the size and color material matching, the actual program customization landing.

3. Craftsmanship Manufacturing: The core components and other materials of hyperbaric chamber series products are manufactured in aerospace material factories.

4. Lifetime maintenance: provide perfect professional after-sales service guarantee

We are committed to providing the highest level of support and services to our customers to ensure their satisfaction with our Hyperbaric Oxygen Chamber product.

Reasons to choose us:

1.Professionalism and authority:

As a manufacturer or brand of hyperbaric chamber, O2EMENG has rich industry experience and technology accumulation, and can provide professional and reliable hyperbaric oxygen therapy equipment.

Its products may have undergone stringent quality control and certification and comply with relevant industry standards to ensure safety and effectiveness in use.

2. Technological advancement:

O2EMENG hyperbaric chambers may adopt the latest technology and design concepts, such as air pressurization or pure oxygen pressurization technology, as well as advanced air supply and exhaust (oxygen) systems, air conditioning systems and control systems. These technologies can ensure that the human body can safely and effectively inhale a high concentration of oxygen in a high-pressure environment for therapeutic purposes.

3. User feedback and reputation:

When choosing hyperbaric chambers, user feedback and reputation are also important reference factors, and O2EMENG hyperbaric chambers may be widely praised for their good treatment effect and user experience.

In addition, O2EMENG Hyperbaric Oxygen Chamber adopts advanced technology and design, which can ensure the safety during the treatment process. At the same time, Hyperbaric Oxygen Therapy, as a non-invasive physical therapy method, has significant therapeutic effects. Compared with other treatment methods, hyperbaric oxygen therapy may have lower cost and higher cost-effectiveness. o2EMENG hyperbaric chamber may further reduce the cost of use and improve economic efficiency by optimizing the design and production processes. Meanwhile this oxygen chamber has the following functions:

(1)Enhance the speed of recovery:

Hyperbaric oxygen chambers can help patients recover quickly from fatigue, illness, etc. O2EMENG Hyperbaric Oxygen Chambers may accelerate the speed of recovery by stimulating aerobic metabolic processes in cells by increasing the level of oxygen in the blood and the distance over which the blood oxygen diffuses.

(2)Immune Enhancement:

Hyperbaric Oxygen Therapy can help the elderly and the physically weak to increase blood oxygen levels and strengthen the body

O2EMENG

HYPERBARIC OXYGEN CHAMBER

- IMPROVE BRAIN BLOOD CIRCULATION
- SPEED UP METABOLISM
- HIGH PRESSURE OXYGEN INHALATION



HARD TYPE HYPERBARIC OXYGEN CHAMBER



Innovative technology, quality assurance

The diagram illustrates a hyperbaric oxygen chamber, a large, dark-colored, rectangular unit with a prominent door on the left side. To the right of the main unit is a smaller, black control cabinet with a digital display and several buttons. A blue curved line connects three callout boxes to the chamber. The first callout, at the top, points to the door and contains an icon of a door and the text 'The cabin door is made of aviation PC material'. The second callout, in the middle, points to the side of the chamber and contains a shield icon and the text 'Multiple safety technologies' followed by 'One-button pressure relief protection system during power outage'. The third callout, at the bottom, points to the control cabinet and contains a robot head icon and the text 'AI intelligent voice control system'. The brand name 'O2EMENG' is visible in blue letters on the front panel of the main chamber unit.

- The cabin door is made of aviation PC material
- Multiple safety technologies
One-button pressure relief protection system during power outage
- AI intelligent voice control system

WHY ARE MORE AND MORE PEOPLE CHOOSING HYPERBARIC OXYGEN CHAMBERS?



In some parts of the body, the capillaries are so narrow that blood cells cannot pass through them, bound oxygen cannot reach them, and they can only rely on dissolved oxygen. The oxygen generator has little change in the dissolved oxygen in the blood, which is also the reason why the ordinary oxygen absorption effect in the plateau area is general.





comfortable lying
adjustment



internal and external
intercom



personalized
multi-color matching



multimedia audio and
video functions



anti-claustrophobia
wide-open view



thoughtful details
anti-pinch hatch

PRODUCT PARAMETER

Product Model	EMY03 Two people hard type hyperbaric oxygen chamber
All-in-one Machine weight	85KG
All-in-one Machine dimensions	40*80*87cm (Length, width and height)
Noise	≤63dB
Voltage	AC:220V(±10%) 50Hz
Rated power	1500W
Oxygen production flow	10L/Min
Oxygen Purity	≥95%±3%
Chamber Capsule Color	White/Grey (Accept Customization)
Total chamber weight	900KG
Chamber size (Length, width and height)	220*150*190cm
Disposition	Chamber,All-in-one Machine,Cooling System,Smart massage seat,Internal and external dual screen control,AI intelligent voice system,Medical grade filter,Negative ion disinfection system,Air dehumidification,Accessories

CONCLUSION



It is more feasible to increase the dissolved oxygen in the blood.
The hyperbaric oxygen chamber can maintain the oxygen partial pressure within a safe and effective range by flexibly adjusting the oxygen concentration or pressure, and gently and effectively improve the oxygen delivery to tissues and organs by increasing dissolved oxygen.



INTEGRATION OF SUPERCHARGING AND OXYGENATION




Overturn the complexity of traditional household





New hyperbaric oxygen equipment



Oxygen equipment
+
air power system
+
Air Cooler

Intelligent all in one machine, no need to operate three devices, easy one-button operation saves time and effort.



Hd LCD screen

AI intelligent voice control system, internal and external dualscreen control

Low tone and low noise

Multiple noise reduction technology is used to solve the noise problem from the two aspects of structure and material, and the noise is controlled within 63 decibels.



Extra large viewable window

Observe the cabin personnel at any time, full light to eliminate tension

Water air conditioner

The cooling and ventilation unit integrates the cooling and ventilation unit with stable airflow and comfortable body, long-lasting and durable, safe and non-toxic



If you want to know more about our products, please contact us!



O 2 E M E N G

After Sale Service

WhatsApp: +86 13760822077

E-mail: kaylayi0612@foxmail.com

Hotline: +86 13760822077

Official website: emoxygen.com



Shenzhen Emeng Health Technology Co., LTD



+86 13760822077



kaylayi0612@foxmail.com



oxygenconcentratormedical.com

Room D, 7th Floor, Guang Long Building, 162 Ping Xin North Road, PingHu Street, LongGang District,
ShenZhen