

Shenzhen Emeng Health Technology Co., LTD oxygenconcentratormedical.com

O2EMENG Advanced Hyperbaric Chamber for Rehabilitation - EMY01 (28% - 32% Oxygen Concentration)

Basic Information

Place of Origin: SHENZHEN EMENG HEALTH

• Brand Name: O2EMENG

• Certification: CE/FDA/FC/ROHS/UKCA/SGS

Model Number: EMY01Minimum Order 1PC

Quantity:

• Price: \$23999

Packaging Details: wooden packingDelivery Time: 15-20 work days

• Payment Terms: D/P, T/T, Western Union, MoneyGram

• Supply Ability: 100PCS



Product Specification

• Properties: Rehabilitation Therapy Supplies

Usage: Medical/Therapeutic

 Chamber Oxygen 28% 32% Concentration:

• Cabin Szie: 185*88*152CM

Model: EMY01

Noise Level: Less Than 63 DB

• Function: Hyperbaric Oxygen Therapy

• Cabin Weight: 1000KGS

• Highlight: medical hyperbaric oxygen chamber,

medical oxygen hyperbaric chamber, therapeutic hyperbaric oxygen chamber



More Images







Product Description

O2EMENG Advanced Hyperbaric Chamber for Rehabilitation - EMY01 (28% - 32% Oxygen Concentration)

Product Description:

This decompression chamber for oxygen therapy is a heavy-duty product, with a cabin weight of 1000KGS. The chamber is spacious and comfortable, providing ample room for patients to relax and receive treatment. The Hyperbaric Oxygen Chamber is easy to operate and comes equipped with a variety of safety features, including pressure sensors, emergency release valves, and alarms.

The Hyperbaric Oxygen Chamber is designed to deliver oxygen at a concentration of 28% 32%, which is higher than the normal atmospheric concentration of 21%. This increased oxygen concentration promotes healing and helps to reduce inflammation in the body.

The chamber is also equipped with a humidifier to provide a comfortable environment for patients during treatment.

The Hyperbaric Oxygen Chamber is rated at 1500W and is designed to operate efficiently and effectively. The chamber is easy to maintain and comes with a variety of accessories to ensure that patients receive the best possible treatment. The Hyperbaric Oxygen Chamber is a powerful tool for medical professionals and patients alike, providing a safe and effective way to deliver oxygen therapy.

Features:

Product Name: Hyperbaric Oxygen Chamber **Properties:** Rehabilitation Therapy Supplies

Usage: Medical/Therapeutic

Function: Hyperbaric Oxygen Therapy

Cabin size: 185*88*152CM Rated power: 1500W

Keywords: High-pressure oxygen chamber, Decompression chamber for oxygen therapy, HBOT chamber

Technical Parameters:

Product Attribute	Value
Model	EMY01
Function	Hyperbaric Oxygen Therapy
Properties	Rehabilitation Therapy Supplies
Usage	Medical/Therapeutic
Chamber Oxygen Concentration	28% 32%
Noise Level	Less Than 63 DB
Cabin Weight	1000KGS
Cabin Size	185*88*152CM
Rated Power	1500W

Applications:

The EMY01 model is a spacious chamber with dimensions of 185*88*152CM, providing ample space for patients to receive treatment comfortably. The cabin weight is 1000KGS, making it a sturdy and durable unit that can withstand frequent use.

The chamber oxygen concentration is set at 28% 32%, which is the optimal range for providing effective hyperbaric oxygen therapy. This therapy is used to treat a variety of medical conditions, including wounds that won't heal, carbon monoxide poisoning, decompression sickness, and more.

The EMY01 model is an ideal choice for rehabilitation therapy supplies and can be used in a variety of settings, including hospitals, clinics, and rehabilitation centers. It is a versatile piece of equipment that can be used for a wide range of applications, including treating sports injuries, post-surgical recovery, and more.

Some common scenarios in which the EMY01 hyperbaric oxygen therapy unit might be used include:

Treating patients with chronic wounds that are difficult to heal using other methods

Helping athletes recover from sports-related injuries

Providing post-surgical care to patients to speed up the healing process

Treating patients with carbon monoxide poisoning or decompression sickness

If you are looking for a reliable and effective compression chamber for hyperbaric oxygen therapy, the O2EMENG EMY01 model is an excellent choice. It is produced by a reputable manufacturer and has been rigorously tested to ensure its effectiveness and safety.

Support and Services:

Our Hyperbaric Oxygen Chamber product comes with comprehensive technical support and services to ensure that you can operate and maintain it with ease. Our team of experts is available to provide you with any assistance you may need, whether it be installation, troubleshooting, or repair.

In addition, we offer some survices, for example:

- **1.Demand customization:** according to your use environment and expectations, recommend the optimal core products and life aesthetics composition scheme.
- **2.Exclusive design:** the advanced technology team and the top design team jointly, from the size and color material matching, the actual program customization landing.
- 3. Craftsmanship Manufacturing: The core components and other materials of hyperbaric chamber series products are manufactured in aerospace material factories.
 - 4. Lifetime maintenance: provide perfect professional after-sales service guarantee

Other details:

- 1.Smart Chair: Ergonomic design with USB-A charging port allows you to easily charge while enjoying hyperbaric oxygen therapy and keep your device fully charged. The frame of the seat is made of 304 stainless steel, which is durable and highlights the quality and style.
- 2. High strength explosion-proof material door: The use of specially developed high strength, high permeability explosion-proof material, combined with excellent sealing performance, to ensure the stability and safety of the cabin environment. The door body design takes into account pressure stability and fast opening function, with Al intelligent voice system, easy to operate and safe.
- 3. Comprehensive data monitoring: Built-in high-precision sensor, real-time monitoring and display of cabin temperature, humidity, pressure and oxygen concentration, so that you know every oxygen therapy process. Accurate data to provide scientific basis for your health management.
- **4. Quiet design:** Pay special attention to noise control, ensure that the cabin environment is quiet and comfortable, the sound level is 63dB, and create a private space away from the noise and relax.
- **5. Stable power supply:** Support 220V/50Hz standard voltage, rated power up to 1500W, to ensure stable operation of equipment, to meet the needs of long-term hyperbaric oxygen therapy.
- 6. Precise oxygen concentration control: The oxygen concentration in the cabin can be accurately adjusted to between 28% and 32%, and the scientific ratio can effectively improve blood oxygen saturation, promote cell metabolism, and accelerate body recovery and repair.
- 7. Comfortable temperature control: Equipped with advanced cooling system, you can freely adjust the cabin temperature according to personal needs, to ensure that you can enjoy the most suitable oxygen therapy environment in the hot summer.

[Applicable audience] :

Accelerate physical recovery and improve performance for athletes and fitness enthusiasts.

Sub-healthy people, improve sleep quality, enhance immunity


Brain workers, relieve brain fatigue, improve work efficiency.

Chronic disease patients, assist rehabilitation treatment, promote health recovery.

Choose our high-end intelligent hyperbaric oxygen chamber, open a new chapter of your health, let every breath become the nourishment of life, enjoy the health and happiness brought by science and technology

POSITIVE EFFECTS OF HYPERBARIC **OXYGEN THERAPY ON THE HUMAN BODY**

It can increase the oxygen concentration by increasing the ambient air pressure, thus improving the efficiency of the body's utilization of oxygen.



Literature sources: Biddlestone J, Bandarra D, Rocha S. The role of hypoxia in inflammatory diseasereview)[J] Int JMol Med,2015,35(4):p.859-69.

The above knowledge science for consumers to learn reference, does not represent the efficacy of the product the actual effect varies from person to person.

PROFESSIONAL REGISTRATION DEPARTMENT



O2 EMENG have an special professional registrationdepartment, which in charge of domestic registrationand oversea registration. We can help our current ofpotential customers to do the registration work. Wehave a good relation with domestic famous institutewho can provide us with authoritative medicaanalysis and authorization report. Quality and Service are always top concern in Olive. We providefull set of service, including date and samples, which are available on request









+00 13/000220// waylayloo12@i0xinali.com vxygenconcentratormedical.com Room D, 7th Floor, Guang Long Building, 162 Ping Xin North Road, PingHu Street, LongGang District,

ShenZhen