

Shenzhen Emeng Health Technology Co., LTD oxygenconcentratormedical.com

CE and ISO Certified Disposable Medical Diagnostic Test Kits for Antigen Rapid Test

Basic Information

Place of Origin: china
Brand Name: DVOT
Certification: CE/ISO
Model Number: 4 IN 1
Minimum Order 800test

Quantity:
• Price: 1 1.2U SD

• Packaging Details: 800T/PER CTN Weight :17.5 KGS CTN SIZE

:57*54*34CM

• Delivery Time: 15-20 working days

• Payment Terms: D/P, T/T, Western Union, MoneyGram

• Supply Ability: 1000000pcs



Product Specification

Certifications: CE ISO

Result Time:
 15 Minutes To Get

Components: Test Cassette, Dropper, Buffer Solution

Test Time: 15 MinutesTime Of Use: DisposableAccuracy: 99%

Type: Rapid Test KitKeyword: Antigen Rapid Test

Highlight: disposable medical diagnostic test kits,

ce clungene antigen rapid test, ce medical diagnostic test kits



More Images









Product Description

CE and ISO Certified Disposable Medical Diagnostic Test Kits for Antigen Rapid Test

Product Description:

DVOT (4-in-1) is including to testing SARS-CoV-2 & Influenza A/B & RSV Antigen Combo Test Kit.

Our Antigen Rapid Test Kit is a type of rapid test kit, meaning that it provides results quickly and efficiently. This is especially important in the case of COVID-19 testing, as time is of the essence in preventing the spread of the virus. With this home testing kit, you can get results in just minutes, allowing you to take the necessary precautions and seek medical attention if needed.

The Antigen Rapid Test Kit is a disposable product, making it easy and convenient to use. Each kit comes with a test cassette, dropper, and buffer solution, providing everything you need to conduct the test at home. The components are designed to be user-friendly, so you don't need any special training or expertise to use the kit.

In addition to its accuracy and convenience, the Antigen Rapid Test Kit is also affordable, making it accessible to more people. With this diagnostic test kit, you can easily monitor your health and keep yourself and your loved ones safe from COVID-19.

Overall, the Antigen Rapid Test Kit is an excellent home testing kit for COVID-19, providing accurate results quickly and affordably. With its disposable components and user-friendly design, it's easy to use and accessible to more people. If you're looking for a reliable way to test for COVID-19 at home, the Antigen Rapid Test Kit is an excellent choice.

Technical Parameters:

Product Attribute	Description
Product Type	Rapid Test Kit
Keyword	Antigen Rapid Test
Test Time	15 Minutes
Result Time	15 Minutes To Get
Accuracy	99%
Certifications	CE,ISO
Components	Test Cassette, Dropper, Buffer Solution
Time Of Use	Disposable

Applications:

The Antigen Rapid Test Kit is primarily used for the detection of SARS-CoV-2 antigen in human nasal swab specimens, which is an important diagnostic tool in identifying Covid-19 infections. It is suitable for use in hospitals, clinics, and laboratories, and is an effective way to quickly test large numbers of people for the virus.

The DVOT 4 IN 1 Antigen Rapid Test Kit is an accurate and reliable way to diagnose Covid 19 infections, and it can also be used in conjunction with the Antibody Diagnostic Kit to determine the stage of the disease. The Antibody Diagnostic Kit is designed to detect the presence of antibodies in the blood, which can help to determine whether a person has been infected with the virus in the past.

The Antigen Rapid Test Kit is ideal for use in a variety of different situations, including in airports, schools, and workplaces, where it can be used to screen large numbers of people quickly and efficiently. It is also suitable for use in emergency situations, where rapid diagnosis is essential.

The DVOT 4 IN 1 Antigen Rapid Test Kit is a valuable tool in the fight against Covid 19, and its accuracy rate of 99% makes it a reliable diagnostic tool for healthcare professionals. It is a cost-effective way to test for Covid 19, and its rapid results mean that patients can receive treatment quickly and efficiently.

In summary, the DVOT 4 IN 1 Antigen Rapid Test Kit is an essential tool for the diagnosis of Covid 19 infections. Its accuracy, speed, and ease of use make it an ideal choice for healthcare professionals, and its versatility means that it can be used in a variety of different settings. Whether in hospitals, clinics, or emergency situations, the Antigen Rapid Test Kit is a valuable diagnostic tool.

Support and Services:

Our Medical Diagnostic Test Kits product comes with comprehensive technical support and services to ensure accurate and reliable results. We also help you get started with the product quickly and efficiently.

Key benefits:

Qualitative detection of SARS-CoV-2, Influenza A/B and RSV in the nose

Cover your winter related respiratory viruses in only one test

Distinguish your positive result(s) at one glance High accuracy, cost-effective and easy to use

COVID: reliable with Alpha, Delta & Omicron variants, BA4, BA5, CH.1.1, and EG.5.

1. Reliable test performance

2.Fast answers in just 15 minutes

3. Supports differential detection of SARS-CoV-2

4.Simple, 4-step process

5. Fully integrated assay, requiring only one application area and one sample

Steps to use:

Step 1: Insert a sterile swab into the patient's nostril. Rotate the swab 5 times over the nasopharyngeal surface and then leave the swab in place for 15 seconds to absorb secretions. Remove the swab from the nasal cavity.

Step 2: Immerse the swab in the sample container, making sure the swab is in full contact with the fluid, but do not apply too much force to avoid splashing. Cover and gently shake the sample container.

Step 3: Place the test instrument on a flat surface and place 4 drops of the immersion sample into the sample hole of the test instrument at a 90° angle.

Step 4: Read test results within 15 to 30 minutes. Do not touch or move the test apparatus until results can be read.









Shenzhen Emeng Health Technology Co., LTD



+86 13760822077



kaylayi0612@foxmail.com



oxygenconcentratormedical.com

 $Room\ D,\ 7th\ Floor,\ Guang\ Long\ Building,\ 162\ Ping\ Xin\ North\ Road,\ PingHu\ Street,\ LongGang\ District,$ ShenZhen