

## Highly Airtight Oxygen health Hyperbaric Chamber 1500W For Medical / Therapeutic

## **Basic Information**

- Place of Origin:
- Brand Name:
- Certification:
- Model Number:
- Minimum Order Quantity:
- Price:
- Packaging Details:
- Delivery Time:

- SHENZHEN EMENG HEALTH **O2EMENG** 
  - CE/FDA/FC/ROHS/UKCA/SGS

- Payment Terms:
- Supply Ability:
- EMY05 1PC \$11000 wooden packing
- 15-20 work days

100PCS

- D/P, T/T, Western Union, MoneyGram

0 2 E M E N G

## **Product Specification**

- Function:
- Tank Weight:
- Properties:
- Tank Szie:
- Chamber Oxygen • Concentration:
- Usage:
- Rated Power:
- Highlight:

Hyperbaric Oxygen Therapy 800KGS **Rehabilitation Therapy Supplies** 220\*120\*190CM

28% 32%

Medical/Therapeutic

1500W

1500w hyperbaric oxygen chamber, 1500w oxygen hyperbaric chamber, medical hyperbaric oxygen chamber



## More Images





## Highly Airtight Oxygen health Hyperbaric Chamber 1500W For Medical / Therapeutic

#### **Product Decrition:**

The Hyperbaric Chamber is a portable hyperbaric chamber that weighs 800KGS, making it easy to transport and store. The chamber is equipped with a powerful 1500W motor that delivers superior performance and reliability. With its advanced features and high-quality construction, this Hyperbaric Chamber is the perfect choice for medical and therapeutic applications.

The Hyperbaric Chamber has a tank size of 220\*120\*190CM, providing ample space for patients to receive their therapy. The chamber is designed to provide a comfortable and relaxing environment for patients while also delivering maximum therapeutic benefits. Whether you are looking to treat a specific medical condition or simply improve your overall health and wellbeing, the Hyperbaric Chamber is the perfect solution.

One of the key benefits of the Hyperbaric Chamber is its versatility. The portable hyperbaric chamber can be used for a wide range of medical and therapeutic applications, including wound healing, sports injuries, chronic fatigue syndrome, autism, and many other conditions. The chamber is also an effective tool for enhancing athletic performance and promoting overall wellness.

Another advantage of the Hyperbaric Chamber is its ease of use. The chamber is designed to be user-friendly, with simple controls and an intuitive interface. This makes it easy for patients to receive their therapy without any hassle or confusion.

### Advantage:

The O2EMMENG EMY05 model hyperbaric chamber has several significant differences and advantages over traditional hyperbaric chambers:

1. Advanced Technology and Control Systems: The EMY05 model utilizes the latest in hyperbaric technology with a precise oxygen concentration and pressure control system. Traditional hyperbaric chambers may use more basic controls that do not provide the same precise regulation and monitoring capabilities.

2. New Hyperbaric Oxygen Equipment: Intelligent all-in-one machine, unlike traditional equipment, does not require the operation of three machines, Oxygen Equipment + Aerodynamic Equipment + Air Cooler, realizing easy one-button operation and saving time and effort.

3. Comfort and User Experience: The EMY05 model is designed with a more spacious and comfortable interior space, taking into account the patient's comfort and convenience during treatment.

4. Durability and Reliability: The EMY05 model undergoes rigorous quality control and testing to ensure stable operation over long periods of time without frequent maintenance.

5. Technological innovations and safety: The EMY05 model integrates new technological innovations that enhance safety and therapeutic efficacy. It allows for more precise control of oxygen delivery to support medical treatments more effectively.

## Support and Services:

The Hyperbaric Oxygen Chamber product comes with comprehensive technical support and services to ensure that you get the most out of your investment. Our team of experienced technicians and engineers are available to provide support and assistance whenever you need it.

Our technical support services include:

**1.Demand customization:** according to your use environment and expectations, recommend the optimal core products and life aesthetics composition scheme.

**2.Exclusive design:** the advanced technology team and the top design team jointly, from the size and color material matching, the actual program customization landing.

3. Craftsmanship Manufacturing: The core components and other materials of hyperbaric chamber series products are manufactured in aerospace material factories.

4. Lifetime maintenance: provide perfect professional after-sales service guarantee

At Hyperbaric Oxygen Chamber, we are committed to providing the highest level of support and service to our customers. Contact us today to learn more about our technical support and services.

# WHY ARE MORE AND MORE PEOPLE CHOOSING HYPERBARIC OXYGEN CHAMBERS?

0

In some parts of the body, the capillaries are so narrow that blood cells cannot pass through them, bound oxygen cannot reach them, and they can only rely on dissolved oxygen. The oxygen generator has little change in the dissolved oxygen in the blood, which is also the reason why the ordinary oxygen absorption effect in the plateau area is general.

O2EMENG



internal and external intercom



comfortable lying adjusment



multimedia audio and video functions



anti-claustrophobia

wide-open view

. .





MEDICAL INSTITUTION





personalized multi-color matching

thoughtful details

anti-pinch hatch



FITNESS INDUSTRY





HOTEL



PLATEAU AREA

INDIVIDUAL FAMILY

WORK PLACE



#### 1. How does hyperbarics work?

Hyperbarics is a technology in which the air pressure in the environment is increased. When a person's body is placed in a more pressure environment, it absorbs more oxygen molecules per volume of compressed air.

The body normally transports oxygen via the hemoglobin of the red blood cells. By increasing the air pressure, oxygen is then driven into the body's fluids, allowing a super-saturation of tissues and organs with oxygen. The increased pressure infuses the body with oxygen, even reaching injuries with damaged circulation.

#### 2. What are the health benefits of hyperbaric oxygen chambers?

Studies have shown that hyperbaric oxygen chambers can treat hypoxic dieases, relieve fatigue, delay aging, prevent cardiovascular disease, improve sleep quality, accelerated wound repair and enhance immunity. We have customers serviced shows hyperbaric oxygen therapy makes very positive results.

#### 3. How long is the treatment?

Usually, we recommend our clients receive the HBOT therapy 5-6 sessions per week, 60-90mins per session. Using frequency or duration can be adjusted according to individual circumstances.

#### 4. Is HBOT safe?

Hyperbaric oxygen therapy is prescribed by a physician and performed under medical supervision. Although there are minor risks like all medical treatments, overall hyperbaric oxygen therapy is extremely safe.

#### 5. Is there anyone who can not use our chamber?

Pregnant ladies, patients who catch a cold or Flu, or those who have rhinitis or otitis media are not recommended to use the chamber.

#### 6. How do you feel inside a hyperbaric oxygen chamber?

When lying in the oxygen chamber, you will feel comfortable and at ease. Also, You may feel a little squeezing of the ears because of the pressure. It is quite normal and can be relieved by chewing or swallowing.

#### 7. After bought the chamber, do I need to by any other things?

No, the price of our chamber includes the chamber capsule, all machines, namely the air compressor, oxygen concentrator, air cooler, and accessories like frame, tubes, etc, so you don't have to buy anything else.

#### 8. Is it complex to install the chamber or change the filter?

Of course not, we will provide you with a video to show you how to install the chamber and change the filter.

We also support online services to help you with all the issues, so don't worry.

#### 9. What is the general time warranty expected for ZHIKANG oxygen chambers?

We will provide one-year warranty and lifelong after-sales services. If any quality problem/fault in material/design is under correct operation within one year, If it's easy to fix, we will send new components for free and teach you how to repair them. (Even after the one-year warranty period, you can also contact us freely and we will be glad to help you.)

# **OUR CERTIFICATE**





