

Hyperbaric Oxygen Chamber for Athletes' Recovery 1.3 ata (adjustable , Leading the New Trend of Oxygen Therapy

Our Product Introduction

for more products please visit us on oxygenconcentratormedical.com

Basic Information

- Place of Origin: SHENZHEN EMENG HEALTH
- Brand Name: O2EMENG
- Certification: CE/FDA/FC/ROHS/UKCA/SGS
- Model Number: EMY04
- Minimum Order Quantity: 1PC
- Price: \$9500
- Packaging Details: wooden packing
- Delivery Time: 15-20 work days
- Payment Terms: D/P, T/T, Western Union, MoneyGram
- Supply Ability: 100PCS




Product Specification

- Tank Size: 135*80*185CM
- Tank Weight: 450KGS
- Machine Size: 80*40*85CM
- Machine Weight: 85KGS
- Rated Power: 1500W
- Voltage: 220V, 50Hz
- Chamber Oxygen Concentration: 28% 32%
- Chamber Temperature: Cooling Adjustable
- Noise Level: Less Than 63 DB
- Highlight: 1500w hyperbaric oxygen chamber, 1500w oxygen hyperbaric chamber



More Images



Product Description

Hyperbaric Oxygen Chamber for Athletes' Recovery, Leading the New Trend of Oxygen Therapy and Wellness

Product Description:

The EMY04 model is O2EMENG's flagship product, integrating advanced technology designed to provide the best user experience. Its new hyperbaric oxygen equipment is an intelligent all-in-one machine that combines oxygen equipment, aerodynamic system and air cooler, easily operated by one button, saving time and labor. It is widely used in healthcare, home as well as beauty salons.

O2EMENG Hyperbaric Oxygen Chamber adopts advanced technology and design, which can ensure the safety during the treatment process. Meanwhile, hyperbaric oxygen therapy, as a non-invasive physical therapy method, has significant therapeutic effects. Compared with other treatment methods, hyperbaric oxygen therapy may have lower cost and higher cost-effectiveness. O2EMENG hyperbaric oxygen chamber may further reduce the cost of use and improve the economic benefits by optimizing the design and production process. At the same time, this oxygen chamber is suitable for people who are easy to be fatigued, work late at night and overtime, insomnia patients, ambiguous skincare people, and sub-healthy people, etc. It can increase the oxygen, blood oxygen saturation, and the diffusion distance of oxygen in the body.

All in all, our hyperbaric chamber is a top-quality product with excellent performance and reliability. Its compact size allows it to be installed in many medical facilities without taking up too much space. It is an essential piece of equipment for any hospital or clinic that provides high-quality rehabilitation therapy supplies. It produces at least 90% oxygen concentration and is one of the best hyperbaric chambers on the market for a great experience. Buy it today and experience the benefits of hyperbaric oxygen therapy.

Reasons to choose us:

1. Professionalism and authority:

As a manufacturer or brand of hyperbaric chamber, O2EMENG has rich industry experience and technology accumulation, and can provide professional and reliable hyperbaric oxygen therapy equipment. Its products may have undergone stringent quality control and certification and comply with relevant industry standards to ensure safety and effectiveness in use.

2. Technological advancement:

O2EMENG hyperbaric chambers may adopt the latest technology and design concepts, such as air pressurization or pure oxygen pressurization technology, as well as advanced air supply and exhaust (oxygen) systems, air conditioning systems and control systems.

These technologies can ensure that the human body can safely and effectively inhale a high concentration of oxygen in a high-pressure environment for therapeutic purposes.

3. User feedback and reputation:

When choosing hyperbaric chambers, user feedback and reputation are also important reference factors, and O2EMENG hyperbaric chambers may be widely praised for their good treatment effect and user experience.

Advantages:

The O2EMENG EMY04 model hyperbaric chamber has several significant differences and advantages over traditional hyperbaric chambers:

1. Advanced Technology and Control Systems: The EMY04 model utilizes the latest in hyperbaric technology with a precise oxygen concentration and pressure control system. Traditional hyperbaric chambers may use more basic controls that do not provide the same precise regulation and monitoring capabilities.

2. New Hyperbaric Oxygen Equipment: Intelligent all-in-one machine, unlike traditional equipment, does not require the operation of three machines, Oxygen Equipment + Aerodynamic Equipment + Air Cooler, realizing easy one-button operation and saving time and effort.

3. Comfort and User Experience: The EMY04 model is designed with a more spacious and comfortable interior space, taking into account the patient's comfort and convenience during treatment.

4. Durability and Reliability: The EMY04 model undergoes rigorous quality control and testing to ensure stable operation over long periods of time without frequent maintenance.

5. Technological innovations and safety: The EMY04 model integrates new technological innovations that

enhance safety and therapeutic efficacy. It allows for more precise control of oxygen delivery to support medical treatments more effectively.

Support and Services:

1.Demand customization: according to your use environment and expectations, recommend the optimal core products and life aesthetics composition scheme.

2.Exclusive design: the advanced technology team and the top design team jointly, from the size and color material matching, the actual program customization landing.

3. Craftsmanship Manufacturing: The core components and other materials of hyperbaric chamber series products are manufactured in aerospace material factories.

4. Lifetime maintenance: provide perfect professional after-sales service guarantee

Our Hyperbaric Oxygen Chamber product comes with comprehensive technical support and services to ensure optimal performance and customer satisfaction. Our team of experienced technicians and engineers are available to assist with installation, maintenance, and troubleshooting of the product.

In addition, we offer training sessions for customers to learn how to properly operate and maintain the Hyperbaric Oxygen Chamber. This includes training on safety protocols and best practices for usage.

Our technical support and services also include regular maintenance and repair services to keep the product functioning at its best. We offer both on-site and remote support options to better serve our customers.

If you have any questions or concerns about our technical support and services for the Hyperbaric Oxygen Chamber, please do not hesitate to contact us.

Supercharged Oxygen




Increase the partial pressure of oxygen in the body.
The hyperbaric chamber can be raised to 30kpa.




Increase blood oxygen saturation.
Oxygen intake is equivalent to more than 2 times that of normal breathing.



Increase the oxygen diffusion distance in the body.
Use the flow of lymph to increase blood flow.




The blood vessel circulation in the human body under normal atmospheric pressure. Various factors impede flow.




Hyperbaric oxygen chamber speeds up human blood vessel circulation, promote blood flow.

Increases the partial pressure of oxygen in the blood, and improves the vascular transport and terminal arrival of oxygen.



Atmospheric Oxygen

At standard atmospheric pressure, blood vessels contain oxygen (atmospheric oxygen)



Supercharged Oxygen

In a hyperbaric oxygen environment, the blood vessels contain oxygen (high concentrations of oxygen)

Israeli Scientists Reverse Aging In Human Trials For The First Time

Human Telomere Lengthening

20%

Number of experimenters

25

Decreased Senescent Cells

37%

Sample age


≥64

Cell Young

25 years old

Experiment period

90 Days



Adapt to the crowd



Body easily fatigued



Work out



Working overtime and staying up late



Insomnia patients



Student population



Beauty-loving skincare crowd



Socializing drunk



Middle-aged and elderly



Sub-healthy population

If you want to know more about our products, please contact us!



O2EMENG


After Sale Service

WhatsApp: +86 13760822077

E-mail: kaylayi0612@foxmail.com

Hotline: +86 13760822077

Official website: emoxygen.com



Shenzhen Emeng Health Technology Co., LTD



+86 13760822077



kaylayi0612@foxmail.com



oxygenconcentratormedical.com

Room D, 7th Floor, Guang Long Building, 162 Ping Xin North Road, PingHu Street, LongGang District,
ShenZhen