



Shenzhen Emeng Health Technology Co., LTD

oxygenconcentratormedical.com

O2EMENG Hard Type 1.3ATA Double Hard Seat Hyperbaric Camera Oxygen Capsules Wholesale HBOT for Sports Recovery

Our Product Introduction

For more products please visit us on oxygenconcentratormedical.com


Basic Information

- Place of Origin: SHENZHEN EMENG HEALTH
- Brand Name: O2EMENG
- Certification: CE/FDA/FC/ROHS/UKCA/SGS
- Model Number: EMY03
- Minimum Order Quantity: 1PC
- Price: \$23500
- Packaging Details: wooden packing
- Delivery Time: 15-20 work days
- Payment Terms: D/P, T/T, Western Union, MoneyGram
- Supply Ability: 100PCS



Product Specification

- Chamber Oxygen Concentration: 28% 32%
- Model: EMY03
- Properties: Rehabilitation Therapy Supplies
- Noise Level: Less Than 63 DB
- Cabin Weight: 900KGS
- Cabin Size: 220*150*190CM
- Function: Hyperbaric Oxygen Therapy
- Working Pressure: 1.3ATA adjustable
- Highlight: **Hyperbaric Camera Oxygen Capsules Wholesale**
,
Double Hard Seat Hyperbaric Oxygen Chamber,
Sports Recovery Hyperbaric Camera Oxygen Capsules



More Images



Product Description

O2EMENG 28%~32% Concentration EMY03 Hyperbaric Oxygen Chamber for Rehabilitation Therapy Supplies

Product Description:

With a chamber oxygen concentration of 28% 32%, this high-pressure oxygen chamber is capable of delivering the right amount of oxygen to your body to promote faster healing and recovery. Hyperbaric oxygen therapy has been proven to be effective in treating a wide range of conditions, including wounds, infections, and even cancer.

The EMY03 model is designed with your comfort in mind. It operates at a noise level of less than 63 DB, which means that you can relax and enjoy your therapy session in peace and quiet. The chamber is spacious and can accommodate individuals of different sizes with ease.

As a rehabilitation therapy supply, this high-pressure oxygen chamber is perfect for individuals who are recovering from an injury or illness. The benefits of hyperbaric oxygen therapy are numerous, from reducing inflammation and swelling to promoting the growth of new blood vessels.

Overall, the EMY03 Hyperbaric Oxygen Chamber is a valuable addition to any rehabilitation therapy program. It is safe, effective, and easy to use. Whether you are looking to improve your overall health and well-being or recover from an injury or illness, this high-pressure oxygen chamber is the perfect solution for you.

Technical Parameters:

Model	EMY03
Noise Level	Less Than 63 DB
Function	Hyperbaric Oxygen Therapy
Usage	Medical/Therapeutic
Chamber Oxygen Concentration	28% 32%
Cabin Weight	900KGS
Working pressure	1.3ATA(adjustable)
Cabin size	220*150*190CM
Properties	Rehabilitation Therapy Supplies

Applications:


Hyperbaric Oxygen Therapy (HBOT) is a well-known treatment for various medical conditions. The EMY03 model is specifically designed to provide a comfortable and effective therapy experience. It is a compression chamber for hyperbaric oxygen therapy that can safely and effectively deliver 100% pure oxygen to patients.

The O2EMENG EMY03 Hyperbaric Oxygen Chamber is suitable for a wide range of medical conditions, including but not limited to:

1. Non-healing wounds and ulcers
2. Burn injuries
3. Carbon monoxide poisoning
4. Crush injuries and compartment syndrome
5. Decompression sickness
6. Gas gangrene
7. Infections
8. Multiple sclerosis
9. Stroke
10. Autism
11. And many more

The O2EMENG EMY03 Hyperbaric Oxygen Chamber is a high-quality medical device with a Cabin Weight of 900KGS . It is designed to be user-friendly, with easy-to-use controls, and has a noise level of less than 63 DB, ensuring a comfortable and quiet therapy experience.

Overall, the O2EMENG EMY03 Hyperbaric Oxygen Chamber is a top choice for anyone in need of hyperbaric oxygen therapy. Whether you are a patient or a healthcare professional, this HBOT chamber is sure to exceed your expectations in terms of quality,



functionality, and ease of use.

Reasons to choose us:

1. Professionalism and authority:

As a manufacturer or brand of hyperbaric chamber, O2EMENG has rich industry experience and technology accumulation, and can provide professional and reliable hyperbaric oxygen therapy equipment.

Its products may have undergone stringent quality control and certification and comply with relevant industry standards to ensure safety and effectiveness in use.

2. Technological advancement:

O2EMENG hyperbaric chambers may adopt the latest technology and design concepts, such as air pressurization or pure oxygen pressurization technology, as well as advanced air supply and exhaust (oxygen) systems, air conditioning systems and control systems.

These technologies can ensure that the human body can safely and effectively inhale a high concentration of oxygen in a high-pressure environment for therapeutic purposes.

3. User feedback and reputation:

When choosing hyperbaric chambers, user feedback and reputation are also important reference factors, and O2EMENG hyperbaric chambers may be widely praised for their good treatment effect and user experience.

Support and Services:

The Hyperbaric Oxygen Chamber product comes with comprehensive technical support and services to ensure that customers have a seamless experience with our product. Our technical support team is available to assist with any queries or issues that may arise during installation, setup, or usage of the product. We also provide regular maintenance and repair services to keep the product in optimal

condition and extend its lifespan. Our services include:

1. Demand customization: according to your use environment and expectations, recommend the optimal core products and life aesthetics composition scheme.

2. Exclusive design: the advanced technology team and the top design team jointly, from the size and color material matching, the actual program customization landing.

3. Craftsmanship Manufacturing: The core components and other materials of hyperbaric chamber series products are manufactured in aerospace material factories.

4. Lifetime maintenance: provide perfect professional after-sales service guarantee

We are committed to providing the highest level of support and services to our customers to ensure their satisfaction with our Hyperbaric Oxygen Chamber product.

APPLICATION SCENARIO



MEDICAL INSTITUTION



OXYGEN BAR CLUB



ELDERLY CARE CENTER



FITNESS INDUSTRY



BEAUTY INDUSTRY



HOTEL



PLATEAU AREA



INDIVIDUAL FAMILY



WORK PLACE

HARD TYPE HYPERBARIC OXYGEN CHAMBER



Innovative technology, quality assurance



The cabin door is made of aviation PC material



Multiple safety technologies
One-button pressure relief protection system during power outage



AI intelligent voice control system



OUR CERTIFICATE



CHAMBER CAPSULE DESCRIPTION



1. Automatic digital display real-time temperature, humidity, oxygen purity and pressure.
2. Rotary timer system, use of time can be pressed in advance Multiple noise reduction technology.
3. Excellent Silence performance, the noise is $\leq 63\text{dB}$. Offer you a peaceful experience environment.
4. Double pressure relief device Spiral step less super smooth rotary valve. Spiral ergonomic hand-held design. Pressure relief is smoother and more convenient. Internal and external two-way integrated polycarbonate valve. One valve, two-way operation.
5. Automatic pressurize relief valves
Two automatic high precision constant pressure valve. Ensure the product is safe and effective by Aluminum alloy steel, high hardness, wear resistance.
6. Communicating system.
Intercom telephone is convenient for users to talk with operator outside.
7. Machine start-up protection system.
The chamber only can be turned on when the chamber door is completely closed, avoiding the danger of the door being left unlocked.
8. Automatic pressurize and depressurize system.
Based on the principle of "Cylinder (engine)" , Pressurize automatically after the door is completely closed.
9. Made by Aluminum alloy steel; With environmental spray paint;
Strong weight bearing capacity, Durable and wear-resistant; non-deformable; long lifespan.
10. In case of unexpected power failure, the door can be opened by emergency safety from inside the chamber.



Shenzhen Emeng Health Technology Co., LTD



+86 13760822077



kaylayi0612@foxmail.com



oxygenconcentratormedical.com

Room D, 7th Floor, Guang Long Building, 162 Ping Xin North Road, PingHu Street, LongGang District, ShenZhen