

## Medical/Therapeutic Medical Oxygen Concentrator with Flow Rate of 10Liter/Min and Noise Level Less Than 63 DB

Our Product Introduction

For products please visit us on [oxygenconcentratormedical.com](http://oxygenconcentratormedical.com)

### Basic Information

- Place of Origin: China
- Brand Name: Tylenolmed
- Certification: CE,ISO13485,SGS ,FCC,RoHs
- Model Number: ND-ZY10SL
- Minimum Order Quantity: Negotiable
- Price: Negotiable
- Packaging Details: OPP packing ,1pcs per Carton. Carton size:540\*370\*695mm,Gross Weight:21kgs.
- Delivery Time: 7-10 Working days
- Payment Terms: T/T, Western Union
- Supply Ability: 1000pcs



### Product Specification

- Voltage: 110V/50Hz; 220V/50Hz
- Chamber Oxygen Concentration: 28% 32%
- Function: Hyperbaric Oxygen Therapy
- Rating: 1500W
- Power Consumption: 1500w
- Cabin Weight: 1000KGS
- Usage: Medical/Therapeutic
- Chamber Colour: White /Black
- Highlight: **10LPM medical oxygen concentrator, quiet therapeutic oxygen machine, high-flow oxygen concentrator**

## Product Description

### Product Description:

The Medical Oxygen Concentrator is a cutting-edge device designed for Hyperbaric Oxygen Therapy, making it an essential tool for clinics and medical facilities. This advanced equipment operates effectively at two different voltages, providing versatility and convenience - 110V/50Hz and 220V/50Hz. With an input power of 700VA, it ensures efficient and reliable performance.

Weighing 1000KGS, the Medical Oxygen Concentrator is built to be sturdy and durable, guaranteeing long-term use and stability in clinical settings. Available in two classic colors - White and Black, this device seamlessly integrates into any environment while exuding a professional and sleek appearance.

As a top-tier Clinic Oxygen Producer, the Medical Oxygen Concentrator stands out for its exceptional quality and functionality. It serves as a vital tool for medical professionals, enabling them to administer Hyperbaric Oxygen Therapy with precision and accuracy. This Clinic Oxygen Producer is engineered to meet the demanding requirements of healthcare facilities, ensuring optimal performance and reliability. Moreover, as a leading Medical Oxygen Generator, this device plays a crucial role in providing a continuous and reliable supply of oxygen for therapeutic purposes. Its innovative design and advanced features make it a standout choice for medical practitioners seeking a dependable and efficient solution for oxygen therapy.

### Features:

Product Name: Medical Oxygen Concentrator

Usage: Medical/Therapeutic

Rating: 1500W

Noise Level: Less Than 63 DB

Input power: 700VA

Product Category: Medical Oxygen Concentrator

### Technical Parameters:

Chamber Oxygen Concentration	28% 32%
Rating	1500W
Product Name	O2EMENG Medical Therapy Machine 1.3ATA Hyperbaric Oxygen Chamber For Healing
Usage	Medical/Therapeutic
Product Category	Medical Oxygen Concentrator
Color	White, Black
Chamber Colour	White /Black
Function	Hyperbaric Oxygen Therapy
Input power	700VA
Voltage	110V/50Hz; 220V/50Hz

### Applications:

Tylenolmed Medical Oxygen Concentrator ND-ZY10SL is a high-quality healthcare oxygen separator that is suitable for a variety of application occasions and scenarios. With certifications including CE, ISO13485, SGS, FCC, and RoHs, this product ensures safety and reliability in medical environments.

Originating from China, the Tylenolmed Oxygen Concentration Machine is available for purchase with a negotiable minimum order quantity and price, making it accessible for different healthcare facilities and professionals. Accepted payment terms include T/T and Western Union, providing flexibility for procurement.

With a supply ability of 1000pcs and a delivery time of 7-10 working days, the Tylenolmed Oxygen Concentrator offers efficiency in meeting the needs of medical settings. The packaging details, including OPP packing and carton dimensions of 540\*370\*695mm, ensure secure transportation of the product.

The chamber color options of white or black allow customization based on preference or facility aesthetics. Operating at 1500W, this oxygen concentration machine consumes 1500w of power, providing reliable performance for hyperbaric oxygen therapy.

Model EMY01 from Tylenolmed is designed to be versatile and effective in various healthcare settings, offering a consistent oxygen supply for patients in need. Whether used in hospitals, clinics, or home care environments, this oxygen concentrator serves as a crucial tool in enhancing respiratory care.

### Support and Services:

Our Product Technical Support and Services for the Medical Oxygen Concentrator include:

- Initial product setup assistance

- Troubleshooting and diagnostic support
- Maintenance guidance and recommendations
- Software updates and upgrades
- User training and education
- Warranty coverage and service agreements

Please feel free to reach out to our support team for any assistance or inquiries regarding your Medical Oxygen Concentrator.

## Packing and Shipping:

### Product Packaging:

The Medical Oxygen Concentrator is carefully packaged to ensure safe delivery to our customers. It is securely wrapped in protective materials to prevent any damage during transportation. The packaging is designed to be durable and reliable, keeping the concentrator safe during shipping.

### Shipping:

Once your order is processed, the Medical Oxygen Concentrator will be shipped to your address using a trusted courier service. We strive to deliver the product to you in a timely manner to ensure you receive it promptly. You will receive a tracking number to monitor the status of your shipment and ensure a smooth delivery process.

## FAQ:

**Q: What is the brand name of this Medical Oxygen Concentrator?**

**A: The brand name is Tyleneolmed.**

**Q: What is the model number of this Medical Oxygen Concentrator?**

**A: The model number is ND-ZY10SL.**

**Q: What certifications does this Medical Oxygen Concentrator have?**

**A: It is certified with CE, ISO13485, SGS, FCC, and RoHs.**

**Q: Where is this Medical Oxygen Concentrator manufactured?**

**A: It is manufactured in China.**

**Q: What are the payment terms accepted for purchasing this Medical Oxygen Concentrator?**

**A: The accepted payment terms are T/T and Western Union.**



## CHAMBER CAPSULE DESCRIPTION



1. Automatic digital display real-time temperature, humidity, oxygen purity and pressure.
2. Rotary timer system, use of time can be pressed in advance Multiple noise reduction technology.
3. Excellent Silence performance, the noise is  $\leq 63\text{dB}$ . Offer you a peaceful experience environment.
4. Double pressure relief device Spiral step less super smooth rotary valve. Spiral ergonomic hand-held design. Pressure relief is smoother and more convenient. Internal and external two-way integrated polycarbonate valve. One valve, two-way operation.
5. Automatic pressurize relief valves  
Two automatic high precision constant pressure valve. Ensure the product is safe and effective by Aluminum alloy steel, high hardness, wear resistance.
6. Communicating system.  
Intercom telephone is convenient for users to talk with operator outside.
7. Machine start-up protection system.  
The chamber only can be turned on when the chamber door is completely closed, avoiding the danger of the door being left unlocked.
8. Automatic pressurize and depressurize system.  
Based on the principle of "Cylinder (engine)" , Pressurize automatically after the door is completely closed.
9. Made by Aluminum alloy steel; With environmental spray paint;  
Strong weight bearing capacity, Durable and wear-resistant; non-deformable; long lifespan.
10. In case of unexpected power failure, the door can be opened by emergency safety from inside the chamber.





## WHY ARE MORE AND MORE PEOPLE CHOOSING HYPERBARIC OXYGEN CHAMBERS?

In some parts of the body, the capillaries are so narrow that blood cells cannot pass through them, bound oxygen cannot reach them, and they can only rely on dissolved oxygen. The oxygen generator has little change in the dissolved oxygen in the blood, which is also the reason why the ordinary oxygen absorption effect in the plateau area is general.



## CONCLUSION



It is more feasible to increase the dissolved oxygen in the blood.  
The hyperbaric oxygen chamber can maintain the oxygen partial pressure within a safe and effective range by flexibly adjusting the oxygen concentration or pressure, and gently and effectively improve the oxygen delivery to tissues and organs by increasing dissolved oxygen.



# HARD TYPE HYPERBARIC OXYGEN CHAMBER



Innovative technology, quality assurance



The cabin door is made of aviation PC material

Multiple safety technologies  
One-button pressure relief protection system during power outage

AI intelligent voice control system

## OUR CERTIFICATE



## SHIPPING AND DELIVERY





## FAQ

### 1. How does hyperbarics work?

Hyperbarics is a technology in which the air pressure in the environment is increased. When a person's body is placed in a more pressure environment, it absorbs more oxygen molecules per volume of compressed air.

The body normally transports oxygen via the hemoglobin of the red blood cells. By increasing the air pressure, oxygen is then driven into the body's fluids, allowing a super-saturation of tissues and organs with oxygen. The increased pressure infuses the body with oxygen, even reaching injuries with damaged circulation.

### 2. What are the health benefits of hyperbaric oxygen chambers?

Studies have shown that hyperbaric oxygen chambers can treat hypoxic diseases, relieve fatigue, delay aging, prevent cardiovascular disease, improve sleep quality, accelerated wound repair and enhance immunity. We have customers serviced shows hyperbaric oxygen therapy makes very positive results.

### 3. How long is the treatment?

Usually, we recommend our clients receive the HBOT therapy 5-6 sessions per week, 60-90mins per session. Using frequency or duration can be adjusted according to individual circumstances.

### 4. Is HBOT safe?

Hyperbaric oxygen therapy is prescribed by a physician and performed under medical supervision. Although there are minor risks like all medical treatments, overall hyperbaric oxygen therapy is extremely safe.

### 5. Is there anyone who can not use our chamber?

Pregnant ladies, patients who catch a cold or Flu, or those who have rhinitis or otitis media are not recommended to use the chamber.

### 6. How do you feel inside a hyperbaric oxygen chamber?

When lying in the oxygen chamber, you will feel comfortable and at ease. Also, You may feel a little squeezing of the ears because of the pressure. It is quite normal and can be relieved by chewing or swallowing.

### 7. After bought the chamber, do I need to by any other things?

No, the price of our chamber includes the chamber capsule, all machines, namely the air compressor, oxygen concentrator, air cooler, and accessories like frame, tubes, etc, so you don't have to buy anything else.

### 8. Is it complex to install the chamber or change the filter?

Of course not, we will provide you with a video to show you how to install the chamber and change the filter.

We also support online services to help you with all the issues, so don't worry.

### 9. What is the general time warranty expected for ZHIKANG oxygen chambers?

We will provide one-year warranty and lifelong after-sales services.

If any quality problem/fault in material/design is under correct operation within one year,

If it's easy to fix, we will send new components for free and teach you how to repair them.

(Even after the one-year warranty period, you can also contact us freely and we will be glad to help you.)



**Shenzhen Emeng Health Technology Co., LTD**



+86 13760822077



kaylayi0612@foxmail.com



oxygenconcentratormedical.com

Room D, 7th Floor, Guang Long Building, 162 Ping Xin North Road, PingHu Street, LongGang District, ShenZhen