



Shenzhen Emeng Health Technology Co., LTD

oxygenconcentratormedical.com

## 28% 32% Chamber Oxygen Concentration Healthcare Oxygen Separator for Oxygen Therapy Equipment

### Our Product Introduction

For products please visit us on [oxygenconcentratormedical.com](http://oxygenconcentratormedical.com)

#### Basic Information

- Place of Origin: China
- Brand Name: Tylenolmed
- Certification: CE,ISO13485,SGS ,FCC,RoHs
- Model Number: ND-ZY10SL
- Minimum Order Quantity: Negotiable
- Price: Negotiable
- Packaging Details: OPP packing ,1pcs per Carton. Carton size:540\*370\*695mm,Gross Weight:21kgs.
- Delivery Time: 7-10 Working days
- Payment Terms: T/T, Western Union
- Supply Ability: 1000pcs



#### Product Specification

- Noise Level: Less Than 63 DB
- Voltage: 110V/50Hz; 220V/50Hz
- Rating: 1500W
- Cabin Weight: 1000KGS
- Product Category: Medical Oxygen Concentrator
- Product Name: O2EMENG Medical Therapy Machine  
1.3ATA Hyperbaric Oxygen Chamber For Healing
- Power Consumption: 1500w
- Usage: Medical/Therapeutic
- Highlight: **medical oxygen concentrator with 28-32% concentration**  
,  
**healthcare oxygen separator for therapy equipment**

## Product Description

### Product Description:

The O2EMENG Medical Therapy Machine 1.3ATA Hyperbaric Oxygen Chamber For Healing is a top-of-the-line Medical Oxygen Concentrator designed to provide high-quality oxygen therapy for medical purposes. This advanced Healthcare Oxygen Separator is equipped with a Flow Rate of 10 liters per minute, ensuring a steady and reliable supply of oxygen for patients in need.

Powered by a 700VA Input power, this Medical Oxygen Generator offers efficient performance and is suitable for use in various healthcare settings such as hospitals, clinics, and home care environments. The Oxygen Concentration Machine is designed to deliver a consistent oxygen flow to support patients with respiratory conditions or those requiring supplemental oxygen therapy.

The O2EMENG Medical Therapy Machine features a sleek and modern design, available in two elegant Chamber Colour options: White or Black. This allows users to choose a style that best fits their preferences or blends seamlessly with the existing decor of their medical facility or home environment.

Despite its powerful capabilities, the O2EMENG Medical Therapy Machine is surprisingly lightweight, with a Cabin Weight of 1000KGS. This makes it easy to move and transport as needed, providing flexibility and convenience for healthcare professionals or caregivers who may need to relocate the device within a facility or between different patient rooms.

Whether used in a clinical setting or at home, the O2EMENG Medical Therapy Machine offers a reliable and efficient solution for oxygen therapy needs. Its high Flow Rate, durable construction, and user-friendly design make it a valuable asset for patients requiring oxygen supplementation for various medical conditions.

### Features:

Product Name: Medical Oxygen Concentrator

Voltage: 110V/50Hz; 220V/50Hz

Usage: Medical/Therapeutic

Noise Level: Less Than 63 DB

Product Category: Medical Oxygen Concentrator

Power Consumption: 1500w

### Technical Parameters:

Chamber Oxygen Concentration	28% 32%
Product Name	O2EMENG Medical Therapy Machine 1.3ATA Hyperbaric Oxygen Chamber For Healing
Voltage	110V/50Hz; 220V/50Hz
Cabin Weight	1000KGS
Product Category	Medical Oxygen Concentrator
Usage	Medical/Therapeutic
Power Consumption	1500w
Noise Level	Less Than 63 DB
Flow Rate	10Liter/Min
Input power	700VA

### Applications:

Tylenolmed Medical Oxygen Concentrator ND-ZY10SL is a top-of-the-line Healthcare Oxygen Separator designed to provide high-quality oxygen therapy in various medical settings. With certifications including CE, ISO13485, SGS, FCC, and RoHs, this Oxygen Concentration Machine guarantees safety and efficiency in delivering pure oxygen to patients.

Originating from China, this Clinic Oxygen Producer offers flexibility in terms of Minimum Order Quantity and Price, making it a cost-effective solution for hospitals, clinics, and healthcare facilities. Payment Terms accepted include T/T and Western Union, ensuring convenience for buyers.

With a Supply Ability of 1000pcs and a short Delivery Time of 7-10 Working days, the Tylenolmed Oxygen Concentrator ND-ZY10SL is readily available to meet the demands of medical institutions. The Packaging Details reveal a sturdy OPP packing with 1pcs per Carton, ensuring safe transport. The Carton size is 540\*370\*695mm with a Gross Weight of 21kgs.

Available in White and Black colors, this device operates quietly with a Noise Level of Less Than 63 DB, ensuring a peaceful environment for patients. The Chamber Oxygen Concentration ranges from 28% to 32%, providing the necessary oxygen levels for therapeutic use.

With a Power Consumption of 1500w, this Oxygen Concentrator is energy-efficient and reliable.

Weighing 1000KGS, the Tylenolmed Medical Oxygen Concentrator ND-ZY10SL is a robust and durable solution for medical oxygen therapy needs. Its advanced features and specifications make it a preferred choice for healthcare professionals seeking a reliable and efficient Oxygen Concentration Machine.

**Support and Services:**

Our Product Technical Support and Services for the Medical Oxygen Concentrator include: -----

- Assistance with setup and installation of the oxygen concentrator
- Troubleshooting guidance for any operational issues
- Maintenance and servicing to ensure optimal performance
- Training on how to use the device safely and effectively

**Packing and Shipping:**

Product Name: Medical Oxygen Concentrator  
Description: This medical oxygen concentrator is designed to provide a continuous flow of oxygen for individuals with respiratory conditions.  
Features:  
Adjustable oxygen flow settings  
Compact and portable design  
Quiet operation  
Package Contents:  
1 Medical Oxygen Concentrator  
1 Power Cord  
1 User Manual  
Shipping Information:  
Shipping Method: Standard Shipping  
Shipping Time: 3-5 business days  
Shipping Cost: Calculated at checkout

**FAQ:**

Q: What is the brand name of this Medical Oxygen Concentrator? -----  
A: The brand name is Tylenolmed.  
Q: What is the model number of this Medical Oxygen Concentrator?  
A: The model number is ND-ZY10SL.  
Q: What certifications does this Medical Oxygen Concentrator have?  
A: It is certified with CE, ISO13485, SGS, FCC, and RoHs.  
Q: Where is this Medical Oxygen Concentrator manufactured?  
A: It is manufactured in China.  
Q: What are the payment terms for purchasing this Medical Oxygen Concentrator?  
A: The payment terms accepted are T/T and Western Union.

O2EMENG

# HYPERBARIC OXYGEN CHAMBER

- IMPROVE BRAIN BLOOD CIRCULATION
- SPEED UP METABOLISM
- HIGH PRESSURE OXYGEN INHALATION



## CHAMBER CAPSULE DESCRIPTION



1. Automatic digital display real-time temperature, humidity, oxygen purity and pressure.
2. Rotary timer system, use of time can be pressed in advance Multiple noise reduction technology.
3. Excellent Silence performance, the noise is  $\leq 63\text{dB}$ . Offer you a peaceful experience environment.
4. Double pressure relief device Spiral step less super smooth rotary valve. Spiral ergonomic hand-held design. Pressure relief is smoother and more convenient. Internal and external two-way integrated polycarbonate valve. One valve, two-way operation.
5. Automatic pressurize relief valves  
Two automatic high precision constant pressure valve. Ensure the product is safe and effective by Aluminum alloy steel, high hardness, wear resistance.
6. Communicating system.  
Intercom telephone is convenient for users to talk with operator outside.
7. Machine start-up protection system.  
The chamber only can be turned on when the chamber door is completely closed, avoiding the danger of the door being left unlocked.
8. Automatic pressurize and depressurize system.  
Based on the principle of "Cylinder (engine)" , Pressurize automatically after the door is completely closed.
9. Made by Aluminum alloy steel; With environmental spray paint;  
Strong weight bearing capacity, Durable and wear-resistant; non-deformable; long lifespan.
10. In case of unexpected power failure, the door can be opened by emergency safety from inside the chamber.





## WHY ARE MORE AND MORE PEOPLE CHOOSING HYPERBARIC OXYGEN CHAMBERS?

In some parts of the body, the capillaries are so narrow that blood cells cannot pass through them, bound oxygen cannot reach them, and they can only rely on dissolved oxygen. The oxygen generator has little change in the dissolved oxygen in the blood, which is also the reason why the ordinary oxygen absorption effect in the plateau area is general.



## CONCLUSION



It is more feasible to increase the dissolved oxygen in the blood. The hyperbaric oxygen chamber can maintain the oxygen partial pressure within a safe and effective range by flexibly adjusting the oxygen concentration or pressure, and gently and effectively improve the oxygen delivery to tissues and organs by increasing dissolved oxygen.



# HARD TYPE HYPERBARIC OXYGEN CHAMBER



Innovative technology, quality assurance



The cabin door is made of aviation PC material

Multiple safety technologies  
One-button pressure relief protection system during power outage

AI intelligent voice control system



## OUR CERTIFICATE



## SHIPPING AND DELIVERY



## FAQ

### 1. How does hyperbarics work?

Hyperbarics is a technology in which the air pressure in the environment is increased. When a person's body is placed in a more pressure environment, it absorbs more oxygen molecules per volume of compressed air.

The body normally transports oxygen via the hemoglobin of the red blood cells. By increasing the air pressure, oxygen is then driven into the body's fluids, allowing a super-saturation of tissues and organs with oxygen. The increased pressure infuses the body with oxygen, even reaching injuries with damaged circulation.

### 2. What are the health benefits of hyperbaric oxygen chambers?

Studies have shown that hyperbaric oxygen chambers can treat hypoxic diseases, relieve fatigue, delay aging, prevent cardiovascular disease, improve sleep quality, accelerated wound repair and enhance immunity. We have customers serviced shows hyperbaric oxygen therapy makes very positive results.

### 3. How long is the treatment?

Usually, we recommend our clients receive the HBOT therapy 5-6 sessions per week, 60-90mins per session. Using frequency or duration can be adjusted according to individual circumstances.

### 4. Is HBOT safe?

Hyperbaric oxygen therapy is prescribed by a physician and performed under medical supervision. Although there are minor risks like all medical treatments, overall hyperbaric oxygen therapy is extremely safe.

### 5. Is there anyone who can not use our chamber?

Pregnant ladies, patients who catch a cold or Flu, or those who have rhinitis or otitis media are not recommended to use the chamber.

### 6. How do you feel inside a hyperbaric oxygen chamber?

When lying in the oxygen chamber, you will feel comfortable and at ease. Also, You may feel a little squeezing of the ears because of the pressure. It is quite normal and can be relieved by chewing or swallowing.

### 7. After bought the chamber, do I need to by any other things?

No, the price of our chamber includes the chamber capsule, all machines, namely the air compressor, oxygen concentrator, air cooler, and accessories like frame, tubes, etc, so you don't have to buy anything else.

### 8. Is it complex to install the chamber or change the filter?

Of course not, we will provide you with a video to show you how to install the chamber and change the filter.

We also support online services to help you with all the issues, so don't worry.

### 9. What is the general time warranty expected for ZHIKANG oxygen chambers?

We will provide one-year warranty and lifelong after-sales services.

If any quality problem/fault in material/design is under correct operation within one year,

If it's easy to fix, we will send new components for free and teach you how to repair them.

(Even after the one-year warranty period, you can also contact us freely and we will be glad to help you.)



**Shenzhen Emeng Health Technology Co., LTD**



+86 13760822077



kaylayi0612@foxmail.com



oxygenconcentratormedical.com

Room D, 7th Floor, Guang Long Building, 162 Ping Xin North Road, PingHu Street, LongGang District, ShenZhen